

Product datasheet for **SC322054**

DLK2 (NM_023932) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DLK2 (NM_023932) Human Untagged Clone
Tag:	Tag Free
Symbol:	DLK2
Synonyms:	DLK-2; EGFL9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_023932, the custom clone sequence may differ by one or more nucleotides

ATGCCAGCGGCTGCCGCTGCCTGCATCTCGTGTGCCTGTTGTGATTCTGGGGCTCCCGTCAAGCTG
 TCCGAGCCGATGACTGCAGCTCCCACTGTGACCTGGCCACGGCTGCTGTGCACCTGACGGCTCCTGCAG
 GTGTGACCCGGGCTGGGAGGGGCTGCACTGTGAGCGCTGTGTGAGGATGCCTGGCTGCCAGCAGGTACC
 TGCCACCAGCCATGGCAGTGCATCTGCCACAGTGGCTGGGCAGGCAAGTTCTGTGACAAAGATGAACATA
 TCTGTACCACGCACTCCCTGCCAGAATGGAGGCCAGTGCATGTATGACGGGGCGGTGAGTACCATTG
 TGTGTGCTTACCAGGCTTCCATGGGCGTGAAGTGCAGGCGCAAGGCTGGACCTGTGAACAGGCAGGCTCC
 CCATGCCGCAATGGCGGCGAGTCCAGGACGACAGGGCTTTGCTCTCACTTCACGTGCCGCTGCTTGG
 TGGGCTTTGTGGGTGCCGCTGTGAGGTAATGTGGATGACTGCCTGATGCGGCCTGTGCTAACGGTGC
 CACCTGCCTTGACGGCATAAACCGCTTCTCCTGCCTCTGTCTGAGGGCTTTGCTGGACGCTTCTGCACC
 ATCAACCTGGATGACTGTGCCAGCCGCCATGCCAGAGAGGGGCCGCTGTGGGACCGTGTCCAGGACT
 TCGACTGCCTCTGCCCCAGTGGCTATGGTGGCAAGACCTGTGAGCTTGTCTTACCTGTCCAGACCCCC
 AACCACAGTGGACACCCCTCTAGGGCCACCTCAGCTGTAGTGGTACCTGCCACGGGGCCAGCCCCCAC
 AGCGCAGGGGCTGGTCTGCTGCGGATCTCAGTGAAGGAGTGGTGGGAGGCAAGAGGCTGGGCTAGGTG
 AGCCTAGCTTGGTGGCCCTGGTGGTGTGGGGCCCTCACTGCTGCCCTGGTCTGGCTACTGTGTTGCT
 GACCTGAGGGCTGGCGCCGGGTGTCTGCCCCCTGGACCTGTTGCTACCTGCCCCACACTATGCT
 CCAGCGTCCAGGACGAGGTGTCAGGTTAGCATGCTGCCAGCAGGGCTCCCTGCCAGCTGACTTGC
 CCCCTGAGCTGGAAAGACCACAGCACTGTGA

Restriction Sites: Please inquire

ACCN: NM_023932



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_023932.2</u> , <u>NP_076421.2</u>
RefSeq Size:	1574 bp
RefSeq ORF:	1152 bp
Locus ID:	65989
UniProt ID:	<u>Q6UY11</u>
Cytogenetics:	6p21.1
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	Regulates adipogenesis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) differs in the 5' UTR, compared to variant 2. Both variants 1 and 2 encode the same isoform (a).