

## Product datasheet for **SC322080**

### OGFRL1 (NM\_024576) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	OGFRL1 (NM_024576) Human Untagged Clone
Tag:	Tag Free
Symbol:	OGFRL1
Synonyms:	dj331H24.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_024576, the custom clone sequence may differ by one or more nucleotides

```
ATGGGCAACCTGCTCGGCGGGGTGAGCTCCGCGAGCCACCACCGTGGAGGACTGCGACTCCACCTGGC
AGACCGACTCGGAGCCGAGCCGAAGAACCAGGGCCGGGCGGCGGCGAGCGAGGGCCCGGGGAGGAGTC
CGAGCAGCCCGCGCAGCCCCGGAGCAAGCCGGGCGGGCGCCGCGCCAGCCCCGCGCCGAGGAC
GCCGAGGCGGGGCGCCGAGCAGGGTGGTATTCCACTGAAGCAACTGCCAAACCAAGAGAAGTTTTT
ATGCTGCCAGGGATTTGTACAAGTACCGACACCAGTACCCAACTTCAAAGATATCCGATATCAAAATGA
CTTGAGCAATCTTCGTTTTATAAGAATAAAATCCATTCAAGCCAGATGGTGTTCACATTGAAGAAGTT
CTAAGTAAATGGAAAGGAGATTATGAAAACTGGAGCACAAACCACACTTACATTCATGGCTTTTCCCC
TGAGAGAACAAAGGCTTGAACCTTATGCCAAAGAACTAATACATATGAAATTGAGGAATCAAAAAAAC
AAAAGAAGCAATTAGAAGATTCCTCGCTTATAAAATGATGCTAGAATTTTTTGAATAAAAGTACTGACT
GATAAACTGAAATGTTGCTCGGGCTGTTAACTGGCAGGAAAGATTTTCAGCATCTGAATGAGTCCCAGC
ACAATATTTAAGAATCACTCGTATTCTTAAAAGCCTTGGTGAGCTTGGATATGAAAGTTTTAAATCTCC
TCTTGAAAAATTTATTCTTCATGAAGCTTGTGGAGAATACTATTCCCAATATTAAGCAGAGTGCTCTA
GAGTATTTGTTTATAACAATTAGAGACAGAAGAGAAAGGAGAAAGCTCCTGCGGTTCGCCAGAAACACT
ACACGCCCTCAGAGAACCTTATCTGGGGACCGCTCGAAAAGAACAGTCGGAGGGAAGCAAGGCCAGAA
AATGCTTCCCCTCTCGCTCCAGTCATAACAGTCAAACCTTCTATGCACAAAAAGCCAAAGGACTCCAAA
AATTCTCCTCAGCTGTTCAATTTAAATAGCAAAACAGCTGAAGACAAAAAAGTGGCACCAAAAGAGGCTG
TGAAGAGACAGACAGGCCAGCCAGCCAGACCCAGCAATGAAGCTGCCAAGCCAAGAAATACAGAGAAGGA
CAGTAATGCTGAGAACATGAATTCTCAACCTGAGAAAACAGTTACTACTCCACAGAAAAAAGGAGAGT
GTATCTCCTGAGAATAACGAAGAAGGTGGAAATGATAACCAAGACAATGAAAATCCTGGAATACAAATT
GCCATGATGTTGTACTAGTACAGTGA
```

Restriction Sites: Please inquire



[View online >](#)

<b>ACCN:</b>	NM_024576
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_024576.3</a></u> , <u><a href="#">NP_078852.3</a></u>
<b>RefSeq Size:</b>	1755 bp
<b>RefSeq ORF:</b>	1356 bp
<b>Locus ID:</b>	79627
<b>UniProt ID:</b>	<u><a href="#">Q5TC84</a></u>
<b>Cytogenetics:</b>	6q13
<b>Domains:</b>	OGFr_N
<b>Protein Families:</b>	Druggable Genome