

Product datasheet for **SC322082**

SATL1 (NM_001012980) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SATL1 (NM_001012980) Human Untagged Clone
Tag:	Tag Free
Symbol:	SATL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for SC322082
 AGCGGACTGGGCCAACCAAGTCCTGAGGCCAACCAACATGAGTCCACCAGGCATGTGGCAA
 CCAGGCGTGCAACAACCAGGCATCAGCCAGCAAGTCCCAGCCACCCAGACATGAGTCAA
 CCAGGCATGAGCCAGCAAGTCCCAGCCAACCAAGGCATAAGGCCAACCAAGACTAGCCAA
 TCATGTAAGAACCAACAGACATGAGCCAACCAAGCGCAAACCAATCAAGTTTATCAGAT
 TCCAACCAAACAGGTATAATCCAGCCAAGCCAAAGCTTACTTGGTATGAACCAAAATGGAC
 ATGAACCAATGGAGTGAAGCCTATATGAAATGAACCAAGTGGACATGAAACAACCAAGC
 ATGAGCCAAGCTGGCATGAGGCAATCAGGTACAAACCTACCAGACATAAACCAACCTGGC
 ATGAAACAACCAGGCACATGGCAATTAGGTAGGAGCCAACCAGGCATGTGGCCACAAAGC
 CTGAGCGAACTAGTCCTGAGTGAAGCAAGCATAAGCCAACCAGGTCCACCAGCAACGAGCC
 CCAAGCCAATCAGGCCCCAGACAATCAAGCAGGAGCCAAGCAGGCACAAACCAATCAGGT
 ATAAGCCAACCAGTGATGTGGCAACTAGACATGAGACAGTCAGGTGGGAGCCAACCAAGC
 ATGAGACAAGTAGGCACCAGCCAATCAGGCACAAGCCAAATAGGCATGAGCCAACCAGGC
 ACATGGCAAACAGGCCTGAGCCAACCAGTCCCGAGGCAACCAACAAGAGTCCACCAGGC
 ATGTGGCAACGAGGCATGTGGCAACCAGGCATGAGCCAGCAAGTCCCAGCCAACAGGC
 ATGAGACAACCAGGCATGAGCCAATCAAGTAAGAACCAACAGGCATGAGCCATCCAGGC
 AGGGGCCAACCAAGGCATATGGGAACCGGGCCGAGTCAGCCAGGCCTGAGCCAACAAAGC
 CTGAACCAATTAGTGCTGAGCCAACCAGGCCTGAGTCAACCAGGCAGGAGCCAACCAAGT
 GTGAGCCAATGGGCATGAGGCAACAAGCATGGATTACTTTCAAATAAGACATGCGAGAG
 GCTGGAGACTGCCAGAAATTTTGCAGTATTAAAGAATTGGCTGCCTGTGAAAACATG
 CTAGATGCAATGGAGTTAACAGCAGCTGATTTACTCAGAGATGGCTTTGGGGACAATCCC
 CTTTTCTACTGCCTGATTGCAGAAGTAAACGATCAACAAAAACCATCAGGCAAACTGACT
 GTTGGATTTGCCATGACTACTTTACATACGACTCATGGACTGGCAAGGTACTTTACCTA
 GAGGACTTTTTATGTCACACAAGCTTACCAAGATAGCCATCACAACCTCAATGTAACATGCAT
 GCACTTTCTTGTGTCATTTGGAACCAGGCTTCTATCACTACTATACTAGTCGAGGGGC
 TTTAGACCTTCTCTGAGGAGGGCTGGCATCTCTCAGGTTTAAACAGAGAAGAACTCCT
 GGACATGGCATGGGAAGAATGAAAGAGGCTTTGTAACCACTCATCCCAACACAGCCTCCA
 ACATTTGACTGTCACCTGAGTCTAACAGCAGTTTACTTTGTTGTACAATACAGCAAGTG
 ATCATCAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire

ACCN: NM_001012980

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001012980.1</u> , <u>NP_001012998.1</u>
RefSeq Size:	1616 bp
RefSeq ORF:	1338 bp
Locus ID:	340562
UniProt ID:	<u>Q86VE3</u>
Cytogenetics:	Xq21.1