

Product datasheet for **SC322128**

ACBD7 (NM_001039844) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ACBD7 (NM_001039844) Human Untagged Clone
Tag: Tag Free
Symbol: ACBD7
Synonyms: bA455B2.2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC (PS100020)
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for SC322128
 TGAGCACCGGACGACAACAGCAGCAGCAACAACGCGGCAGCCGCCACCATGGCCCTGCAG
 GCTGATTTTGACAGGGCTGCAGAAGATGTGAGGAAGCTGAAAGCAAGACCAGATGATGGA
 GAACTGAAAGAAGCTCTATGGGCTTTACAACAAGCAATAGTTGGAGACATTAATATTGCG
 TGTCAGGAATGCTAGATTTGAAAGGCAAAGCCAAATGGGAAGCATGGAACCTCAAAAAA
 GGGTTGTCGACGGAAGATGCGACGAGTGCCTATATTTCTAAAGCAAAGGAGCTGATAGAA
 AAATACGGAATTTAGAATACAGCATATGAGGAATTTTTCTTTTGAAGACTTCCAATGC
 TATCATGACCTAACATTTAGAGGGAGAGGCATACTGTTAACTTGATGTATCATGTATATT
 TTTGCTATTAGCATGAGCTGTAGTACTTAGAAGTATACTGAAACTACATAGTTAACTCAA
 TTGTAAGTCTTAAACCAGGTGTCTTTAAAAGTCTTTTTAAGAAAATGTGGCACATATA
 CACCATGGAATACTATGGAGCCATAAAAAGGGATGAGTTCATGTCCTTTGCAGGGACATG
 GATGAAGCTGGAACCATCATTCTCAGCAAATATCACAAGATCAGAAAACCAAGCATGT
 TCTCACTCATAAGTGGGAGTTGAACAATGAGAACACATGGACACAGGGAAGGGAACATCA
 CAACTGGGGCCTGTCGGGGACTAGGGGTTAGGGGAGGGATAACATTAGGAGACATACC
 TGATGTAGGTGACCGGTTGATGGGTACAGCAAACCACCATGGCAGATGTATACCTATGTA
 ACAAAGTGCACGTTCTGCATATGTATCCAGAACTTAAAGTATATAAAAAAAATTTCTT
 TTAGAAAAGTATTTTTATTTTTATAGAGTTAGGGAGTACAGTTGCAGTTTTTGTGTTGTT
 TGTTTGTTTTTCTTGAGGAGTCTCACTCCAGCCTGGGCGACAAAACAAGACTCCGTCTC
 AAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire
ACCN: NM_001039844



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001039844.2, NP_001034933.1</u>
RefSeq Size:	3383 bp
RefSeq ORF:	267 bp
Locus ID:	414149
UniProt ID:	<u>Q8N6N7</u>
Cytogenetics:	10p13
Gene Summary:	Binds medium- and long-chain acyl-CoA esters.[UniProtKB/Swiss-Prot Function]