

Product datasheet for SC322232

TMEM106C (NM_024056) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM106C (NM_024056) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM106C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC322232

```
GACATGACACCAGTGGCATATCACGGCCATGGGGTCTCAGCATTCCGCTGCTGCTCGCCC
CTCCTCCTGCAGGCGAAAGCAAGAAGATGACAGGGACGGTTTGCTGGCTGAACGAGAGCA
GGAAGAAGCCATTGCTCAGTTCCCATATGTGGAATTCACCGGGAGAGATAGCATCACCTG
TCTCACGTGCCAGGGGACAGGCTACATTCCAACAGAGCAAGTAAATGAGTTGGTGGCTTT
GATCCCACACAGTGATCAGAGATTGCGCCCTCAGCGAACTAAGCAATATGTCCTCCTGTC
CATCCTGCTTTGTCTCCTGGCATCTGGTTTGGTGGTTTTCTTCTGTTTCCGCATTCACT
CCTTGTGGATGATGACGGCATCAAAGTGGTGAAAGTCACATTTAATAAGCAAGACTCCCT
TGTAATTTCCACCATCATGGCCACCCTGAAAATCAGGAACTCAAACCTTCTACACGGTGGC
AGTGACCAGCCTGCCAGCCAGATTCACTACATGAACACAGTGGTCACTACATATGTGAC
TACTAACGTCTCCCTTATCCACCTCGGAGTGAGCAACTGGTGAATTTTACCGGGAAGGC
CGAGATGGGAGGACCGTTTTCTATGTGTACTTCTTCTGCACGGTACCTGAGATCCTGGT
GCACAACATAGTGATCTTCATGCGAACTTCAGTGAAGATTTTATACATTGGCCTCATGAC
CCAGAGCTCCTTGGAGACACATCACTATGTGGATTGTGGAGGAAATCCACAGCTATTTA
ACAACCTGCTATTGGTTCTTCCACACAGCGCCTGTAGAAGAGAGCACAGCATATGTTCCCA
AGGCCTGAGTTCTGGACCTACCCCCACGTGGTGAAGCAGAGGAGGAATTGGTCACTTA
ACTCCCAGCAAACATCCTCCTGCCACTTAGGAGGAAACACCTCCCTATGGTACCATTAT
GTTTCTCAGAACCAGCAGAATCAGTGCCTAGCCTGTGCCAGCAAATAGTTGGCACTCAA
TAAAGATTTGCAGAATTTAATACAGATCTTTTCAGCTGTTCTTAGGGCATTATAAATGGA
AATCATAAACGTGGTTCTAGGTTATCAAACCATGGAGTGATGTGGAGCTAGGATTGTGAGT
GACCTGCAGGCCATTATCAGTGCCTCATCTGTGAGAAGTGGCAGCAGAGAGGGACCATC
CAAATACCTAAGAGAAAACAGACCTAGTCAGGATATGAATTTGTTTCAGCTGTTCCCAAA
GGCCTGGGAGCTTTTTGAAAAGAAAAGAAAAGTGTGTTGGCTTTTTTTTTTTTTTAGAA
AGTTAGAATTGTTTTTACCAAGAGTCTATGTGGGGCTTGATTACCCTTCATCCATTGGC
TGGAACATGGATTGGGGATTTGATAGAAAAATAAACCTGCTTTTTGATTCAAAAA
AAAAAAAA
```

Restriction Sites: Please inquire



[View online »](#)

ACCN:	NM_024056
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024056.2 , NP_076961.1
RefSeq Size:	1449 bp
RefSeq ORF:	753 bp
Locus ID:	79022
UniProt ID:	Q9BVX2
Cytogenetics:	12q13.11
Protein Families:	Transmembrane