

Product datasheet for SC322314

AGA (NM_000027) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGA (NM_000027) Human Untagged Clone
Tag:	Tag Free
Symbol:	AGA
Synonyms:	AGU; ASRG; GA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for SC322314
GATCGCTGGCTGCCGGGACTTTTCTCGCGCTGGTCTCTTCGGTGGTCAGGGATGGCGCGG
AAGTCGAACTTGCCGTGTGCTTCTCGTGCCGTTTCTGCTCTGCCAGGCCCTAGTGCCTGC
TCCAGCCCTCTGCCCTGGTCGTCAACACTTGGCCCTTTAAGAATGCAACCGAAGCAGCG
TGGAGGGCATTAGCATCTGGAGGCTCTGCCCTGGATGCAGTGGAGAGCGGCTGTGCCATG
TGTGAGAGAGAGCAGTGTGACGGCTCTGTAGGCTTTGGAGGAAGTCTGATGAACCTGGA
GAAACCACACTAGATGCCATGATCATGGATGGCACTACTATGGATGTAGGAGCAGTAGGA
GATCTCAGACGAATTAATAATGCTATTGGTGTGGCAGGAAAGTACTGGAACATAACA
CACACACTTTTAGTAGGAGAGTCAGCCACCACATTTGCTCAAAGTATGGGGTTTATCAAT
GAAGACTTATCTACCAGTCTTCAAGCTCTTCAATTCAGATTGGCTTGCTCGGAATTGC
CAGCCAAATTATTGGAGGAATGTTATACCAGATCCCTCAAATACTGCGGACCCTACAAA
CCACCTGGTATCTTAAAGCAGGATATTCTATCCATAAAGAAACAGAAGATGATCGTGGT
CATGACACTATTGGCATGGTTGTAATCCATAAGACAGGACATATTGCTGCTGGTACATCT
ACAAATGGTATAAAAATCAAATACATGGCCGTGTAGGAGACTCACCAATACCTGGAGCT
GGAGCCTATGCTGACGATACTGCAGGGGCGAGCCGAGCCACTGGGAATGGTGATATATTG
ATGCGCTTCTGCCAAGCTACCAAGCTGTAGAATACATGAGAAGAGGAGAAGATCCAACC
ATAGCTTGCCAAAAGTGATTTCAAGAATCCAGAAGCATTTTCCAGAATCTTTGGGGCT
GTTATATGTGCCAATGTGACTGGAAGTTACGGTGTGCTTGAATAAACTTTCAACATTT
ACTCAGTTTAGTTTCATGGTTTATAATTCCGAAAAAATCAGCCAAGTGGGAAAAAGTG
GACTGCATCTAATCCATCTTTACTGTCAACATCTGTATTTAAAAAAAAAAAAAAAAAAAA

Restriction Sites:	Please inquire
ACCN:	NM_000027



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_000027.2 , NP_000018.1
RefSeq Size:	2041 bp
RefSeq ORF:	1041 bp
Locus ID:	175
UniProt ID:	P20933
Cytogenetics:	4q34.3
Domains:	Asparaginase_2
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Lysosome, Other glycan degradation
Gene Summary:	<p>This gene encodes a member of the N-terminal nucleophile (Ntn) hydrolase family of proteins. The encoded preproprotein is proteolytically processed to generate alpha and beta chains that comprise the mature enzyme. This enzyme is involved in the catabolism of N-linked oligosaccharides of glycoproteins. It cleaves asparagine from N-acetylglucosamines as one of the final steps in the lysosomal breakdown of glycoproteins. Mutations in this gene are associated with the lysosomal storage disease aspartylglycosaminuria that results in progressive neurodegeneration. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is subject to proteolytic processing. [provided by RefSeq, Nov 2015]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).</p>