

Product datasheet for **SC322322**

SPINT2 (NM_021102) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPINT2 (NM_021102) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPINT2
Synonyms:	DIAR3; HAI-2; HAI2; Kop; PB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for SC322322
 CGACCCACGCGTCCGGCGACCTCCGCGCCTGGGAGGTGTAGCGCGGCTCTGAACGCGCT
 GAGGGCCGTTGAGTGTGCGAGGCGGCGAGGGCGCGAGTGAGGAGCAGACCCAGGCATCGC
 GCGCCGAGAAGGCCGGGCGTCCCCACACTGAAGTCCGGAAAGGCGACTTCCGGGGGCTT
 TGGCACCTGGCGGACCTCCCGGAGCGTCGGCACCTGAACGCGAGGCGCTCCATTGCGCG
 TGCGCGTTGAGGGGCTTCCCGCACCTGATCGCGAGACCCCAACGGCTGGTGGCGTCGCT
 GCGCGTCTCGGCTGAGCTGGCCATGGCGCAGCTGTGCGGGCTGAGGCGGAGCCGGGCGTT
 TCTCGCCCTGCTGGGATCGCTGCTCCTCTCTGGGGTCTTGGCGGCGACCGAGAACGCAG
 CATCCACGACTTCTGCCTGGTGTGCGAAGGTGGTGGGCAGATGCCGGGCCTCCATGCCTAG
 GTGGTGTGACAATGTCACTGACGGATCCTGCCAGCTGTTTGTGTATGGGGCTGTGACGG
 AAACAGCAATAATTACCTGACCAAGGAGGAGTGCCTCAAGAAATGTCCACTGTACAGAG
 GAATGCCACGGGTGACCTGGCCACCAGCAGGAATGCAGCGGATTCTCTGTCCCAAGTGC
 TCCAGAAAGCAGGATTCTGAAGACCACTCCAGCGATATGTTCAACTATGAAGAATACTG
 CACCGCCAACGCAGTCACTGGGCTTGGCGTGCATCCTTCCACGCTGGTACTTTGACGT
 GGAGAGGAACCTGCAATAACTTCATCTATGGAGGCTGCCGGGCAATAAGAACAGCTA
 CCGCTCTGAGGAGGCTGCATGCTCCGCTGCTTCCGCCAGCAGGAGAATCTCCCTGCC
 CCTTGGCTCAAAGGTGGTCTTCTGGCGGGCTGTTCTGTGATGGTGTGATCCTCTTCT
 GGGAGCCTCCATGGTCTACCTGATCCGGGTGGCACGGAGGAACCAGGAGCGTGCCCTGCG
 CACCGTCTGGAGTCCGGAGATGACAAGGAGCAGCTGGTGAAGAACACATATGTCCTGTG
 ACCGCCCTGTGCGCAAGAGGACTGGGGAAGGGAGGGGAGACTATGTGTGAGCTTTTTTTA
 AATAGAGGGATTGACTCGGATTTGAGTGATCATTAGGGCTGAGGTCTGTTTCTCTGGGAG
 GTAGGACGGCTGCTTCTGGTCTGGCAGGATGGGTTTGTCTTGGAAATCCTCTAGGAGG
 CTCCTCCTCGCATGGCCTGCAGTCTGGCAGCAGCCCGAGTTGTTTCTCTCGCTGATCGAT
 TTTCTTCCAGGTAGATTTTTCTTTGCTTATGTTGAATTCATTGCCTTTTTCTCAT
 CACAGAAGTGATGTTGGAATCGTTTCTTTGTTGTCTGATTTATGGTTTTTTAAGTAT
 AAACAAAAGTTTTTTATTAGCATTCTGAAAAGAAAGTAAAATGTACAAGTTTAAATA
 AAAGGGGCTTCCCTTTAGAATAAATTCAGCATGTGCTTTCAAAAAAAAAAAAAAAAAAA
 AA

Restriction Sites: Please inquire

ACCN: NM_021102

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021102.2](#), [NP_066925.1](#)

RefSeq Size: 1571 bp

RefSeq ORF: 759 bp

Locus ID: 10653

UniProt ID: [O43291](#)

Cytogenetics: 19q13.2

Domains: KU

Protein Families: Transmembrane

Gene Summary: This gene encodes a transmembrane protein with two extracellular Kunitz domains that inhibits a variety of serine proteases. The protein inhibits HGF activator which prevents the formation of active hepatocyte growth factor. This gene is a putative tumor suppressor, and mutations in this gene result in congenital sodium diarrhea. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]
Transcript Variant: This variant (a) represents the longer transcript and encodes the longer isoform (a).