

Product datasheet for SC322325

ETFA (NM_000126) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ETFA (NM_000126) Human Untagged Clone
Tag: Tag Free
Symbol: ETFA
Synonyms: EMA; GA2; MADD
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC (PS100020)
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for SC322325
CGAGGTTGCGGCGGAAGCGGAGACCATGTTCCGAGCGCGGCTCCGGGGCAGCTCCGGCG
GGCGGCCTCATTGCTACGATTTTCAGAGTACCCTGGTAATAGCTGAGCATGCAAATGATTC
CCTAGCACCCATTACTTTAAATACCATTACTGCAGCCACACGCCTTGGAGGTGAAGTGTC
CTGCTTAGTAGCTGGAACCAAATGTGACAAGGTGGCACAAGATCTCTGTAAGTAGCAGG
CATAGCAAAGTCTGGTGGCTCAGCATGATGTGTACAAAGGCCTACTTCCAGAGGAACT
GACACCATTGATTTGGCAACTCAGAAGCAGTTCAATTACACACACATCTGTGCTGGAGC
ATCTGCCTTCGAAAGAACCTTTTGGCCAGAGTAGCAGCCAACTTGGAGTTGCCCGAT
TTCTGACATCATTGCAATCAAGTCACTGACACATTTGTGAGAACTATTTATGCAGGAAA
TGCTCTATGTACAGTGAAGTGTGATGAGAAAGTAAAGTGTCTGTCCGTGGAACATC
CTTTGATGCTGCAGCAACAAGTGGCGGTAGTGCCAGTTCAGAAAAGGCATCAAGTACTTC
ACCAGTGGAAATATCAGAGTGGCTTGACCAGAAATTAACAAAAAGTGATCGACCAGAGCT
AACAGGTGCCAAAGTGGTGGTATCTGGTGGTCCGAGGCTTGAAGAGTGGAGAGAACTTAA
GTTGTTATATGACTTGGCAGATCAACTACATGCTGCAGTTGGTGGTCTCCCGTGGTGTGT
TGATGCTGGCTTTGTTCCCAATGACATGCAAGTTGGACAGACGGGAAAAATAGTAGCACC
AGAACTTTATATTGCTGTTGGAATATCTGGAGCCATCCAACATTTAGCTGGGATGAAAGA
CAGCAAGACAATTGTGGCAATTAATAAAGACCCAGAAGCTCCAATTTTCCAAGTGGCAGA
TTATGGAATAGTTGCAGATTTATTTAAGGTAGTTCCCTGAAATGACTGAGATATTGAAGAA
AAAATGAATCAGGATCATGCCTTAAAAAGAAAACCTTTTGTAAAGTATTCCTGAAATC
ACAGATATTTGTTGGTATTATAACAATCATTGGAAAGCATGGAGAGCTACATTTTCATAAT
TTGAGGAAAAATTTCTAACAGATGCCAGAATGCTTGTATGGGATTGCTGTGTTTCCTT
TTAATTATTTGTTGGTCCAAACAATTATTGTTTGAACCTTTTTAATTCTGTACTAAAAATC
TATAATAAAGCTTTTCCACAGCTTTAAACTAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Please inquire
ACCN: NM_000126



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_000126.2</u> , <u>NP_000117.1</u>
RefSeq Size:	1402 bp
RefSeq ORF:	1002 bp
Locus ID:	2108
UniProt ID:	<u>P13804</u>
Cytogenetics:	15q24.2-q24.3
Domains:	ETF_alpha
Protein Families:	Druggable Genome
Gene Summary:	<p>ETF A participates in catalyzing the initial step of the mitochondrial fatty acid beta-oxidation. It shuttles electrons between primary flavoprotein dehydrogenases and the membrane-bound electron transfer flavoprotein ubiquinone oxidoreductase. Defects in electron-transfer-flavoprotein have been implicated in type II glutaricaciduria in which multiple acyl-CoA dehydrogenase deficiencies result in large excretion of glutaric, lactic, ethylmalonic, butyric, isobutyric, 2-methyl-butyrilic, and isovaleric acids. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).</p>