

Product datasheet for **SC322330**

MFSD3 (NM_138431) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MFSD3 (NM_138431) Human Untagged Clone
Tag:	Tag Free
Symbol:	MFSD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC322330

```

CGAGAGACGTGTCCAGACCCTAGCCTGTACGACGCTGACTCTGCCCGGTCCCAGAACCA
AAGCCATGCCGGGTGTGGCCTCTGACCCACGCGGAGGGGACCTCGCCTTGCGGGACCC
CACCTGGAACCCGACCTCCCAGCCTCGCAGCCGGCCTGAGCCGCCATGCGCGGGAAGTTG
CTGCCGCTGGCCGGCCTTACCTGGTGCAGGGCCTGCCCTACGGGCTCCAGTCCGGCCTC
CTGCCAGTGCTGCTGCGTGCCGGCGGCCTCTCGCTGACGCGCGTGGGGCTGGCCAAGTT
CTGTACGCTCCGTGGTGCTCAAGCTGGCTTGGGCCCCGCTGGTGGACGCGCAGGGCTCG
GCGAGGGCCTGGGTGACGCGCAGCAGCGGGCCTGGGCTGGTGTGGGCTGCTTGCC
GGGCTGCCCCCTCTGGAGCTGGCCAGGCCGGGCTGCCCGCCGCTGTGGCGGGGTGCTG
CTGTTGTTGAACCTGGGTGCCGCATGCAGGATGTGGCCCTGGACGCGCTGGCTGTGCAG
CTGCTGGAGCCGCCGAACCTGGGCGCGGCAATACCGTGCAGGTGGTCCGCTACAAGCTG
GGGCGCGCTAGCTGGGGCGCGCTGCTGGCGCTGCTGCCACCTTCTCGTGGCCGCAA
CTCTTTCTGCTCCTGGCTGCCACCTACTGGCTGGCCGCGGCCCTGGCCTGGGCTGCACCA
GCCCTGCGGCGGCTCCCACAGCAGCCCCCTTCCGAGCAGCGTCCCACACCGCGCACCTT
CTGCGGGACGTGCTAGCCGTGCCGGGACCGTGTGGACGCGCAGGCTTGTGCTCACCTAC
AAGCTGGGTGAGCAGGTGCCAGCAGCCTGTTTCTTCTCTCTGCTGGACCACGGCGTT
TCTGCTCCCGAGTTGGGACTGTGAATGGTGTGGGTGCTGTGGTCTGCTCCATCGCTGGC
TCCTCCCTGGGTGGGACCTTGTGGCCAAGCACTGGAACTGCTGCCTCTGTTGAGGTG
GTGCTGCGCTTCCGCTCGGGGGCCTAGCCTGTCAGACTGCCTTGGTCTTCCACCTGGAC
ACCTTGGGGGCCAGCATGGACGCTGGCACAATCTTGAGAGGGTCAGCCTTGTGACCTA
TGTCTGCAGCACTTCTTGGGAGGCTGGTCACCACAGTCACTTCACTGGGATGATGCGC
TGCAGCCAGCTGGCCCCAGGGCCCTGCAGGCCACACTACAGCCTTCTGGCCACGCTG
GAGCTGCTGGGAAGCTGCTGCTGGGCACTCTGGCCGAGGCTGGCTGATGGGTTGGGG
CCACATCCCTGCTTCTTGTCTCTGCTATCCTCTCTGCTTCCCGTTCTGTACCTGGAC
CTAGCACCCAGCACCTTCTTTGAGCTGAGTGGCTGGAGTGGTCAATAAAGCCACATGTG
CCTGTGAAAAAAAAAAAAAAAAAAAAA

```

Restriction Sites: Please inquire



[View online »](#)

ACCN:	NM_138431
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_138431.1, NP_612440.1</u>
RefSeq Size:	1468 bp
RefSeq ORF:	1239 bp
Locus ID:	113655
UniProt ID:	<u>Q96ES6</u>
Cytogenetics:	8q24.3
Protein Families:	Transmembrane