

## Product datasheet for **SC322408**

### PBLD (NM\_022129) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** PBLD (NM\_022129) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** PBLD  
**Synonyms:** HEL-S-306; MAWBP; MAWDBP  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC (PS100020)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**Fully Sequenced ORF:** >OriGene sequence for SC322408  
 GCCTCCTTTGTCCTGGATTTGGGAGTTGAGCACCTTCGTCGCCATTGGCTTTCTCCCC  
 AGCTCCAGCCTCTCATCTTGGGAATCTGCGTCAGAAGTCACTCGCAGTCCCGTCAGCC  
 CAGAAGAAGACGTAAGCAGGCTACCAGCAATTTTGAGAACTTGCAAAAACAGCTTGCAA  
 GGAAAATGAAGCTTCTATTTTCATAGCAGATGCATTCACAGCAAGAGCATTTCGTGGGA  
 ATCCTGCTGCTGTTGCTCCTAGAAAATGAATTGGATGAAGACATGCATCAGAAAATTG  
 CAAGGGAGATGAACCTCTGAAACTGCTTTTATCCGAAAACGCACCCGACAGACAACT  
 TTGCACAAAGTTCCTGCTTTGGACTGAGATGGTTTACACCAGCGAGTGAGGTCCCACTCT  
 GTGGCCATGCCACCCTGGCTTCTGCAGCTGTGCTGTTTACAAAAATAAAAAACATGAATA  
 GCACGCTCACGTTTGTCACTCTGAGTGGAGAATAAGGGCCAGACGAGCAGAGGACGGCA  
 TCGTCCTGGACTTGCCTCTTATCCAGCCACCCCAAGGACTTCCATGAAGTAGAGGACT  
 TGATAAAGACTGCCATAGGCAACACACTGGTCCAGGACATCTGTTATTCTCCAGATACCC  
 AAAAGCTCCTCGTCCGCCTCAGTGACGTTTACAACAGGTCGTTTCTGGAGAACCTGAAAG  
 TGAACACGGAGAACTCTGCTGCAAGTTGAAAACACAGGGAAGGTGAAAGGGCTTATTCTTA  
 CCCTTAAAGGAGAGCCTGGTGGGACAGCCCAAGCATTTGACTTTTACTCAAGATATTTTG  
 CACCGTGGGTTGGTGTGGCTGAAGACCCAGTGACAGGGTCTGCACACGCTGTTCTCAGCA  
 GCTACTGGTCCCAGCATCTGGGGAAGAAAGAAATGCATGCTTTTTCAGTGTCCCACCGAG  
 GAGGAGAGCTGGGAATTTCCCTTCGTCCAGACGGAAGGGTTGACATTAGAGGAGGTGCAG  
 CTGTTGTTTTAGAGGGCACACTGACAGCCTAGAGGTGGTTATGCTGTGACGCTGCTGTCT  
 CTAAACCACCAAGTATTTTCTGCTTAAAAAGAAATGTAAGGGGCTGCCTTTAGCAAATGTG  
 CGTAGTAGTCTACTTAATCCTCATGTATAAAAAAAAAAAAAA

**Restriction Sites:** Please inquire  
**ACCN:** NM\_022129



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_022129.3</a></u> , <u><a href="#">NP_071412.2</a></u>
<b>RefSeq Size:</b>	2587 bp
<b>RefSeq ORF:</b>	867 bp
<b>Locus ID:</b>	64081
<b>UniProt ID:</b>	<u><a href="#">P30039</a></u>
<b>Cytogenetics:</b>	10q21.3
<b>Domains:</b>	PhzC-PhzF