

Product datasheet for SC322416

DNAJC12 (NM_021800) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DNAJC12 (NM_021800) Human Untagged Clone
Tag:	Tag Free
Symbol:	DNAJC12
Synonyms:	HPANBH4; JDP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC322416 GGAACATGCTGCTCTAATTCAGATTTAAAGAGTTTCTTCTGTTAATTCGAAGCTCACTG TGCCTCTTGTTTCCGAGGGAAGAAGGACTGATTAAGTCATCTAAATGGATGCAATACTGA ATTACAGGTGAGAAGATACTGAAGATTACTACACATTACTGGGATGTGATGAACATATCTT CGGTTGAACAAATCCTGGCAGAATTTAAAGTCAGAGCTCTGGAATGTCACCCAGACAAGC ATCCTGAAAACCCAAAGCTGTGGAGACTTTTCAGAACTGCAGAAGGCAAAGGAGATTC TGACCAATGAAGAGAGTCGAGCCCGCTATGACCACTGGCGAAGGAGCCAGATGTCGATGC CATTCCAGCAGTGGGAAGCTTTGAATGACTCAGTGAAGACGTCAATGCACTGGGTTGTCA GAGGTAAAAAGACCTGATGCTGGAAGAATCTGACAAGACTCATACCAAGATGGAAA ATGAGGAATGTAATGAGCAAAGAGAAAGAAAGAAAGAGGAGCTGGCTTCAACCGCAGAGA AAACGGAGCAGAAAGAACCCAGCCCTAGAGAAGTCAGTCTCCCGCAAAATTCAGATT CTTGAGGTTTTGCAGATGTGAATGGTTGGCACCTTCGTTTCCGCTGGTCCAAGGATGCTC CCTCAGAACTCCTGAGGAAGTTCAGAACTATGAAATATGAAATATCTCTGCTTCAAAA ATGAGGAAGAGCAAGACTGTCCCTATGCTGCCAATGCAGTCTTTGTTTATGTCTTAA AAATGTCATGTTTATGTCATGCTGTGAATTGCTGAGTACTAATTGATTCTCCATCCTT GAATCAGTTCTCATAATGCTTTTTAAATAAGAAAAATTCAGAAGATGAATTTCTTCCAAT ATTTGAATAAATTAAGCTCTTAGATACAGAGTAGATTGTATTATGCTTTTTCTTATT AATACTACTTATAGAAATCCATTAAGCAATCTCTGTACAGTGATTTAAATATTTCA TTGACATACTGTGATCTCTACAATGAAATGTGAAATTAATTGTCTGAACCCCGTGGGGAG AAATAAATAATTTCCCAAGTTCAAAAAAAAAAAAAAAAAA
Restriction Sites:	Please inquire
ACCN:	NM_021800


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_021800.2 , NP_068572.1
RefSeq Size:	1218 bp
RefSeq ORF:	597 bp
Locus ID:	56521
UniProt ID:	Q9UKB3
Cytogenetics:	10q21.3
Domains:	DnaJ
Gene Summary:	This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (1) represents the longer transcript, and encodes the longer isoform (a).