

Product datasheet for **SC322438**

CCDC74B (NM_207310) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC74B (NM_207310) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC74B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC322438

```

GGCCACTGAGCCGGGTGCAGTGGCAGCGGAGAGTACCTGGCGATGGCGATATGAGCGG
TGCGGGGTGGCGCTGGGACGCGGCCCCAGCTCGCCGACCCGGGCTCTCGGCGCCG
GCGCCAGCGCCCCCTCTGTGGGCGTCCAGTCCTTGAGGCCGAGAGCCCGAGCTCAGGCA
GAGCGACCCGAGAAACGGAACCTGGACCTGGAGAAAAGCCTGCAGTTCTGCAGCAGCA
GCACTCGGAGATGCTGGCCAAGCTCCATGAGGAGATCGAGCATCTGAAGCGGAAAAACA
GGATCTCCATTACAAGCTCATAATGAATCAGACATCACAGAAGAAAGATGGCCCCCAGG
AAACCACCTTTCCAGGGCCTCTGCTCCCTTGGGCGCTCGCTGGGTCTGCATCAACGGAGT
GTGGGTAGAGCCGGGAGGACCCAGCCCTGCCAGGCTGAAGGAGGGCTCCTCACGGACACA
CAGGCCAGGAGGCAAGCGTGGGCGTCTTGCGGGCGGTAGCGCCGACACTGTGCGCTCTCC
TGCAGACAGCCTCTCCATGTCAAGCTTCCAGTCTGTCAAGTCCATCTCTAATTCAGGCAA
GGCCAGGCCCCAGCCCGGCTCCTTCAACAAGCAAGATTCAAAAGCTGACGTCTCCAGAA
GGCGGACCTGGAAGAGGAGCCCCCTACTTCAACAAGCAAGCTGGACAAAGTTCTGGGGT
ACAAGGGCAGGCCAGAAAGGAGAAAGCAGAGGCCTCTAATGCAGGAGCTGCCTGTATGGG
GAACAGCCAGCACCAGGGCAGGCAGATGGGGGCGGGGCACACCCCCAATGATCCTGCC
CCTTCCCTGCGAAAGCCACCACACTTAGGCAGTGCGAAGTGCTCATCCGCGAGCTGTG
GAATACCAACCTCCTGCAGACCAAGAGCTGCGGCACCTCAAGTCCCTCCTGGAAGGGAG
CCAGAGGCCCCAGGCAGCCCCGAGGAAGCTAGCTTTCCAGGGACCAAGAAGCCACGCA
TTTCCCCAAGGTCTCCACCAAGAGCCTCTCCAAGAAATGCCTGAGCCACCTGTGGCGGA
GGTGCCATCCTGCCGCACTGAAGCAGACCCCGAAGAACAACTTTGCCGAGAGGCAGAA
GAGGCTGCAGGCAATGCAGAAACGGCGCCTGCATCGCTCAGTGCTTTGAGCCACCCCAAT
CTGGTCAGTGCCAGGCCACCAACCTGCAGCTGGAGACTGGCTCTCTATAGCATTTCCTG
ATACTTCGCTACTTTTAGGCCTGGCTAAATTCCAAGACAGATAACACTCAAGATAGATA
AAGTACTTGATCTCCAACTGACAACTGTTTATTTTCTAGCTGTTATTTTGCTATTTGG
CATTTACATAAAAGCACACGATGAAGCAGGTATCGCCTTACCTGTTGAACTGAAATAA
AGCTTGTTTATTTCAAAAAAAAAAAAAAAAAA

```

Restriction Sites: Please inquire



[View online »](#)

ACCN:	NM_207310
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207310.1, NP_997193.1</u>
RefSeq Size:	1471 bp
RefSeq ORF:	1143 bp
Locus ID:	91409
UniProt ID:	<u>Q96LY2</u>
Cytogenetics:	2q21.1