

## Product datasheet for SC322444

## IKB beta (NFKBIB) (NM 002503) Human Untagged Clone

**Product data:** 

**Product Type: Expression Plasmids** 

**Product Name:** IKB beta (NFKBIB) (NM\_002503) Human Untagged Clone

Tag: Tag Free IKB beta Symbol:

IKBB; TRIP9 Synonyms: **Mammalian Cell** 

Selection:

Neomycin

Vector: pCMV6-AC (PS100020) E. coli Selection: Ampicillin (100 ug/mL)

**Fully Sequenced ORF:** >OriGene sequence for SC322444

> GCGACTGCGGGGGCCCCTGAGGCGGCGGGGCCATGGCTGGGGTCGCGTGCTTGGGAAA AGCTGCCGACGCAGATGAATGGTGCGACAGCGGCCTGGGCTCCCTGGGTCCGGACGCAGC GGCCCCGGAGGACCTGGGTTGGGCGCGGAGTTGGGCCCGGGGCTGTCGTGGGCTCCCCT CGTCTTCGGCTACGTCACTGAGGATGGGGACACGGCACTGCACTTGGCTGTGATTCATCA GCAGAATGACCTAGGCCAGACAGCCCTGCACCTGGCAGCCATCCTGGGGGAGACATCCAC GGCGCTGCACCTGGCCTGTGGGGGGCACACGCCTGTGCCCGTGCCCTGCTTCAGCC CCGCCCCGGCGCCCCAGGGAAGCCCCCGACACCTACCTCGCTCAGGGCCCTGACCGTAC AGAGGAGAGTGAGGACTGGAAGCTGCAGCTGGAGGCTGAAAACTACGAGGGCCACAC CCCACTCCACGTGGCCGTTATCCACAAAGATGTGGAGATGGTCCGGCTGCTCCGAGATGC TGGAGCTGACCTTGACAAACCGGAGCCCACGTGCGGCCGGAGCCCCCTTCATTTGGCAGT GGAGGCCCAGGCAGCCGATGTGCTGGAGCTTCTCCTGAGGGCAGGCGCGAACCCTGCTGC CCGCATGTACGGTGGCCGCACCCCACTCGGCAGTGCCATGCTCCGGCCCAACCCCATCCT CGCCCGCCTCCTCCGTGCACACGGAGCCCCTGAGCCCGAGGGCGAGGACGAGAAATCCGG CCCCTGCAGCAGCAGTAGCGACAGCGACAGCGGAGACGAGGGCGATGAATACGACGACAT TGTGGTTCACAGCAGCCGAGCCAAACCCGGCTGCCTCCCACCCCAGCCTCAAAACCTCT TCCTGACGACCCCCGCCCCGTGTGATTTGTTTCATTGTTAATATAATTTCCAGTTTAATA

**Restriction Sites:** Please inquire ACCN: NM 002503



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

## IKB beta (NFKBIB) (NM\_002503) Human Untagged Clone - SC322444

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a

point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative

RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning

into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 002503.3, NP 002494.2

RefSeq Size: 1198 bp RefSeq ORF: 1071 bp Locus ID: 4793

 UniProt ID:
 Q15653

 Cytogenetics:
 19q13.2

Domains: ANK

**Protein Families:** Stem cell - Pluripotency, Transcription Factors

**Protein Pathways:** Adipocytokine signaling pathway, B cell receptor signaling pathway, Chemokine signaling

pathway, Cytosolic DNA-sensing pathway, Neurotrophin signaling pathway, NOD-like receptor signaling pathway, RIG-I-like receptor signaling pathway, T cell receptor signaling pathway

**Gene Summary:** The protein encoded by this gene belongs to the NF-kappa-B inhibitor family, which inhibit

NF-kappa-B by complexing with, and trapping it in the cytoplasm. Phosphorylation of serine residues on these proteins by kinases marks them for destruction via the ubiquitination pathway, thereby allowing activation of the NF-kappa-B, which translocates to the nucleus to function as a transcription factor. Alternatively spliced transcript variants have been found for

this gene.[provided by RefSeq, Jul 2011]

Transcript Variant: This variant (1) represents the predominant transcript and encodes the

longest isoform (1).