

Product datasheet for SC322612

C11ORF77 (HARBI1) (NM_173811) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C11ORF77 (HARBI1) (NM_173811) Human Untagged Clone
Tag:	Tag Free
Symbol:	C11ORF77
Synonyms:	C11orf77
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC322612 <pre> GGGAAGTGTGGCTGGTAACCTGGCAGCCGCGAGAGAGAGAAGATTATAAATGGCAGAGC CATTTAACTGTGGTTAGCTGTTGGATTCTGATACTTTCTTAAAAATACGTTCTTGCACCA ACATCTTCATTGGGAACAGTTTCAGAAAAAGCAAAGAGGAGAGCCAAACATTCACATTTACC ATGGCTATACCAATAACAGTGCTTGACTGTGACCTCTTGCTATATGGCCGTGGTCACCGG ACATTGGACCGTTTTAAGCTGGATGATGTGACTGATGAATACTTGATGTCCATGTATGGG TTTCCACGGCAGTTTCATTATTACTTGGTGGAGCTCTTGGGGGCGAATCTTTCTAGGCCT ACTCAGCGATCCAGGGCTATTAGCCCAGAGACACAGGTCCTTGACGATTGGGTTTTTAT ACCTCAGGTTCTTCCAGACTCGGATGGGAGATGCCATTGGAATCAGTCAGGCGTCTATG AGTCGTTGTGTTGCCAATGTCACTGAAGCACTTGTGAAAGGGCCTCACAGTTCATTGCG TTTCCAGCTGATGAAGCCTCCATTCAGGCTCTGAAGGATGAATTCTATGGGTTGGCAGGG ATGCCAGGGGTGATGGGGGTGGTTGACTGTATCCATGTGGCCATCAAGGCACCAATGCT GAAGACCTCTCCTATGTGAACCGAAAAGGCCTGCATTCTTTAACTGCCTGATGGTGTGT GACATTAGAGGGACACTAATGACCGTGGAGACAACTGGCCCGGCAGCCTACAGGACTGT GCTGTGCTGCAGCAGTCTTCCCTCAGTAGTCAGTTTGAAGCGGTATGCACAAAGATAGC TGGCTTCTGGGTGACAGTTCTTCTTTCTTGAACCTGGCTCATGACCCCACTTCACATT CCTGAACTCCAGCAGAATATCGCTATAACATGGCCCATCTGCAACTCACAGTGTGATT GAGAAGACTTTCCGAACCTCTGCTCCCGATTCCGCTGCCTGGATGGATCCAAGGGGGCA CTGCAGTACTCACCAGAGAAATCCAGCCATATCATCTTGGCCTGTTGTGTCCTCCACAAC ATCTCCCTGGAGCATGGGATGGATGTTTGGTCCTCTCCAATGACAGGACCATGGAACAG CCCCCGGAAGAGGAGTATGAGCACATGGAGTCCCTGGACTTAGAGGCTGACCGTATTGCT CAGGAGCTAATGCTCACTCATTTTAGCTAATGTAGAAGGTGGAGAGGAGGGATACTTCCC AGGAGTTGTGACAGACTTTCCTCCTCATCACCTTTTACACAGTTCATCATCTAGCATGA CTGAGTATACAGATACTGTGCATAAACTGACATTTAATATGTGTGTTTTGGTAAGGTTGG GGCTATGCCAGAATATCTTGATTCAATTTGCATATGCATTAATTTAACTGAAACCAAGACA AA AA </pre>



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_173811
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173811.2</u> , <u>NP_776172.1</u>
RefSeq Size:	1966 bp
RefSeq ORF:	1050 bp
Locus ID:	283254
UniProt ID:	<u>Q96MB7</u>
Cytogenetics:	11p11.2
Gene Summary:	Transposase-derived protein that may have nuclease activity (Potential). Does not have transposase activity.[UniProtKB/Swiss-Prot Function]