

Product datasheet for **SC322672**

TAB1 (NM_006116) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAB1 (NM_006116) Human Untagged Clone
Tag:	Tag Free
Symbol:	TAB1
Synonyms:	3'-Tab1; MAP3K7IP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for SC322672
 CTCTAGCCGCCCACTCGCTGCTCTGCGGGGAGGCGGGCGCTCCCGCAGGGGTTCTCCA
 AGATGGCGGCGCAGAGGAGGAGCTTGCTGCAGAGTGAGCAGCAGCCAAGCTGGACAGATG
 ACCTGCCTCTCTGCCACCTCTCTGGGGTTGGCTCAGCCTCCAACCGCAGCTACTCTGCTG
 ATGGCAAGGGCACTGAGAGCCACCCGCCAGAGGACAGCTGGCTCAAGTTCAGGAGTGAGA
 ACAACTGCTTCTGTATGGGGTCTTCAACGGCTATGATGGCAACCGAGTGACCAACTTCG
 TGGCCCAGCGGCTGTCCGCAGAGCTCCTGCTGGGCCAGCTGAATGCCGAGCACGCCGAGG
 CCGATGTGCGGCGTGTGCTGCTGCAGGCCTTCGATGTGGTGGAGAGGAGCTTCTGGAGT
 CCATTGACGACGCTTGGCTGAGAAGGCAAGCCTCCAGTCGCAATTGCCAGAGGGAGTCC
 CTCAGCACCAGCTGCCTCCTCAGTATCAGAAGATCCTTGAGAGACTCAAGACGTTAGAGA
 GGGAAATTTCCGGAGGGGCCATGGCCGTTGTGGCGGTCTTCTCAACAACAAGCTCTACG
 TCGCCAATGTCCGTACAAACCGTCACTTTTATGCAAATCGACAGTGGATGGGTTGCAGG
 TGACACAGCTGAACGTGGACCACACCACAGAGAACGAGGATGAGCTCTTCCGCTTTTCG
 AGCTGGGCTTGGATGCTGGAAAGATCAAGCAGGTGGGGATCATCTGTGGGCAGGAGAGCA
 CCCGGCGGATCGGGGATTACAAGTTAAATATGGCTACACGGACATTGACCTTCTCAGCG
 CTGCCAAGTCAAACCAATCATCGCAGAGCCAGAAATCCATGGGGCACAGCCGCTGGATG
 GGGTGACGGGCTTCTTGGTGTGATGTGCGAGGGGTTGTACAAGGCCCTAGAGGCAGCCC
 ATGGGCCCTGGGCAGGCCAACAGGAGATTGCTGCGATGATTGACACTGAGTTTGCCAAGC
 AGACCTCCCTGGACGCACTGGCCAGGCCGTCGTGGACCGGGTGAAGCGCATCCACAGCG
 ACACCTTCGCCAGTGGTGGGAGCGTGCCAGTTCTGCCCGGCACGAGGACATGACCC
 TGCTAGTGAGGAACCTTGGCTACCCGCTGGGCGAAATGAGCCAGCCACACCCGAGCCAG
 CCCCAGCTGCAGGAGGACGAGTGTACCCTGTGTCTGTGCCATACTCCAGCGCCAGAGCA
 CCAGCAAGACCAGCGTGACCCTCTCCCTTGTGATGCCCTCCAGGGCCAGATGGTCAACG
 GGGCTCACAGTGTTCACCCCTGGACGAAGCCACCCCCACCTCACCAACCAAAGCCGA
 CCTTAACCCTGCAGTCCACCAACACGCACACGCAGAGCAGCAGCTCCAGCTCTGACGGAG
 GCCTCTTCCGCTCCCGGCCCGCCACTCGCTCCCGCTGGCGAGGACGGTCGTGTTGAGC
 CCTATGTGGACTTGTGAGTTTTACCGCTCTGGAGCGTGGACCATGGCGAGCAGAGCG



[View online »](#)

```
TGGTGACAGCACCGTAGGGCAGCCGGAGAATGCAGCCCAAGCAGGGCCTGGCATGGGGCA
GGACAGGGTCCAGCCTTTTCTAACATCTGCCTGTGCCACAACGGCCAGCAGGTGCCCA
TCCTCTGCCCACAGCAGACTCTGTCCCATGGCTCTCCGGCAGTAGAGTGTGTAGTGCA
GACTGGACCTGTGGTTTACCTTGTACCACCCGGGAAGCTGAAGGCCACTTCTCCCA
GATGGCCTCAGCCAGGACCATCGCCCTTCTCAGAGCAGAGGGCCAGGTATAGAAACCGC
AGTGGGCTGCAAGCCGCCCCGAGCCTCCCAGCAGCCTCCTACAGAGCAGGAAGAGCCGCT
CCTGTGAACCCCTGAGTGTTCAGGCCCAGCAGACCCTGCTGTCCCAAGCCACCCCTCCT
CCCACCATCACCTCCCTCACCTCGGGACAGTAGCCCTCCACTTCTCCAGCCTCTCAGCCC
TGTGCTCCTGCATCCAGAGTGGAACCCAGGCTGGTGTCCGCATCTGTCCCTGGGCCCCAC
CCCTGGACCTGCCTTGGTTGTGCATCTGTTGTAACGTTTCCAGGAGGACCAGGGCAGCAT
CTGGGGCTGGGATGGCCACAGAAGGGGCAGGCCAGGTGAAAGGAGCCAGGGGGAAGTG
GTCTAAGAGACCTGAACTGCCAGAGGATGGCGGCTGGGCTTCCCAGAGCCAGGCGTG
CGGGAGAGGTGAGGACTGGCCCCGGTGGGCTGAGGCAGGGGCCGCTGTCGTCAGGCTGA
GCCAGGCTGAGCTGGTGCCTGCCTTGCATTTTCTTCTGGTGCTGAAGACCATAGGCT
GGCAGGCAGCTGAGATGAACTGTCTTTACCACTGATGAGGGCCCTGCGCGGCTGAGGGT
AGCAAGCAGGGGTTGTGAGTCAGGCTGGGGACTTGTGTTGAAAGAAAGAGGAGTGAAAA
TGGTCCAGGAGGGAAGAGGTTCTTTGAGACACAGTACCCTGGGAGGCATAGGAGAGGG
TCGGGCCAGCCCAGCCCAGGGCCTGAGTTAGACTATTTCCACATGTTCTCTGCCTTCAG
TGGGGAGGGGTGCCACCAGGGCTGTCCGCCAGGATTGCCACTCCTGTTTCAGAGGAAGC
AGGCCGAGAGACTTGACCTTGGCCAAGCCACACAATCAGTGGGGCAGCCAGAGCTCAGA
CCTGAGCCATTTTGTGATATCCAGGACCCCGGATTCTCCAGCCCTCCCATCTCCC
AGTCTCCCTGCCCCCATGCCCCAGACCCGGCCACCAGGGACTAGCCGCTGTCGCACAGC
CTCTGGGGTGTGGTCTCTGCGAAGTCAAAGGCTGACAGCTGTGGCCTGGGAATCC
ATTTTCTGCGGCAGAGCAGGGCCTGGTGTGGAACAGGGAGCTGTGGGAAGCCACAGCA
GAAATGGAAGAAAAACAGGTCTCAGCCCAGGGTCTCGCTCACTCCCTCACTCCCACTT
TGAAGCCATCTCTGTTCTGCAGGTGAGAGGATTTAAAGTCAGTCACAAAGGCTTGGGAAC
AAAAGGAATTCTTCCAAGAATGCCTCTGTGGTGTGTTTGGTCTCTAGAAACAGGGTC
ACTTTTTTAATGTAGTAAAGAAGTAATAAATGGTGGCATTACACATTGTGATATTAATA
AAAAAAAAAAAAA
```

Restriction Sites:

Please inquire

ACCN:

NM_006116

OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_006116.2, NP_006107.1</u>
RefSeq Size:	3240 bp
RefSeq ORF:	1515 bp
Locus ID:	10454
UniProt ID:	<u>Q15750</u>
Cytogenetics:	22q13.1
Domains:	PP2C
Protein Families:	Druggable Genome
Protein Pathways:	MAPK signaling pathway, NOD-like receptor signaling pathway, Toll-like receptor signaling pathway
Gene Summary:	<p>The protein encoded by this gene was identified as a regulator of the MAP kinase kinase kinase MAP3K7/TAK1, which is known to mediate various intracellular signaling pathways, such as those induced by TGF beta, interleukin 1, and WNT-1. This protein interacts and thus activates TAK1 kinase. It has been shown that the C-terminal portion of this protein is sufficient for binding and activation of TAK1, while a portion of the N-terminus acts as a dominant-negative inhibitor of TGF beta, suggesting that this protein may function as a mediator between TGF beta receptors and TAK1. This protein can also interact with and activate the mitogen-activated protein kinase 14 (MAPK14/p38alpha), and thus represents an alternative activation pathway, in addition to the MAPKK pathways, which contributes to the biological responses of MAPK14 to various stimuli. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (alpha) encodes the longer isoform (alpha).</p>