


```

CCCGGCGGGCCCCCGCCCCCTCCCCCTCCCCCTTCCCCTTCCCCTTCCCCTCCCAGCG
CGCCCCGCGCCCCGCGGCCCTCGGCGAGCAGCTCGGCTCCCCCAGCGCTCCCCGGGC
CAAAGATATGGCAATGGTAGTTAGCAGCTGGCGAGATCCGCAGGACGACGTGGCCGGGG
CAACCCCGGGCGCCCCAACCCCGCAGCGCAGGCGGCCCGCGGGCGGGCGGGCGGCCGG
CGAGCAGCAGCAGCAGGCGGGCTCGGGCGCGCCGCACACGCCGCAGACCCGGGCCAGCC
CGGAGCGCCCCGCCACCCCGGCACGGCGGGGGACAAGGGCCAGGGCCCCCGCGTTCGGG
CCAGAGCCAGCAGCACATCGAGTGCCTGGTGTGCGGGGACAAGTCGAGCGGCAAGACTA
CGGCCAATTACCTGCGAGGGCTGCAAAAAGTTTCTTCAAGAGGAGCGTCCGCAGGAATT
AACTTACACATGCCGTGCCAACAGGAACTGTCCCATCGACCAGCACACCACGCAACCAAGT
CCAATACTGCCGCCTCAAGAAGTGCCTCAAAGTGGGCATGAGGCGGGAAGCGTTTACGG
AGGAAGAATGCCTCAACCCAGCCCAATCCAGGCCAGTACGCACTACCAACGGGGACCC
CCTCAACGGCCACTGCTACCTGTCCGGTACATCTCGTGTGCTGCGCGCCGAGCCCTA
CCCCACGTGCGCTACGGCAGCCAGTGCATGCAGCCCAACAACATTATGGGCATCGAGAA
CATCTGCGAGCTGGCCGCGCCCTGCTCTTCAGCGCCGTCGAGTGGGCCCGCAACATCCC
TTTCTCCGGATCTGCAGATCACCGACCAGGTGTCCCTGCTACGCCTCACCTGGAGCGA
GCTGTTTGTGCTCAACGCGGCCAGTGCTCTATGCCGCTGCACGTGGCGCCGTTGCTGGC
CGCCGCGGGCCTGCATGCCTCGCCCATGTCTGCCGACCGCGTCTGTCGCTTTCATGGACCA
CATCCGCATCTTCCAGGAGCAGGTGGAGAAGCTCAAGGCGCTACACGTGCACTCAGCCGA
GTACAGCTGCCTCAAAGCCATCGTGTGTTACGTGAGCAGCCTGTGGCCTGTGGATGC
GGCCACATCGAGAGCCTGCAGGAGAAGTCGAGTGGCAGTGGAGGAGTACGTGAGGAG
CCAGTACCCCAACCAGCCAGCCGTTTTGGCAAAGTGTGCTGCGACTGCCCTCGCTGCG
CACCGTCTCTCCTCCGTCATCGAGCAGCTCTTCTTCGTCGTTTGGTAGTAAAACCCC
CATCGAAACTCTCATCCCGGATATGTTACTGTCTGGGAGCAGCTTCAACTGGCCTTACAT
GTCCATCCAGTGTCTCTAGACCTTGGGCGCTTCCACCTGCCCGTCCCCCTAGAGACTC
AGAGGACCCACCTGGGCCAAGGACTCCAAAGCCGCGGGGACACCGGGAAGTGCAGCGGGC
CAGGCAGGCTGGGTGGGAGGAGGAGGGCCGAGACAGGAGCAGCCACCCAGCAGAAATA
CAATCCGAGCTACAAAGCATGGGAAAAAGAGACTTTTTAGGATCAGATCTGTGAGCAGC
TTGGCGAGGAAAAACAAAAA
    
```

Restriction Sites:

Please inquire

ACCN:

NM_005654

OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:

This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_005654.4 , NP_005645.1
RefSeq Size:	3210 bp
RefSeq ORF:	1272 bp
Locus ID:	7025
UniProt ID:	P10589
Cytogenetics:	5q15
Domains:	HOLI, zf-C4
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Gene Summary:	The protein encoded by this gene is a nuclear hormone receptor and transcriptional regulator. The encoded protein acts as a homodimer and binds to 5'-AGGTCA-3' repeats. Defects in this gene are a cause of Bosch-Boonstra optic atrophy syndrome (BBOAS). [provided by RefSeq, Apr 2014]