

## Product datasheet for **SC322767**

### FGFR1 (NM\_023105) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	FGFR1 (NM_023105) Human Untagged Clone
Tag:	Tag Free
Symbol:	FGFR1
Synonyms:	bFGF-R-1; BFGFR; CD331; CEK; ECCL; FGFBR; FGFR-1; FLG; FLT-2; FLT2; HBGFR; HH2; HRTFDS; KAL2; N-SAM; OGD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for SC322767

```

GAAAAGGAGGATCGAGCTCACTGTGGAGTATCCATGGAGATGTGGAGCCTTGTCCACCAAC
CTCTAACTGCAGAACTGGGATGTGGAGCTGGAAGTGCCTCCTCTCTGGGCTGTGCTGGT
CACAGCCACACTCTGCACCGTAGGCCGTCCCCGACCTTGCCCTGAACAAGATGCTCTCCC
CTCCTCGGAGGATGATGATGATGATGATGACTCCTCTTCAGAGGAGAAAGAAACAGATAA
CACCAAACCAAACCCCGTAGCTCCATATTGGACATCCCCAGAAAAGATGAAAAGAAATT
GCATGCAGTGCCGGCTGCCAAGACAGTGAAGTTCAAATGCCCTTCCAGTGGGACCCCAAA
CCCCACACTGCGCTGGTTGAAAAATGGCAAAGAATTCAAACCTGACCACAGAATTGGAGG
CTACAAGGTCCGTTATGCCACCTGGAGCATCATAATGGACTCTGTGGTCCCTCTGACAA
GGGCAACTACACCTGCATTGTGGAGAATGAGTACGGCAGCATCAACCACACATACCAGCT
GGATGTCGTGGAGCGGTCCCCTCACGGCCCATCCTGCAAGCAGGGTTGCCCGCCAAACA
AACAGTGGCCCTGGGTAGCAACGTGGAGTTCATGTGTAAGGTGTACAGTGACCCGCAGCC
GCACATCCAGTGGCTAAAGCACATCGAGGTGAATGGGAGCAAGATTGGCCAGACAACT
GCCTTATGTCCAGATCTTGAAGACTGCTGGAGTTAATACCACCGACAAAGAGATGGAGGT
GCTTCACTTAAGAAATGTCTCCTTTGAGGACGCAGGGGAGTATACGTGCTTGGCGGGTAA
CTCTATCGGACTCTCCATCACTCTGCATGGTTGACCGTTCTGGAAGCCCTGGAAGAGAG
GCCGGCAGTGATGACCTCGCCCCGTACCTGGAGATCATCATCTATTGCACAGGGGCCCTT
CCTCATCTCCTGCATGGTGGGTCGGTCATCGTCTACAAGATGAAGAGTGGTACCAAGAA
GAGTGACTTCCACAGCCAGATGGCTGTGCACAAGCTGGCCAAGAGCATCCCTCTGCGCAG
ACAGGTAACAGTGTCTGCTGACTCCAGTGCATCCATGAACTCTGGGGTCTTCTGGTTCC
GCCATCACGGCTCTCCTCCAGTGGGACTCCCATGCTAGCAGGGGTCTCTGAGTATGAGCT
TCCCGAAGACCCTCGCTGGGAGCTGCCTCGGGACAGACTGGTCTTAGGCAAACCCCTGGG
AGAGGGCTGCTTTGGGCAAGTGGTGTGGCAGAGGCTATCGGGCTGGACAAGGACAAACC
CAACCGTGTGACCAAAGTGGCTGTGAAGATGTTGAAGTCGGACGCAACAGAGAAAGACTT
GTCAGACCTGATCTCAGAAATGGAGATGATGAAGATGATCGGGAAGCATAAGAATATCAT
CAACCTGCTGGGGCCTGCACGCAGGATGGTCCCTTGATGTCATCGTGGAGTATGCCTC

```



[View online »](#)

```

CAAGGGCAACCTGCGGGAGTACCTGCAGGCCCGGAGGCCCCAGGGCTGGAATACTGCTA
CAACCCAGCCACAACCCAGAGGAGCAGCTCTCCTCCAAGGACCTGGTGTCTGCGCCTA
CCAGGTGGCCCGAGGCATGGAGTATCTGGCCTCCAAGAAGTGCATACACCGAGACCTGGC
AGCCAGGAATGTCTGGTGACAGAGGACAATGTGATGAAGATAGCAGACTTTGGCCTCGC
ACGGGACATTACCACATCGACTACTATAAAAAGACAACCAACGGCCGACTGCCTGTGAA
GTGGATGGCACCCGAGGCATTATTTGACCGGATCTACACCCACAGAGTGTGTGGTTC
TTTCGGGGTGTCTCTGTGGGAGATCTTCACTCTGGGCGGCTCCCCATACCCCGGTGCC
TGTGGAGGAACCTTTCAAGCTGCTGAAGGAGGGTCACCGCATGGACAAGCCAGTAACTG
CACCAACGAGCTGTACATGATGATGCGGGACTGCTGGCATGCAAGTGCCTCACAGAGACC
CACCTTCAAGCAGCTGGTGAAGACCTGGACCGCATCGTGGCCTTGACCTCCAACAGGA
GTACCTGGACCTGTCCATGCCCCGGACCACTACTCCCCAGCTTTCCGACACCCGGAG
CTCTACGTGCTCCTCAGGGGAGGATTCCGTCTTCTCTCATGAGCCGCTGCCCGAGGAGCC
CTGCTGCCCCGACACCCAGCCAGCTTGCCAATGGCGGACTCAAACGCCGCTGACTGCC
ACCCACACGCCCTCCCCAGACTCCACCGTCAGCTGTAACCCTCACCCACAGCCCCGCTG
GGCCACACCTGTCCGTCCCTGTCCCCTTTCTGCTGGCAGGAGCCGGTGCCTACCAG
GGCCTTCTGTGTGGCCTGCCTTACCCCACTCAGCTCACCTCTCCCTCCACCTCCTCT
CCACCTGTGGTGTGAGAGGTGCAAAGAGGCAGATCTTTGCTGCCAGCCACTTCATCCCCTC
CCAGATGTTGGACCAACACCCCTCCCTGCCACCAAGGCACTGCCTGGAGGGCAGGGAGTGG
GAGCCAATGAACAGGCATGCAAGTGAGAGCTTCTGAGCTTTCTCCTGTCCGGTTTGGTCT
GTTTTGCCTTACCCATAAGCCCTCGCACTCTGGTGGCAGGTGCCTTGTCTCAGGGCT
ACAGCAGTAGGGAGGTGAGTGTTCGTGCCTCGATTGAAGGTGACCTCTGCCAGATAG
GTGGTGCCAGTGGCTTATTAATCCGATACTAGTTTGTCTTGACCAATGCCTGGTA
CCAGAGGATGGTGAAGCCGAGGTCAGTTGGGGCAGTGTGTGGCCCTGGGGCCAGCC
CCAAACTGGGGGCTCTGTATATAGCTATGAAGAAAACACAAAGTGTATAAATCTGAGTAT
ATATTTACATGTCTTTTTAAAAGGGTCGTTACCAAGAGATTTACCCATCGGGTAAGATGCT
CCTGGTGGCTGGGAGGCATCAGTTGCTATATATTAACCAAAAAAGAAAAAAGGAAA
ATGTTTTTAAAAGGTCATATATTTTTGCTACTTTTGTGTTTTATTTTTTAAATTAT
GTTCTAAACCTATTTTCAGTTTAGGTCCCTCAATAAAAATTGCTGCTGCTCAAAAAAAA
AAAAAAAAAAAAAAAA
    
```

**Restriction Sites:**

Please inquire

**ACCN:**

NM\_023105

**Insert Size:**

3200 bp

**OTI Disclaimer:**

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_023105.2</a> , <a href="#">NP_075593.1</a>
<b>RefSeq Size:</b>	5650 bp
<b>RefSeq ORF:</b>	2202 bp
<b>Locus ID:</b>	2260
<b>UniProt ID:</b>	<a href="#">P11362</a>
<b>Cytogenetics:</b>	8p11.23
<b>Domains:</b>	pkinase, TyrKc, S_TKc, ig, IGc2, IG
<b>Protein Families:</b>	Druggable Genome, Protein Kinase, Transmembrane
<b>Protein Pathways:</b>	Adherens junction, MAPK signaling pathway, Melanoma, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton

**Gene Summary:**

The protein encoded by this gene is a member of the fibroblast growth factor receptor (FGFR) family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds both acidic and basic fibroblast growth factors and is involved in limb induction. Mutations in this gene have been associated with Pfeiffer syndrome, Jackson-Weiss syndrome, Antley-Bixler syndrome, osteoglophonic dysplasia, and autosomal dominant Kallmann syndrome 2. Chromosomal aberrations involving this gene are associated with stem cell myeloproliferative disorder and stem cell leukemia lymphoma syndrome. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (3) lacks an alternate in-frame exon, compared to variant 1. The resulting isoform (3), also known as isoform Beta A1, II, H2, and the 2-Ig Domain+2 AA insert form, lacks the first Ig domain, compared to isoform 1. Both variants 3 and 13 encode the same isoform.