

Product datasheet for **SC322791**

CAPON (NOS1AP) (NM_001126060) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CAPON (NOS1AP) (NM_001126060) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAPON
Synonyms:	6330408P19Rik; CAPON; NPHS22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001128077, the custom clone sequence may differ by one or more nucleotides

[ATGAGCTACTATGGCACTACTACGGAGGACTGGGCTATGGCTATGACTGTAAATATAGTTATACCTCTG
GCTTTGGTGCCTTTAGAATCCTGGACTGTGGCTACAGATGTGGCTGTGGTGGGGTATGGATTTGA](#)

Restriction Sites:	Please inquire
ACCN:	NM_001126060
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001126060.1](#), [NP_001119532.2](#)

RefSeq Size: 5559 bp

RefSeq ORF: 636 bp

Locus ID: 9722

UniProt ID: [O75052](#)

Cytogenetics: 1q23.3

Gene Summary: This gene encodes a cytosolic protein that binds to the signaling molecule, neuronal nitric oxide synthase (nNOS). This protein has a C-terminal PDZ-binding domain that mediates interactions with nNOS and an N-terminal phosphotyrosine binding (PTB) domain that binds to the small monomeric G protein, Dexas1. Studies of the related mouse and rat proteins have shown that this protein functions as an adapter protein linking nNOS to specific targets, such as Dexas1 and the synapsins. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2009]
Transcript Variant: This variant (2) consists of the two terminal exons of variant 1 and results in a protein (isoform 2) that is significantly shorter than isoform 1(see PMID: 16146415).