

## Product datasheet for **SC322802**

### ZNF673 (KRBOX4) (NM\_001129900) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF673 (KRBOX4) (NM_001129900) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF673
Synonyms:	ZNF673
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	<p>&gt;NCBI ORF sequence for NM_001129900, the custom clone sequence may differ by one or more nucleotides</p> <p>ATGGCCATGTCCAGGAATCATTGACCTTCAAGGACGTGTTTGTGGACTTCACCTGGAG  GAGTGGCAGCAACTGGACTCTGCCAGAAGAACCTCTACAGGGATGTCATGCTTGAGAAC  TACAGCCACCTGGTGTCCGTGGGTATCTAGTTGCGAAGCCAGATGTGATCTTCAGGTTG  GGACCAGGTGAAGAGTCTGGATGGCAGATGGGGGACCCCGGTACGGACCTGTGCAGAT  GTATCCATTTTGCCTCATGTATTTGAAGTGCTGTTAT</p>
Restriction Sites:	Please inquire
ACCN:	NM_001129900
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>



[View online »](#)

**RefSeq:** NM\_001129900.1, NP\_001123372.1

**RefSeq Size:** 2487 bp

**RefSeq ORF:** 282 bp

**Locus ID:** 55634

**UniProt ID:** Q5JUW0

**Cytogenetics:** Xp11.3

**Protein Families:** Transcription Factors

**Gene Summary:** This encodes a zinc finger protein with an N-terminal KRAB (Kruppel-associated) domain found in transcriptional repressors. This gene is located in a region of the X chromosome thought to be involved in nonsyndromic X-linked cognitive disability. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]

Transcript Variant: This variant (4) uses a different splice site and contains a different segment in the 3' coding region that shifts the reading frame, compared to variant 1. The resulting protein (isoform 4) has a distinct C-terminus, compared to isoform 1. Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.