

Product datasheet for **SC322880**

HP1 alpha (CBX5) (NM_001127322) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HP1 alpha (CBX5) (NM_001127322) Human Untagged Clone
Tag:	Tag Free
Symbol:	CBX5
Synonyms:	HEL25; HP1; HP1A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC322880 representing NM_001127322. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGGGAAAGAAAACCAAGCGGACAGCTGACAGTTCTTCTTCAGAGGATGAGGAGGAGTATGTTGTGGAG
AAGGTGCTAGACAGGCGCGTGGTTAAGGGACAAGTGGAAATCTACTGAAGTGAAAGGCTTTTCTGAG
GAGCACAATACTTGGAACCTGAGAAAACTTGGATTGCCCTGAGCTAATTTCTGAATTTATGAAAAAG
TATAAGAAGATGAAGGAGGTGAAAATAATAAACCCAGGGAGAAGTCAGAAAGTAACAAGAGGAAATCC
AATTTCTCAAACAGTGCCGATGACATCAAATCTAAAAAAGAGAGAGCAGCAATGATATCGCTCGG
GGCTTTGAGAGAGACTGGAACCAGAAAAGATCATTGGGGCAACAGATTCCTGTGGTGATTTAATGTTC
CTAATGAAATGGAAAGACACAGATGAAGCTGACCTGGTTCTTGCAAAGAAGCTAATGTGAAATGTCCA
CAAATTGTGATAGCATTTTATGAAGAGAGACTGACATGGCATGCATATCCTGAGGATGCGAAAAACAA
GAGAAAGAAACAGCAAAGAGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
Plasmid Map:	<input type="checkbox"/>
ACCN:	NM_001127322
Insert Size:	576 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001127322.1
RefSeq Size:	11723 bp
RefSeq ORF:	576 bp
Locus ID:	23468
UniProt ID:	P45973
Cytogenetics:	12q13.13
MW:	22.2 kDa
Gene Summary:	<p>This gene encodes a highly conserved nonhistone protein, which is a member of the heterochromatin protein family. The protein is enriched in the heterochromatin and associated with centromeres. The protein has a single N-terminal chromodomain which can bind to histone proteins via methylated lysine residues, and a C-terminal chromo shadow-domain (CSD) which is responsible for the homodimerization and interaction with a number of chromatin-associated nonhistone proteins. The encoded product is involved in the formation of functional kinetochore through interaction with essential kinetochore proteins. The gene has a pseudogene located on chromosome 3. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2 and 3 encode the same protein. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence and transcripts to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>