

## Product datasheet for SC322914

### SYCE1L (NM\_001129979) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SYCE1L (NM_001129979) Human Untagged Clone
Tag:	Tag Free
Symbol:	SYCE1L
Synonyms:	MRP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC322914 representing NM_001129979. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGGATCGCC
ATGGCGGGGAAGCTGAAACCTCTGAATGTGGAGGCGCCAGAAGCTACTGAGGAGGCTGAAGGGCAAGCC
AAGTCTTTGAAGACTGAAGACTTGTGGCAATGGTGATAAAGCTGCAGAAAGAGGGAAAGCCTGGAGCCA
CAGATAGAGGACCTGATTAGCCGATTAATGATCTTCAGCAAGCAAGAAGAAATCCAGTGAGGAACTG
AGAGAGACCCACAGTCTCTGGGAGGCCCTGCATAGGGAATTAGACTCCTGAATGGAGAGAAAGTGCAC
CTAGAGGAGGTCTTGGCAAAAAGCAAGAGGCACTGAGGATCCTCCAGATGCACTGCCAAGAGAAGGAA
AGCGAGGCTCAGAGGTTGGATGTCAGAGGACAGCTGGAGGATCTGATGGGCCAGCACAAAGACCTCTGG
GAATCCACATGCTGGAGCAGCGACTGGCCCGGAGATCCGTGCCCTGGAGAGAAGCAAGGAGCAGCTG
CTCTCGGAGAGGAGGCTGGTGCGCGCAAGCTGCGGGAGGTGGAGCGCGGCTGCACTCGCCGCTGAG
GTTCGAGGGCGCCATGGCGGTGAATGACGGGCTGAAGGCGGAGCTGGAGATATTCGGGGAGCAGGTCGGG
AGCGCCCCGAGGTCGGGGCCGGCGAGGGAGAGCCGGACCTGAGCTCCCCCGCGCTCGCGACGAGGAG
GATCCCGAGCCCGGTGGCTGCCCTGACGCCCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001129979
Insert Size:	729 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001129979.1</a>
<b>RefSeq Size:</b>	877 bp
<b>RefSeq ORF:</b>	729 bp
<b>Locus ID:</b>	100130958
<b>UniProt ID:</b>	<a href="#">A8MT33</a>
<b>Cytogenetics:</b>	16q23.1
<b>MW:</b>	27.4 kDa
<b>Gene Summary:</b>	May be involved in meiosis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.