

Product datasheet for **SC322916**

ING4 (NM_001127583) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ING4 (NM_001127583) Human Untagged Clone
Tag:	Tag Free
Symbol:	ING4
Synonyms:	my036; p29ING4
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001127583, the custom clone sequence may differ by one or more nucleotides ATGGCTGCGGGGATGATTTGGAACATTATCTGGACAGTATTGAAAACCTTCCCTTTGAA TTACAGAGAAACTTTCAGCTCATGAGGGACCTAGACCAAAGAACAGAGGACCTGAAGGCT GAAATTGACAAGTTGGCCACTGAGTATATGAGTAGTGCCCGCAGCCTGAGCTCCGAGGAA AAATTGGCCCTTCTCAACAGATCCAGGAAGCCTATGGCAAGTGAAGGAATTTGGTGAC GACAAGGTGCAGCTTGCCATGCAGACCTATGAGATGGTGGACAAACACATTCGGCGGCTG GACACAGACCTGGCCGTTTTGAGGCTGATCTCAAGGAGAAACAGATTGAGTCAAGTGAC TATGACAGCTCTCCAGCAAAGAAGGCCGACTCAAAGGAGAAGAAAGCTGCTCGTGCT CGTTCCAAAGGGAAAACTCGGATGAAGAAGCCCCAAGACTGCCCAGAAGAAGTTAAAG CTCGTGCCGACAAGTCCCTGAGTATGGGATGCCCTCAGTGACCTTTGGCAGTGTCACCCC TCTGATGTGTTGGATATGCCTGTGGATCCCAACGAACCCACCTATTGCCTTTGTCACCAG GTCTCCTATGGAGAGATGATTGGCTGTGACAACCCTGATTGTTCCATTGAGTGGTTCCAT TTTGCCTGTGTGGGGCTGACAACCAAGCCTCGGGGAAATGGTTTTGCCACGCTGCTCC CAAGAACGGAAGAAGAAA
Restriction Sites:	Please inquire
ACCN:	NM_001127583
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001127583.1, NP_001121055.1</u>
RefSeq Size:	1452 bp
RefSeq ORF:	741 bp
Locus ID:	51147
UniProt ID:	<u>Q9UNL4</u>
Cytogenetics:	12p13.31
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	<p>This gene encodes a tumor suppressor protein that contains a PHD-finger, which is a common motif in proteins involved in chromatin remodeling. This protein can bind TP53 and EP300/p300, a component of the histone acetyl transferase complex, suggesting its involvement in the TP53-dependent regulatory pathway. Multiple alternatively spliced transcript variants have been observed that encode distinct proteins. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3) uses a different splice site in the coding region, compared to variant 9. The resulting protein (isoform 3) is shorter by 3 aa compared to isoform 9. Other names for this variant are ING4_v3 and '9 bp skip variant'.</p>