

Product datasheet for SC322951

SYNPR (NM_001130003) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	SYNPR (NM_001130003) Human Untagged Clone
Tag:	Tag Free
Symbol:	SYNPR
Synonyms:	SPO
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC322951 representing NM_001130003. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGACCTGTGAGTCAGCTGGCCTCTGCGGGCACCTCCGGGTGCTGAAGGAGCCCTTGCCTTCTG
CGAGCCCTGGAATTGCTTTTGAATCTTTGCATTTGCAACATGCGGTGGCTATTCTGGAGGCTGCGG
CTGAGTGTGGACTGCGTCAACAAGACAGAAAGTAACTCAGCATCGACATAGCGTTTGCCTACCCATTC
AGGTTGCACCAGGTGACGTTTGAGGTGCCACCTGCGAGGAAAGGAACGGCAGAAGCTGGCATTGATT
GGTGACTCCTCGTCTTCAGCAGAGTTCTTCGTCACTGTTGCTGTCTTCGCCTTCTCTACTCTTTGGCT
GCCACTGTCGTTTACATTTTCTCCAGAACAAATACCGGAAAACAACCGGGGCCACTCATTGACTTC
ATTGTCACTGTAGTCTTTTCGTTCTTGTGGTGGTGGTTCATCAGCTTGGGCAAAGGACTGTCTGAC
GTCAAAGTTGCAACGGATCCCAAGGAAGTATTGCTACTAATGTCAGCTTGCAAACAGCCATCCAACAAA
TGCATGGCTATCCACAGCCCTGTTATGTCAAGCTTAAACACTTCTGTGGTCTTTGGATTCTTGAACCTT
ATTCTCTGGGCTGGAAACATATGGTTTGTTTTCAAGGAGACCGGCTGGCATTCTTCGGGACAGAGATAT
CTTTCAGATCCAATGGAGAAGCACTCCAGCAGCTATAATCAAGGTGGTTACAACCAAGACAGCTATGGA
TCAAGCAGTGGGTACAGTCAGCAGGCGAGTTTGGGGCCAACCTCAGATGAGTTTGGCCAACAGCCTACT
GGCCCCACTTCTTTACCAATCAGATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001130003

Insert Size: 858 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130003.1](#)

RefSeq Size: 2594 bp

RefSeq ORF: 858 bp

Locus ID:	132204
UniProt ID:	<u>Q8TBG9</u>
Cytogenetics:	3p14.2
Protein Families:	Secreted Protein, Transmembrane
MW:	31.3 kDa
Gene Summary:	<p>Intrinsic membrane protein of small synaptic vesicles. Probable vesicular channel protein (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>