

## Product datasheet for **SC322956**

### Vgl4 (VGLL4) (NM\_001128219) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Vgl4 (VGLL4) (NM_001128219) Human Untagged Clone
Tag:	Tag Free
Symbol:	Vgl4
Synonyms:	VGL-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC322956 representing NM_001128219. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTATTTATGAAGATGGACCTGTTGAACTATCAGTACTTGGACAAGATGAACAACAATATCGGCATT
CTGTGCTACGAAGGCGAAGCTGCTCTCAGGGGAGAAGCCAGAATACAGACCCTGCCGGTGGCCTCTGCC
CTCAGCAGTCACCGCACCGGCCCTCCCCAATCAGCCCCAGCAAGAGGAAGTTCAGCATGGAGCCAGGT
GACGAGGACCTAGACTGTGACAACGACCACGTCTCCAAAATGAGTCGCATCTTCAACCCCATCTGAAC
AAGACTGCCAATGGAGACTGCCGACAGACCCCGGGAGCGGAGCCGAGCCCATCGAGCGCGCTGTG
GCCCCACCATGAGCCTGCACGGCAGCCACCTGTACACCTCCCTCCCGAGCCTTGGCCTGGAGCAGCCC
CTCGACTGACCAAGAACAGCCTGGACGCCAGCAGGCCAGCCGGCCTCTCGCCCACACTGACCCCGGGG
GAGCGGCAGCAGAACCGGCCCTCCGTGATCACCTGTGCCTCGGCTGGCGCCCGCAACTGCAACCTCTCG
CACTGCCCCATCGCGCACAGCGGCTGTGCCGCGCCCGGGCCTGCCAGCTACCGGAGGCCACCGAGCGCT
GCCACCACCTGTGACCCCGTGGTGGAGGAGCATTTCCGCAGGAGCCTGGGCAAGAATTACAAGGAGCCC
GAGCCGGCACCCAACCTCCGTGTCCATCACGGGCTCCGTGGACGACCACTTGGCAAGCTCTGGGTGAC
ACGTGGCTCCAGATCAAAGCGGCCAAGGACGGAGCATCCAGCAGCCCTGAGTCCGCTCTCGCAGGGGC
CAGCCCGCCAGCCCTCTGCCACATGGTCAGCCACAGTCACTCCCCCTGTGGTCTCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001128219
Insert Size:	891 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001128219.2</a></u>
<b>RefSeq Size:</b>	3845 bp
<b>RefSeq ORF:</b>	891 bp
<b>Locus ID:</b>	9686
<b>UniProt ID:</b>	<u><a href="#">Q14135</a></u>
<b>Cytogenetics:</b>	3p25.3-p25.2
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	31.9 kDa
<b>Gene Summary:</b>	May act as a specific coactivator for the mammalian TEFs.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (a).