

Product datasheet for **SC322959**

FAM71F2 (NM_001128926) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM71F2 (NM_001128926) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM71F2
Synonyms:	FAM137B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC322959 representing NM_001128926. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGAGTAAATCAGGGGCTCCACCTGAGGTCAGGGAACAGGCCCTGGAGTGGAACTGGGGTGGAA
AATGGCCTTCTTTGTCAACTGATTCATTCTCCAGAATTCAATTTGTTCTCCAACCTCGGTGGTGTGAA
AGCAACTTTATCCAGGTCACTAAGCCCGGGAACGGAGAGATGTCTGTGAAGGGTCTGCCACCGTGATC
CTCGGGGTCACTCCTCGGTGCCCTCCCTGCCACTCCCAATGTCCTCCTGATGGCCAATGTTACCTGG
CCCCAGGGTCCATTTACCACCTGGAGCACACCTGGTGACGCCCCAGTCATCAACCTCAGCAGGCTTCTC
CCCCTGAAGTACGTGGAGCTACGAATCTACGACCGGCTCCAGCGCATCCTGAGGGTTAGGACAGTGACC
GAAAAGATCTACTATCTGAAGCTCCACGAAAAACCCAGAGATTGTGTTTCAATTCTGGGTCCGCTTG
GTGAAAATTCTGCAGAAAGGCCTGTCCATCACCACCAAGACCCGAGAATCAAATTCACCTCACTGCCTG
GTGCCCAAGATGCCCACCAACTCCACAGAAACAACCTGAAAACAGCCTCCTGTATCCCCCAGCCC
AGCGAGCCCTCGTGCTGCTGGCGGCTGAGCAGACAGTGGCAGTTTCTCACAGCTCTCAGGAAAGCCC
CAGCTCACAGCAGACAGGAACAATGACTGTCATTGAAATAGACAACTGCAGCAGCTACAAGATACCC
TCTCCAGTGGCATCTCAATCAATTTGAATATACCCATGAGAGCTGCCTTGAGCCACAGCCTCTGGGAA
CAAGAGGACTGGAATGAGCACCTTCTACAAGTTTATATAGCCAGTTACCTAGGAGAGCACTTCTTGGGA
GCTTGA
ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001128926
Insert Size:	903 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001128926.3</u>
RefSeq Size:	5748 bp
RefSeq ORF:	903 bp
Locus ID:	346653
UniProt ID:	<u>Q6NXP2</u>
Cytogenetics:	7q32.1
MW:	33.5 kDa