

Product datasheet for **SC323019**

PTGER3 (NM_001126044) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: PTGER3 (NM_001126044) Human Untagged Clone
Tag: Tag Free
Symbol: PTGER3
Synonyms: EP3; EP3-I; EP3-II; EP3-III; EP3-IV; EP3-VI; EP3e; Inc003875; PGE2-R
Vector: pCMV6 series

Fully Sequenced ORF: >NCBI ORF sequence for NM_001126044, the custom clone sequence may differ by one or more nucleotides

```
ATGAAGGAGACCCGGGGCTACGGAGGGGATGCCCCCTTCTGCACCCGCTCAACCACTCC
TACACAGGCATGTGGGCGCCGAGCGTTCGCCGAGGCGGGGCAACCTCACGCGCCCT
CCAGGGTCTGGCGAGGATTGCGGATCGGTGTCCGTGGCCCTCCCGATCACCATGTGCTC
ACTGGTTTCGTGGGCAACGCACTGGCCATGCTGCTCGTGTGCGCAGCTACCGGCGCCGG
GAGAGCAAGCGCAAGAAGTCCTTCTGCTGTGCATCGGCTGGCTGGCGCTCACCGACCTG
GTCGGGCAGCTTCTCACCAACCCGGTCTGTCATCGTGTACCTGTCCAAGCAGCGTTGG
GAGCACATCGACCCGTCGGGGCGGCTCTGCACCTTTTTCGGGCTGACCATGACTGTTTTT
GGGCTCTCCTCGTTGTTATCGCCAGCGCCATGGCCGTCGAGCGGGCGCTGGCCATCAGG
GCGCCGCACTGGTATGCGAGCCACATGAAGACGCGTGCCACCCGCGCTGTGCTGCTCGGC
GTGTGGCTGGCCGTGCTCGCCTTCGCCCTGCTGCCGGTGTGGGCGTGGGCCAGTACACC
GTCCAGTGGCCCGGACGTGGTGTTCATCAGCACCGGGCGAGGGGGCAACGGGACTAGC
TCTTCGCATAACTGGGGCAACCTTTTCTCGCCTCTGCCTTTGCCTTCTCGGGCTCTTG
GCGCTGACAGTCACTTTTCTGCAACCTGGCCACCATTAAGGCCCTGGTGTCCCCTGC
CGGGCAAGGCCACGGCATCTCAGTCCAGTGCCAGTGGGGCCGATCACGACCGAGACG
GCCATTCAGCTTATGGGGATCATGTGCGTGTGCGTCTGCTGGTCTCCGCTCCTGATA
ATGATGTTGAAAATGATCTTCAATCAGACATCAGTTGAGCACTGCAAGACACACACGGAG
AAGCAGAAAGAATGCAACTTCTTCTTAATAGCTGTTGCGCTGGCTTCACTGAACCAGATC
TTGGATCCTTGGGTTTACCTGCTGTTAAGAAAGATCCTTCTTCGAAAGTTTTGCCAGATC
AGGTACCACACAAACAATATGCATCCAGCTCCACCTCCTTACCCTGCCAGTGTTCCTCA
ACCTTGATGTGGAGCGACCATTGGAAAGA
```

Restriction Sites: Please inquire

ACCN: NM_001126044

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001126044.1 , NP_001119516.1
RefSeq Size:	2055 bp
RefSeq ORF:	1173 bp
Locus ID:	5733
UniProt ID:	P43115
Cytogenetics:	1p31.1
Protein Families:	Druggable Genome, GPCR, Transcription Factors, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction
Gene Summary:	<p>The protein encoded by this gene is a member of the G-protein coupled receptor family. This protein is one of four receptors identified for prostaglandin E2 (PGE2). This receptor may have many biological functions, which involve digestion, nervous system, kidney reabsorption, and uterine contraction activities. Studies of the mouse counterpart suggest that this receptor may also mediate adrenocorticotrophic hormone response as well as fever generation in response to exogenous and endogenous stimuli. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009]</p> <p>Transcript Variant: This variant (11) has multiple differences compared to variant 1. The resulting protein (isoform 4) has a distinct and shorter C-terminus, as compared to isoform 1. Transcript variants 4, 9 and 11 encode the same protein. Another name for this transcript is EP3 subtype 1c.</p>