

Product datasheet for SC323030

SPRED2 (NM_001128210) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPRED2 (NM_001128210) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPRED2
Synonyms:	Spred-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC323030 representing NM_001128210. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCATCGCCTGGCAGTGACAGCTATATTGTGCGTGTCAAGGCTGTGGTTATGACCAGAGATGACTCC
AGCGGGGGATGGTTCCACAGGAAGGAGCGGGATCAGTCGCGTCGGGGTCTGTAAGGTCATGCACCCC
GAAGGCAATGGACGAAGCGGCTTCTCATCCATGGTGAACGACAGAAAGACAACTGGTGGTATTGGAA
TGCTATGTAAGAAAGGACTTGGTCTACACCAAAGCCAATCCAACGTTTCATCACTGGAAGGTCGATAAT
AGGAAGTTTGGACTTACTTTCAAAGCCCTGCTGATGCCCGAGCCTTGGACAGGGGAGTAAGGAAAGCA
ATCGAAGACCTTATAGAAGTTCAACAAGTCACTTCCACCATCCATAATGAAGCTGAGCTTGGCGAT
GATGACGTTTTTACAACAGCTACAGACAGTTCTTCTAATTCTCTCAGAAGAGAGCAACCTACTCGG
ACAATCTCCTCTCCACATCCTGTGAGCACCGGAGGATTTATACCCTGGGCCACCTCCACGACTCATA
CCCACAGACCACTATCACCTCGATCAGCCGATGCCAAGGCCCTACCGCCAGGTGAGCTTCCCGGACGAC
GACGAGGAGATCGTGCGCATCAACCCCGGGAGAAGATCTGGATGACGGGTACGAGGATTACCGGCAC
GCACCCGTCAGGGGAAGTACCCGACCCTCGGAGGACGCGGACTCCTCCTACGTGCGCTTCGCCAAG
GGCGAGGTCCCAAGCATGACTACAACCTACCCCTACGTGGACTCCTCAGACTTTGGCCTAGGCGAGGAC
CCCAAAGCCCGCGGGGGCAGCGTGATCAAGACGACGCCCTCCCGGGCAAGTCGCGGCGGGGAAGGAG
GACGGAGAGCGCTCGCGGTGCGTACTGCAGGGACATGTTCAACCACGAGGAGAACC CGCGGGCCAC
TGCCAGGACGCGCCCGACTCCGTGAGAATTCATCCCGGGTGAGCTGCATGTGGTGC CGCGGACAGC
ATGCTCTATCACTGTATGTCGGACCCCGAGGGAGACTATACAGACCCCTGCTCGTGCGATACTAGCGAC
GAGAAGTTTTGCCTCCGGTGGATGGCTCTATTGCCTTGTCTTTCCTGGCCCCCTGTATGTGCTGTAC
CTGCCCTTCGGGCTGCTACCACTGCGGAGTGATGTGCAGGTGCTGTGGCGGGAAGCACAAAGCGGCC
GCGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001128210

Insert Size: 1248 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001128210.1](#)

RefSeq Size: 4104 bp

RefSeq ORF: 1248 bp

Locus ID: 200734

UniProt ID: [Q7Z698](#)

Cytogenetics: 2p14

Protein Pathways: Jak-STAT signaling pathway

MW: 47 kDa

Gene Summary: SPRED2 is a member of the Sprouty (see SPRY1; MIM 602465)/SPRED family of proteins that regulate growth factor-induced activation of the MAP kinase cascade (see MAPK1; MIM 176948) (Nonami et al., 2004 [PubMed 15465815]).[supplied by OMIM, Mar 2008]
Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (b) has a shorter and distinct N-terminus compared to isoform a.