

## Product datasheet for **SC323032**

### Dermokine (DMKN) (NM\_001126056) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Dermokine (DMKN) (NM_001126056) Human Untagged Clone
Tag:	Tag Free
Symbol:	DMKN
Synonyms:	UNQ729; ZD52F10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC323032 representing NM\_001126056.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGTTCAGGGGCCCTGGCCTGCCTCTGCTGGCCCTCTGCTGGGCAGTGGGGAGGCTGGCCCC
CTGCAGAGCGGAGAGAAAGCACTGGGACAAATATTGGGGAGGCCCTTGGACATGGCCTGGGAGACGCC
CTGAGCGAAGGGGTGGGAAAGGCCATTGGCAAAGAGGCCGGAGGGGCAGCTGGCTCTAAAGTCAGTGAG
GCCCTTGGCCAAGGGACCAGAGAAGCAGTTGGCACTGGAGTCAGGCAGGTTCCAGGCTTTGGCGTAGCA
GATGCTTTGGGAACAGGGTTCGGGAAGCAGCCATGCTCTGGGAAACTGGGCACGAGATTGGCAGA
CAGGCAGAAGATGTCATTGACACGGAGCAGATGCTGTCCGGGCTCTGGCAGGGGGTGCCTGGCCAC
AATGGTGTGGGAACTTCTGGAGGCCATGGCATCTTTGGCTCTCAAGGTGGCCTTGGAGGCCAGGGC
CAGGGCAATCCTGGAGGTCTGGGACTCCGTGGTCCACGGATACCCGGAAACTCAGCAGGCAGCTTT
GGAATGAATCCTCAGGGAGCTCCCTGGGTCAAGGAGGCAATGGAGGGCCACAAACTTTGGACCAAC
ACTCAGGGAGCTGTGCCAGCCTGGCTATGGTTCAGTGAGGCCAGCAACCAGAATGAAGGTTGCACG
AATCCCCACCATCTGGCTCAGGTGGAGGCTCCAGCAACTCTGGGGATCCAGCACCGGCTCCTCCTCC
GGCAACCACGGTGGGAGCGGGGAGGAAATGGACATAAACCCGGGTGTGAAAAGCCAGGGAAATGAAGCC
CGCGGGAGCGGGAACTCTGGGATTAGGGCTTCAGAGGACAGGGAGTTTCCAGCAACATGAGGAAATA
AGCAAAGAGGGCAATCGCTCCTTGGAGGCTCTGGAGACAATTATAGGGGGCAAGGGTCGAGCTGGGGC
AGTGGAGGAGGTGACGCTGTTGGTGGAGTCAATACTGTGAACCTCTGAGACGCTCCTGGGATGTTTAACT
TTTGACACTTTCTGGAAGAATTTAAATCCAAGCTGGGTTTCATCAACTGGGATGCCATAAACAAGAAC
CAGGTCGGCCCCCAGCACCCGAGCCCTCTCTACTTCAGCCGACTCTGGGAGGATTTCAAACAGAAC
ACTCCTTTCTCAACTGGAAAGCAATTATTGAGGGTGGGACGCGTCATCACTGCAGAAACGTGCAGGC
AGAGCCGATCAGAACTACAATTACAACAGCATGCGTATCCCACTGCCTATGGTGGGAAGTACTCAGTC
AAGACCCCTGCAAAGGGGGAGTCTCACCTTCTCCTCGGCTTCCCGGTGCAACCTGGCCTGCTGCAG
TGGGTGAAGTTTTGGTAG
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

**Restriction Sites:**

Sgfl-RsrII

**Plasmid Map:**



**ACCN:** NM\_001126056

**Insert Size:** 1398 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001126056.2</a></u>
<b>RefSeq Size:</b>	1913 bp
<b>RefSeq ORF:</b>	1398 bp
<b>Locus ID:</b>	93099
<b>UniProt ID:</b>	<u><a href="#">Q6E0U4</a></u>
<b>Cytogenetics:</b>	19q13.12
<b>MW:</b>	46.9 kDa
<b>Gene Summary:</b>	<p>This gene is upregulated in inflammatory diseases, and it was first observed as expressed in the differentiated layers of skin. The most interesting aspect of this gene is the differential use of promoters and terminators to generate isoforms with unique cellular distributions and domain components. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jun 2010]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon and contains three additional alternate in-frame exons in the coding region, compared to variant 2. This results in a shorter protein (isoform 3, also referred to as isoform beta-2), compared to isoform 2.</p>