

Product datasheet for SC323047

LIPM (NM_001128215) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LIPM (NM_001128215) Human Untagged Clone
Tag:	Tag Free
Symbol:	LIPM
Synonyms:	bA304I5.1; LIPL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC323047 representing NM_001128215. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGTTGGAAACCTTGTC AAGACAGTGGATTGTCTCACACAGAATGGAAATGTGGCTTCTGATTCTGGTG
GCGTATATGTTCCAGAGAAATGTGAATTCAGTACATATGCCAACTAAAGCTGTGGACCCAGAAGCATTTC
ATGAATATTAGTGAATCATCCAACATCAAGGCTATCCCTGTGAGGAATATGAAGTCGCAACTGAAGAT
GGGTATATCCTTTCTGTTAACAGGATTCCTCGAGGCCTAGTGCAACCTAAGAAGACAGGTTCCAGGCCT
GTGGTGTACTGCAGCATGGCCTAGTTGGAGGTGCTAGCAACTGGATTTCCAACCTGCCCAACAATAGC
CTGGGCTTCATTCTGGCAGATGCTGGTTTTGACGTGTGGATGGGGAACAGCAGGGGAAACGCCTGGTCT
CGAAAACAAGACACTCTCCATAGACCAAGATGAGTTCTGGGCTTTCAGTTATGATGAGATGGCTAGG
TTTGACCTTCTGCAGTGATAAACTTTATTTTGCAGAAAACGGGCCAGGAAAAGATCTATTATGTCCGGC
TATTCACAGGGCACCCACCATGGGCTTTATTGCATTTTCCACCATGCCAGAGCTGGCTCAGAAAATCAAA
ATGTATTTTGCTTTAGCACCCATAGCCACTGTTAAGCATGCAAAAAGCCCCGGGACCAAAATTTTGTTG
CTGCCAGATATGATGATCAAGGGATTGTTTGGCAAAAAGAAATTTCTGTATCAGACCAGATTTCTCAGA
CAACTTGTATTACCTTTGTGGCCAGGTGATTCTTGATCAGATTTGTAGTAATATCATGTTACTTCTG
GGTGGATTCAACCAACAATATGAACATGAGCCGAGCAAGTGTATATGCTGCCACACTTTGACTGG
ACATCTGTGCAAAATATTCTACTGGAGCCAGGCAGTGAATTTGGTGAATCCGGGCATTTGACTGG
GGGAGTGAGACCAAAAATCTGGAAAAATGCAATCAGCCAACCTCTGTAAGGTACAGAGTCAGAGATATG
ACGGTCCCTACAGCAATGTGGACAGGAGGTCAGGACTGGCTTTCAATCCAGAAGACGTGAAAATGCTG
CTCTCTGAGGTGACCAACCTCATCTACCAATAAGAATATTCCTGAATGGGCTCACGTGGATTTTCATCTGG
GGTTTGGATGCTCCTCACCGTATGTACAATGAAATCATCCATCTGATGCAGCAGGAGGAGACCAACCTT
TCCAGGGACGGTGTGAGGCCGATTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001128215

Insert Size: 1272 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001128215.1](#)

RefSeq Size: 1484 bp

RefSeq ORF: 1272 bp

Locus ID: 340654

UniProt ID: [Q5VYY2](#)

Cytogenetics: 10q23.31

Protein Families: Druggable Genome

MW: 48.2 kDa

Gene Summary: Plays a highly specific role in the last step of keratinocyte differentiation. May have an essential function in lipid metabolism of the most differentiated epidermal layers.
[UniProtKB/Swiss-Prot Function]