

Product datasheet for **SC323054**

TBCEL (NM_001130047) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TBCEL (NM_001130047) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBCEL
Synonyms:	EI; LRRC35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC323054 representing NM_001130047. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGATCAACCTAGTGAAGAAGTTTCATGCAAGTATTATGTGAAAAATATAGTCCTGAAAAATTTTCCT
TATCGCGTGGCCGGGGATGGGAGTCCATGTCCAGCCACACCTCAGGGCTCTCCTATGAAAGATCGC
CTCAACCTCCCAAGTGTACTAGTGTGAACAGCTGTGGAATAACCTGTGCAGGAGATGAAAAAGAAATT
GCTGCTTCTGCGCTCATGTGTCGGAAGTATGATCTTTCTGACAACAACTCGAAGACTGGCATGAGGTC
AGTAAAAATTGTGCAAAATGTTCCCTCAGTTGGAGTTTCTAAACCTGAGTTCCAACCCTCTGAATTTGTCG
GTTTTAGAAAGAACATGTGCTGGGTCCTTCTCTGGGGTTCGCAAATTGCTCTCAACAACAGCAAAGCT
TCTTGGGAGACGGTCCACATGATACTACAGGAGTTACCAGATTTGGAGGAGCTTCTCCTGTGCCTTAAT
GACTATGAAACAGTGTCTTGTCTTCTATTTGCTGTCATTCTCTTAAGCTACTACATATAACAGACAAT
AACCTCCAAGACTGGACTGAAATACGAAAGTTAGGAGTTATGTTTCTTCACTGGATACCCTCGTCTG
GCCAACAAATCATTTGAATGCTATTGAGGAGCCTGATGATTCATTGGCCAGGTTGTTTCTAATCTTCGA
TCCATCAGCCTCCACAAGTCAGTTTGCAGTCTGGGAAGACATTGATAAACTAAATTCATTTCCAAA
CTGGAAGAAGTGAGATTGTTAGGAATCCTCTTCTGCAGCCATATACCACCGAGGAGCGAAGGAAATTG
GTAATAGCCAGATTGCCATCAGTTTCAAACCTAATGGCAGCGTTGTTACTGATGGTGAACGAGAAGAT
TCTGAGAGATTTTTTATTCGTTACTATGTGGATGTTCCACAGGAAGAAGTCCATTTCAGGTATCATGAA
CTGATCACTAAATATGGGAAGTTGGAGCCTTGGCAGAAGTGGACCTAAGACCCAGAGCAGTGCAAAA
GTAGAAGTCCACTTTAACGATCAGTGGAAGAAATGAGCATTTCGCTGGACCAAACAGTGGCAGAATA
AAGAAACAGTTAAAACTCTAGTACAATTACCCACAAGCAACATGCTTCTACTATTTTACCATGAA
GCACCCTTTGGCCAGAGGAAATGAAGTACAGCTCTCGGGCATTGCATTCCTTTGGCATTAGGGATGGA
GATAAAATTTACGTGGAATCCAAAACAAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001130047

Insert Size: 1275 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130047.1](#)

RefSeq Size: 3914 bp

RefSeq ORF: 1275 bp

Locus ID: 219899

UniProt ID: [Q5QJ74](#)

Cytogenetics: 11q23.3

Protein Families: Druggable Genome

MW: 48.2 kDa

Gene Summary: Acts as a regulator of tubulin stability.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.