

Product datasheet for SC323060

C13orf31 (LACC1) (NM_001128303) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C13orf31 (LACC1) (NM_001128303) Human Untagged Clone
Tag:	Tag Free
Symbol:	LACC1
Synonyms:	C13orf31; FAMIN; JUVAR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC323060 representing NM_001128303. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTGAAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGAAGCTGTTTTGATTGATCTTTTTGGTTTGAATTTGAACTCTCAAAAAACTGCCATCAGACA
TTACTGAAGACTTTGAATGCTGTCCAATACCACCATGCTGCCAAGGCCAAGTTTCTGTATAATGTGT
TGCAGTAACATCAGCTATGAAAGGGATGGAGAACAAGATAATTGTGAAATAGAAACAAGCAATGGATTA
TCAGCTCTCTTGGAAAGATTTGAGATTGTTAGCTGTCCCAGCATGGCTGCCACTTTGTATACCATTAAA
CAGAAAAATTGATAAAAAATCTGAGCAGCATTAAAGTAATTGTACCCAGGCACAGGAAGACATTAATG
AAAGCTTTTATTGATCAACTCTTCACTGATGTTTACAATTTTGAATTTGAAGATTTGCAAGTGACTTTT
AGGGGAGGGCTTTTTAAACAGTCCATTGAAATAAACGTAATCACAGCTCAAGAACTAAGAGGAATTCAG
AATGAAATAGAAACATTTTTGAGAAGTCTGCCAGCACTGAGAGGAAAATTAACATTTATCACTTCTTCT
TTGATCCCAGATATTTTCATACATGGATTTACTACAAGAACAGGTGGGATATCTTATATACCAACTCTT
AGCTCATTCAATCTTTCAGTAGTTCCAAACGGAGAGATCCCAAGGTAGTGGTTCAAGAAAATCTGCGT
AGGTTGGCGAATGCTGCAGGATTTAATGTGGAGAAATTTACCGAATAAAGACTCATCATTCCAATGAC
ATCTGGATTATGGGAAGAAAGGAGCCTGACTCTTATGATGGAATAACCACAAATCAGAGAGGAGTCAACA
ATAGCAGCTCTTGGTGCAGACTGTATACCGATAGTTTTGCAGATCCAGTCAAAAAAGCATGTGGGGT
GCTCAGCTGGTTGGAAGGACTTTGTTGGGTGTTGCTATGGCTACAGTGAATGCTATGATAGCAGAA
TATGGCTGCAGTTTGGAAAGACATTTGTTGTTGACTTGGACCTTCAGTAGGACCTTGCTGTTTTACTCTT
CCAAGGAATCAGCAGAGGCATTTTCATAATCTTCACTCCTGCATGTGTACAACTATTTGATTCACCAAT
CCCTGTATCGACATCCGTAAGCCACAAGGATTCTTCTAGAACAGGGAGGAATCTTCCACAGAATATT
CAGGACCAGAACCAAGATCTCAACCTCTGTACATCTTGCCATCCTGACAAGTTTTCTCCCATGTCCGA
GATGGCCTTAATTTTGGTACACAGATTGGCTTCATATCAATTAAGAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001128303

Insert Size: 1293 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001128303.1](#)

RefSeq Size: 4288 bp

RefSeq ORF: 1293 bp

Locus ID: 144811

UniProt ID: [Q8IV20](#)

Cytogenetics: 13q14.11

MW: 47.8 kDa

Gene Summary: This gene encodes an oxidoreductase that promotes fatty-acid oxidation, with concomitant inflammasome activation, mitochondrial and NADPH-oxidase-dependent reactive oxygen species production, and bactericidal activity of macrophages. The encoded protein forms a complex with fatty acid synthase on peroxisomes and is thought to be modulated by peroxisome proliferator-activated receptor signaling events. Naturally occurring mutations in this gene are associated with inflammatory bowel disease, Behcet's disease, leprosy, ulcerative colitis, early-onset Crohn's disease, and systemic juvenile idiopathic arthritis. [provided by RefSeq, Apr 2017]

Transcript Variant: This variant (1) encodes the longest isoform (1). Variants 1-7 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.