

## Product datasheet for **SC323061**

### GCNT7 (NM\_080615) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	GCNT7 (NM_080615) Human Untagged Clone
Tag:	Tag Free
Symbol:	GCNT7
Synonyms:	C20orf105; dj1153D9.2; gcnt
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC323061 representing NM_080615. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCCAGCTTCGAGCCACAAAGTCTGGACTTGTCTGTGAGGGCAGTCATTTGCATCTTCATTTTTCTT
TACTTAAGGAATCCAACCTCTGCAGAATCAGAAGAAGAACCTGCCAACCAGAAGTAGTAGAATGTGGC
TTTTACCCAGATGAACTGTGTTCCGCTTTATTTGAAGGGAAAGGGCGGCCCCCAAAATTGCAAAATTT
TGCAAAACCCCTCATAAATCTGAAATACATGCTCATTTACACACACCAGGAACTGCTCCAGGATTTCT
CGGGGGCTGCATTTATAACCAGACCCCTGTCTGCAGAAGAGGGCGATTTCTCTTTGGCATATATTATA
ACTATTCATAAGGAGCTGGCCATGTTTGTGCAGCTTCTCAGAGCTATTTATGTACCTCAAATGTTTAT
TGTATTCATGTTGATGAAAAGGCCCAATGAAGTATAAGACTGCTGTGCAAACCTTGGTTAACTGTTTT
GAAAACGTTTTTATTTCTCCAAGACAGAGAAGGTGGCTTATGCTGGCTTTACAAGACTACAGGCAGAT
ATTAATTTGATGAAAGTTCTAGTGCATTCTAAATTTCAATGGAATATGTCATCAATCTTTGTGGACAG
GATTTTCCCATCAAAACCAACAGAGAAATCATACACTACATCAGAAGCAAATGGAGTGATAAAAAATATT
ACTCCTGGAGTAAATCCAACCATTCACATTAATCCAAGACAAGTCAAAGTCATCTCGAATTCGTTCCCT
AAAGGAAGTATCTATGCACCTCCAATAACAGATTCAAAGACAAACCACCCATAACTTAACCATTTAT
TTTGAAGTGTCTACTATGTACTTACAAGGAAGTTTGTAGAGTTCATACTGACTGACATCCATGCAAAA
GACATGCTTCAGTGGTCCAAGACATCCGAGCCAGAGCAACTACTGGGTGACCCCTGAACCGACTA
AAAGATGCTCCAGGTGCTACACCAACGCTGGCTGGGAAGGAAATGTTTCGAGCCATTAAGCGGAAAAGT
GAGGAAGGAAATGTTTCATGATGGATGTAAGGCCGCTACGTCGAAGACATCTGCGTATACGGACCAGGA
GACCTGCCGTGGCTCATTAGTGCCTTCTCTGTTTGCCAACAAATTCGAACCCTCGACAGACCCGCTT
GTGTTTACCTGCTTAGAGCGGGCCACAGACTTCAAGTGTGAGGCAGGCAGAAGTTCCTATAGAGCCA
CACTGGCATTTCACAGCAGAGTCATTTCAATATGAGACTGAACCGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_080615
<b>Insert Size:</b>	1293 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_080615.1</u>
<b>RefSeq Size:</b>	2064 bp
<b>RefSeq ORF:</b>	1293 bp
<b>Locus ID:</b>	140687
<b>UniProt ID:</b>	<u>Q6ZNI0</u>
<b>Cytogenetics:</b>	20q13.31
<b>MW:</b>	49.3 kDa
<b>Gene Summary:</b>	Glycosyltransferase.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1, coding) differs at a polymorphic nucleotide compared to the reference genome sequence and represents the protein coding allele.