

## Product datasheet for **SC323101**

### TC2N (NM\_001128596) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TC2N (NM_001128596) Human Untagged Clone
Tag:	Tag Free
Symbol:	TC2N
Synonyms:	C2CD1; C14orf47; MTAC2D1; Tac2-N
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC323101 representing NM\_001128596.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAACAGAATTTATAAAGAGTTGCTGTGGAGGATGTTTCTATGGTGAACAGAAAAACACAACCTT
TCTGTGGAAAGAGATTTTAAAGCAGCAGTCCCAAATAGTCAAAATGCTACTATCTCTGTACCTCCATTG
ACTTCTGTTTCTGTAAAGCCTCAGCTTGGCTGTACTGAGGATTATTTGCTTCCAAATTACCATCTGAT
GGCAAAGAAGTACCATTTGTGGTGCCCAAGTTTAAATTATCTTACATTCAACCCAGGACACAAGAACT
CCTTCACATCTGGAAGAAGTGAAGGATCTGCCAGAGCATCTTTTGGAGATCGAAAGGTAGAAGTTC
AGTTCATCCCAGCACGGACCTAGCTATGATGTGTATAACCCATTCTATATGTATCAGCACATTTACCT
GATTTGAGTCGACGCTTCTCCCGTTTCAAGAAGTACGAGACTGTATGGATCGGTTTGTGATTTAAGG
ACGAACAACTTCCCGTTCCCGTGGGCTAAGCAAATCTATGTTTGTATCTACAACTCATCTCAGCGA
TTCATCCAGAGACATGATTCATTGTCCAGTGTACCCAGTAGTTCTTCTCAAGGAAAAATCTCAGGGG
AGTAACAGAAGCCTGGATACAATTACTCTATCAGGAGATGAAAGGACTTTGGGAGACTGAATGTGAAA
TTGTTTTATAATCTTCCAGTAGAACAGATCTGGATCACAGTTTTACAGTGCAGAGATTTAAGTTGGCCC
TCTAGTTATGGAGACACTCCTACTGTTTCTATAAAGGAATACTTACATTGCCAAACCAGTGCATTTT
AAATCTTCAGCCAAGGAAGTTCCAACGCTATTGAATTTATGAAACGTTTGTATTTGCTATTAACCTT
CAAAATCTACAACTGTAAAGACTTGTATTTAAGATTCAAACCCAGACTCCCAGGAAGAAAACCATTTGGA
GAATGCTCAATGTCACTCAGAACCCTTAGCACACAGGAAATGGATTACTCTTTGGATATAACACCACCT
TCAAAAATTTCTGTTTCCCATGCAGAACTTGAATTGGGGACTTGTTCGAGTAAATAGCAGAATT
CAGTTACAAATCTTGAGGCACGGTACCTTCCAAGCTCATCAACACCTCTGACTTTGAGTTTTTCGTG
AAGGTGGGAATGTTAGCTCGGGAGAGTTGATTTATAAGAAAAAGACGCTTACTGAAGCCTCCAAAT
GGAAGAGTCAAGTGGGGAGAGACTATGATTTTTTCCACTTATACAGAGTAAAAAGAAATTGTTTTCTC
ATTAAGCTTTACAGTCGAAGCTCTGTAAGAAGAAAACTTTGTGGGCCAGATTTGGATAAGTGAAGAC
AGTAATAACATTGAAGCAGTGAACCAGTGGAAAGAGACAGTAATAAATCCAGAAAAGTTGTTATCAGG
TGGCACAATTAATCCATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001128596

<b>Insert Size:</b>	1473 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001128596.2</a></u>
<b>RefSeq Size:</b>	5153 bp
<b>RefSeq ORF:</b>	1473 bp
<b>Locus ID:</b>	123036
<b>UniProt ID:</b>	<u><a href="#">Q8N9U0</a></u>
<b>Cytogenetics:</b>	14q32.12
<b>MW:</b>	55.3 kDa