

Product datasheet for **SC323102**

ZNF222 (NM_001129996) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF222 (NM_001129996) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF222
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC323102 representing NM_001129996.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGATCGATTCAGGAGAAAAGAAGCCTGGGCGGAGAGCAGAGGAGGCAGTGACCTTCAAGGATGTGGCT
GTGATCTTCACTGAGGAGGAGCTGGGGCTGCTGGACCTGCCAGAGGAAGCTGTACCGAGATGTGATG
CTTGAGAACTTCAGGAACCTGCTCTCAGTGGGGCATCAACCATTCCATGGAGATACTTTCCACTTCCTA
AGGGAAGAAAAGTTTGGGTGATGGGGACAACAAGCCAAAGAGAAGGAATTTGGGAGGCAAGATCCAA
ACTGAGATGGAGACTGTTCCAGAAGCAGGAACACATGAAGAATTTTCTGCAAGCAAATTTGGGAACAA
ATTGCAAGTGACTTAACCAAGTCTCAAGATACCACCATAAGTAACCTCAGTTATTTGAACAAGATGAC
AACCCCTCCAGATTAAGCAAGACTATCTACAGTTCACACAAGAGAAAAACCTTTCCAGGGTGAAAT
TGTAACAGTTCTTCAGTGATGTTTCCTTCTTTGATCTTCCTCAGCAGTTATATTCAGGAGAGAAGTCT
CATACCTGTGATGAGTGTGAAAAAGCTTCTGTTACATCTCAGCCCTTCATATTCATCAGAGGGTCCAC
ATGGGAGTGAAATGCTATAAGTGTGATGTGTGGTAAGGAATTTAGTCAGAGCTCACGTCTGAAACT
CATCAAAGAGTCCACACTGGAGAGAAACCATTCAAATGTGAGCAGTGTGGGAAAGGCTTCAGATGTAGA
TCAGCACTTAAAGTTCATTGCAATTACACATGAGAGAGAAACCTTATAATTGTGAGAAATGTGGGAAG
GCTTTCATGCACAATTTCCAGCTTCAGAAACATCAGAAATTCATACTGGGAGAGGCCATTCAAATGT
GAAATATGTGGTAAGAGCTTCTGTCTTAGGTCAAGTCTTAATAGGCATTGCATGGTCCACACAGCAGAG
AACTGTACAAATCTGAAAAGTATGGAAGAGGTTTCATTGATAGGCTAGATTTGCATAAGCATCAGATG
ATTCATATGGGACAGAAACCATATAATTGTAAAGAATGTGGGAAGAGCTTCAAATGGTCCCTCATATCTT
TTGGTCCATCAACGAGTCCCACTGGAGAAAAGCCATACAAATGTGAGGAGTGTGGGAAGGGCTACATT
AGTAAGTCAGGTCTTGACTTCCACCATAGAACCACACGGGAGAGAGATCTTATAACTGTGATAACTGC
GGGAAGAGCTTTAGACATGCTTCTAGTATTTTGAATCATAAGAACTCCACTGCCAAAGAAAGCCATTG
AAATGTGAAGACTGTGGAAGAGGCTTGTATGCCGGTCATACTGTAAAGACCAACAAAGAGACCACAGT
GGAGAAAACCCATCCAAATGTGAGGACTGTGGGAAGCGCTACAAGAGGCGCTGAATCTGGATATAATT
TTATCATTATTTTAAATGACATATA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001129996

Insert Size: 1476 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001129996.1</u>
RefSeq Size:	1680 bp
RefSeq ORF:	1476 bp
Locus ID:	7673
UniProt ID:	<u>Q9UK12</u>
Cytogenetics:	19q13.31
Protein Families:	Transcription Factors
MW:	57.1 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and it encodes the longer protein (isoform 1).