

## Product datasheet for **SC323105**

### ATF7 (NM\_001130059) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ATF7 (NM_001130059) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATF7
Synonyms:	ATFA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC323105 representing NM\_001130059.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAGACGACAGACCGTTTGTGTGCAATGCCCGGGCTGTGGACAGAGATTTACAAACGAGGACCAC
CTGGCAGTTCATAAACACAAGCATGAGATGACATTGAAATTTGCCAGCCCGAACTGACTCAGTCATC
ATTGCAGATCAAACGCCTACTCCAAGTATTCCTGAAGAACTGTGAGGAGGTGGGACTCTCAATGAA
CTAGCTAGCTCCTTTGAACATGAATTCAAGAAAGCTGCAGATGAGGATGAGAAAAAGCAAGAAGCAGG
ACTGTTGCCAAAAAACTGGTGGCTGCTGCTGGGCCCTTGACATGTCTCTGCCTTCCACACCAGACATC
AAAATCAAAGAAGAAGAGCCAGTGGAGGTAGACTCATCCCACCTGATAGCCCTGCCTCTAGTCCCTGT
TCCCCACCACTGAAGGAGAAGGAGGTTACCCCAAAGCCTGTTCTGATCTTACCCCCACACCCACCATT
GTACGCTCTGGCTCCCTGCCTCTCCACTGGGCTATGATCCACTTCATCCAACCTTCCCTCCCAACC
TCTGTCATCACACAGGCTCCACCATCAAACAGGCAAATGGGGTCTCCCACTGGCTCCCTCCCTCTGTG
ATGCATCTTGCTAATGGACAGACCATGCCTGTGTTGCCAGGGCCTCCAGTACAGATGCCGTCTGTTATA
TCGCTGGCCAGACCTGTGTCCATGGTGCCCAACATTCCCTGGTATCCCTGGCCACCAGTTAACAGTAGT
GGTCCATTTCTCCCTCTGGCCACCCTATACCATCAGAAGCCAAGATGAGACTGAAAGCCACCCTAACT
CACCAAGTCTCCTCAATCAATGGTGGTGTGGAATGGTGGTGGTACTGCCAGCACCATGGTGACAGCC
CGCCAGAGCAGAGCCAGATTCTCATCCAGCACCTGATGCCCATCCCTGCCAGCCACAGGTCTCA
CCAGCTCAGCCCCACCCTAGTACTGGGGGGCGACGGCGGCACAGTAGATGAAGATCCAGATGAGCGA
CGGCAGCGCTTTCTGGAGCGCAACCGGGCTGCAGCCTCCCGTCCCGCCAAAAGCGAAAGCTGTGGGTG
TCCTCCCTAGAGAAGAAGGCCAAGAAGTCACTTCTCAGAACATTCAGCTGAGTAATGAAGTCACATTA
CTACGCAATGAGGTGCCCAGTTGAAACAGCTACTGTTAGCTCATAAAGACTGCCAGTCACTGCACTA
CAGAAAAAGACTCAAGGCTATTTAGAAAAGCCCAAGGAAAGCTCAGAGCCAACGGGTTCTCCAGCCCT
GTGATTACGACAGCTCAGCAACAGCCCTAGCAATGGCCTCAGTGTTCGCTCTGCAGCTGAAGCTGTG
GCCACCTCGTCTCACTCAGATGGCCAGCCAAAGGACAGAAGTGAAGTATGCCGATACAATCGCATGTA
ATCATGACCCACAGTCCAGTCTGCGGGCAGA TGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:**

NM\_001130059

<b>Insert Size:</b>	1485 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001130059.1</a></u>
<b>RefSeq Size:</b>	6600 bp
<b>RefSeq ORF:</b>	1485 bp
<b>Locus ID:</b>	11016
<b>Cytogenetics:</b>	12q13.13
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	53 kDa
<b>Gene Summary:</b>	<p>Plays important functions in early cell signaling. Binds the cAMP response element (CRE) (consensus: 5'-GTGACGT[AG][AG]-3'), a sequence present in many viral and cellular promoters. Activator of the NF-ELAM1/delta-A site of the E-selectin promoter. Has no intrinsic transcriptional activity, but activates transcription on formation of JUN or FOS heterodimers. Also can bind TRE promoter sequences when heterodimerized with members of the JUN family.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1, also known as ATFa3).</p>