

## Product datasheet for **SC323117**

### **CYPIVF11 (CYP4F11) (NM\_001128932) Human Untagged Clone**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids                                    |
| Product Name:             | CYPIVF11 (CYP4F11) (NM_001128932) Human Untagged Clone |
| Tag:                      | Tag Free   |
| Symbol:                   | CYP4F11  |
| Synonyms:                 | CYPIVF11   |
| Mammalian Cell Selection: | Neomycin   |
| Vector:                   | pCMV6-Entry (PS100001)                                 |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                   |



[View online »](#)


**Fully Sequenced ORF:** >SC323117 representing NM\_001128932.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCGAGCTGAGCCTGTCTGGCTGGCCCTCGGGCCGTGGCAGCATCCCCGTGGCTGCTTCTGCTG
CTGGTTGGGGCTCCTGGCTCCTGGCCCGCTCTGGCCTGGACCTACACCTTCTATGACAACTGCCGC
CGCTCCAGTGTTTTCTCAACCCCGAAACAGAAGTGGTTTTGGGGACACCAGGCCTGGTCACTCCC
ACGGAAGAGGGCATGAAGACATTGACCCAGCTGGTGACCACATATCCCCAGGGCTTTAAGTTGTGGCTG
GGTCTACCTTCCCCCTCCTCATTTTATGCCACCCTGACATTATCCGGCCTATCACCAGTGCCTCAGCT
GCTGTGCGACCCAAGGATATGATTTTCTATGGCTTCTGAAGCCCTGGCTGGGGGATGGGCTCCTGCTG
AGTGGTGGTGACAAGTGGAGCCGCCACCGTCGGATGTTGACGCCTGCCTTCCATTTCAACATCTTGAAG
CCTTATATGAAGATTTTCAACAAGAGTGTGAACATCATGCACGACAAGTGGCAGCGCCTGGCCTCAGAG
GGCAGCGCCAGACTGGACATGTTTGAACACATCAGCCTCATGACCTTGGACAGTCTGCAGAAATGTGTC
TTCAGCTTTGAAAGCAATTGTCAGGAGAAGCCAGTGAATATATTGCCGCATCTTGGAGCTCAGTGCC
TTTGTAGAAAAGAGAAACCAGCAGATTCTTGCACACGGACTTCTGTATTATCTCACTCCTGATGGG
CAGCGCTTCCGCAGGGCCTGCCACCTGGTGACAGACTTACAGATGCCGTATCCAGGAGCGGGCCTGC
ACCTCCCCACTCAGGATTTGATGATTTCTCAAGAACAAGGCAAAGTCCAAGACTTTAGACTTCATT
GATGTGCTTCTGCTGAGCAAGGATGAAGATGGGAAGGAATTGTCTGATGAGGACATAAGAGCAGAAGCT
GACACCTTCATGTTTGGGGCCATGACACTACAGCCAGTGGTCTCTCCTGGTCTATACCACCTTGCA
AAGCACCAGAATACCAGGAACAGTGGCCGCAAGAAGTGAAGAGCTTCTGAAGGACCGTGAACCTATA
GAGATTGAATGGGACGACCTGGCCAGCTGCCCTTCTGACCATGTGCATTAAGGAGAGCCTGCGGTTG
CATCCCCCAGTCCCGTCATCTCCGATGTTGCACGCAGGACTTGTGCTCCAGAGCGCCGCTCAGC
CCCAAAGGCATTGTCTGCCTCATCAATATTATCGGGATCCATTACAACCAACTGTGTGGCCAGACCCT
GAGGTCTACGACCCCTTCCGTTTTCGACCAAGAGAACATCAAGGAGAGGTACCTCTGGCTTTTATCCC
TTCTCGGCAGGGCCAGAAACTGCATCGGGCAGGCTTCGCCATGGCTGAGATGAAGGTGGTCTGGCG
CTCACCTGTGCACTTCCGCATCTGCCGACCCACACTGAACCCCGCAGGAAACCCGAGCTGATATTG
CGCGCAGAGGGTGGACTTGGCTGCGGGTGGAGCCCTGGGTGCGAACTCACAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001128932

|                               |   |
|-------------------------------|---|
| <b>Insert Size:</b>           | 1575 bp   |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>   |
| <b>RefSeq:</b>                | <a href="#">NM_001128932.1</a>  |
| <b>RefSeq Size:</b>           | 2973 bp   |
| <b>RefSeq ORF:</b>            | 1575 bp   |
| <b>Locus ID:</b>              | 57834   |
| <b>UniProt ID:</b>            | <a href="#">Q9HBI6</a>  |
| <b>Cytogenetics:</b>          | 19p13.12  |
| <b>Protein Families:</b>      | Druggable Genome, P450, Transmembrane   |
| <b>MW:</b>                    | 60.1 kDa  |
| <b>Gene Summary:</b>          | <p>This gene, CYP4F11, encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This gene is part of a cluster of cytochrome P450 genes on chromosome 19. Another member of this family, CYP4F2, is approximately 16 kb away. Alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) lacks an internal segment in the 5' UTR, as compared to variant 1. Both variants 1 and 2 encode the same protein.</p> |