

Product datasheet for **SC323122**

ZNF837 (NM_138466) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF837 (NM_138466) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF837
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC323122 representing NM_138466.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGCTCCAGCCCAGAAGGCTGGGCAGGGAGGACTCCCCAAGGCCGATGCCAGGGTGCCTCCGGG
GCCCGGAGAGAAGAGGCCCGAGGAGCCGAGGCCCTCGAAGAGGACCGAGCTGGGAGCCGCCCACTCAA
AAGGGGGACCTGCGTGGAGCGGGCGGTAGGACGACTCCCCAGGCGGGGGCTCCCGGGCTGCAGC
CTCGGGGTGAGCCCCGGCCGGGACCCGGCACAGCGCCGGGACCAGACCCCTCGTGCGGGAGCCGTGC
GGCCCCACCTCTTCTCAGAACCCTGAGCTGTTATCCCCGAGGGGCTGCAAGCCCGGAGGGCCCTGT
AGGAGCCAGCGCTGGCGGGGACTGCAGCAGGAACCTCTGCCTGGCGTGGCATCGCGGGCGCCCGCT
GGGAGACGCCACCCGTGTGTACCCCTGTCCGGAGCGGATCCAGAACCACCCCGGACTCAACTGTGT
GAGGTCCACACGACTGTTGGCCGTGCCAACCAGGACTGGCGCTCCGACCTGCCGAGGACCCCAAAG
CCGACCTCCCGGGGAGGAACCCCTTGGTGGAGCAGCCCGGGCTTGCCTGTGCGCGAGGCCTTTGCG
TGGAGGGCCCTGCGGATCCCCAGGAGAGGCTGCAGGCGACGGAGGAGCCCGTCCGTGTGCCCGTGC
GGGAAGCGCTTCCGCCCAACCAGCAGCAGCAGCGGGCAAGAGTCCCCAGTGTGTCTGAGTGCGGC
CAAACCTCGCGACCTCGCCCTATTGTCCCCGACCCCCCGGCCAGCGACTGTACGCTTGCAGCAGTGC
GGCAAGGCCTTACGCGCACCTCCAGCCTGCTGCAGCACCAGCGCATCCACAGGGCGAGCGGCCCTAC
GAGTGCGCCGAGTGCGGCAAGGCCTTCTGCGCTGCTCCGGCCTGTACCGCCACCAGAAGACGACTCG
GCCGAGCGCCACCGCGTGGCCCCGTCTGGCGCGGCGCGCTTCCGGCTGGGGTGCCCGCCCTGCGGG
GACTACAGCGAGCGGAGTCCCCGACGGGGTGGGAGCCGGGAGAACCGTACGAGTGCGCCGACTGC
GCCAAGGCCTTCGGGCTGTTCTCGCACCTCGTGGAGCACCAGCGCGTGCACACCGCGAGAAGCCCTAC
GCGTGCCCCGAGTGGGCAAGGCCTTCAACCAGCGCTCGAACCTGAGCCGGCACCAGCGCACTCACAGC
AGCGCCAAGCCCTACGCGTGCCCACTGTGCGAAAAGGCCTTCAAGGGCCGCTCGGGCCTGGTGCAACAC
CAGCGCGCGCACACCGCGAGCGGCCCTATGGCTGCTCCGAGTGCAGAAAGACCTTCCGCGGCTGCTCC
GAGCTGCGCCAGCAGCGCCTGCACTCGGGCGAGAAGCCCTACATCTGCCGCGACTGCGGCAAGGCC
TTCGTGCGCAACTGCAGCCTGGTGCACCCTGCGCACGACACGGGCGAGCGGCCCTACGCGTGCAGGA
GATTGCGGCCGCGCTTACGCCAACGCTCAACCTCAACGAGCACCAGGAGCGGACGGGGCGCGCC
GCGCCTGA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_138466

Insert Size: 1596 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_138466.1</u>
RefSeq Size:	1948 bp
RefSeq ORF:	1596 bp
Locus ID:	116412
UniProt ID:	<u>Q96EG3</u>
Cytogenetics:	19q13.43
MW:	58.1 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) uses an alternate splice site, compared to variant 1, and is protein-coding.