

Product datasheet for **SC323132**

OSGIN2 (NM_001126111) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OSGIN2 (NM_001126111) Human Untagged Clone
Tag:	Tag Free
Symbol:	OSGIN2
Synonyms:	C8orf1; hT41
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC323132 representing NM_001126111.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCCGTGTGGTGTGCCGCTGCTCCCTGGCCGGTCATTTAGAAAATATAGTGACACTGAAACTGAA
GGAGAGATTTTAAATTCCTTAGTGCAACTTTTGGTGACAACCTGGGGCGAAAAGTTAAAGCGATGCCA
TTAGTTGAAGAAACTCTTTACTGGAAGATTCGTCACTGACTTTTCTGTGGTAATAATAGGAAATGGA
CCCTCAGGAATATGCCTTTCTTATATGTTATCAGGCTACAGACCGTATTTATCATCAGAAGCAATACAC
CCAAATACAATCTTAAATAGTAAATAGAAGAAGCAAGACATCTTCCATTGTTGATCAGGACTTAGAA
TACTTGTCTGAGGGCCTTGAGGGCCGATCATCCAATCCAGTTGCAGTACTTTTCGATACACTTCTTCAT
CCAGATGCTGACTTTGGGTATGATTATCCATCCGTTTTGCATTGGAAATTAGAGCAACATCATTATATC
CCTCACGTAGTTCTTGGTAAAGTCCACCTGGTGGGGCTTGGCATAATATGGAAGGCTCCATGTTGACA
ATCAGCTTTGGAAGTTGGATGAACTACCTGGACTTAAATTTAAGGACTGGGTATCAAGTAAACGAAGG
AGCCTAAAAGGGGATCGAGTTATGCCAGAGGAAATAGCTCGCTACTATAAACATTATGTAAGTCAATG
GGTCTTCAGAAGAATTCAGAGAGAATACTTACATAACTTCCGTATCAAGACTCTACAGAGATCAAGAT
GATGATGATATCAAGACAGAGATATTTCAACAAAGCATTACAGATAGAGAAGTCAAACCTTTATCAAG
AGAAACTGGGAAATAGGGGTTATCAGCGAATAGCTGATGTTTCTCATGTTCCCTTCTGCCTCTTTGCT
GAGAATGTAGCGCTGGCAACTGGAACGCTGGATTCTCCTGCCATCTGGAATGGAAGGGGAAGATTTT
CCTTTTGTGTTTCATTCAATGCCTGAATTTGGAGCTGCTATAACAAAGGAAAGTTGCGTGGCAAAGTG
GATCCAGTGTTAATTGTAGTTCTGGGCTTACTGCCGCTGACGCAGTACTGTGTGCTTACAACAGTAAT
ATCCCTGTGATTCATGTGTTTCGCAGACGAGTAAGTATCAATGATGTTACTCAGTCATATTTCTGAGACTCAAT
AAGCTGTATCCTGAATATCATAAAGTCTATCATATGATGTACTCAGTCATATTTCTGAGACTCAAT
CTTTTATCTGATTATACCAGCTTTCCCGAGCACCGTGTGCTTTCCCTTTAAGTCGGACATGAAATGTGT
CTCAAAGCGTTTCTGGATTGAAGAAAATATTTAAGCTGTCTGCAGCAGTAGTATTGATAGGTTCTCAT
CCTAATCTGTCTTTTCTGAAGGATCAAGGGTGTACCTAGGCCATAAGTCAAGCCAGCCAATCACATGT
AAGGGTAACTCTGTGAAATAGATACATATACCTATGAGTGTATTAAGAAGCCAACCTTTTTCATTG
GGTCTTTGGTTGGAGACAATTTTGTTCGATTTTAAAGGGAGGGGCGCTGGGTGTTACACGCTGTTTA
GCTACAAGACAGAAGAAAAGCATTGTTTGTGAAAGAGGAGGAGATGGGATAGCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_001126111
Insert Size:	1650 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001126111.2</u>
RefSeq Size:	4281 bp
RefSeq ORF:	1650 bp
Locus ID:	734
UniProt ID:	<u>Q9Y236</u>
Cytogenetics:	8q21.3
MW:	61.7 kDa
Gene Summary:	<p>May be involved in meiosis or the maturation of germ cells.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) uses an alternate 5'UTR that includes a translation start site that is upstream and in-frame with that used for variant 2. The resulting protein (isoform 1) is longer and has a distinct N-terminus compared to isoform 2.</p>