

Product datasheet for **SC323156**

SPIRE1 (NM_001128627) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPIRE1 (NM_001128627) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPIRE1
Synonyms:	Spir-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC323156 representing NM_001128627.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAACAGAGGTCATTGAATCTTTGGGAATTATTATTTATAAAGCACTGGACTATGGTTGAAGGAG
AATGAAGAAAGGGAATTAAGCCCTCCCCTAGAGCAGCTTATCGATCACATGGCCAACACGGTGGAAGCT
GACGGTAGCAATGATGAGGGCTATGAGGCTGCAGAAGAAGGCCTGGGAGATGAAGATGAAAAGAGAAAA
ATCTCAGCTATTCGGTCATATAGAGATGTCATGAAGTTGTGTGCTGCTCATCTCCCTACTGAATCAGAT
GCACCAATCATTATCAGGCAGTATGTCGTGCACTGTTTGCAGAAACAATGGAGCTCCATACATTTCTG
ACCAAAATTAAGAGTGCAGAAAGAGAATCTTAAGAAGATTCAAGAAATGAAAAGAGCGATGAATCTAGC
ACAGACTTGGAAGAGCTGAAAACGCTGACTGGCAGCATTCTGGGTACAGGTGATGAGGGATTTGAGG
AATGGGGTAAAACCTAAGAAGGTCCAAGAGCGGCAGTACAACCCTTTGCCATTGAATATCAGCTCACC
CCTTATGAGATGTTAATGGATGACATTCGCTGCAAAAGATACACCTTGCAGAAAAGTGATGGTGAATGGT
GATATTCCTCGGTTAAAAAGAGTGCTCATGAAATCATCCTCGACTTCATCAGATCCAGACCTCCT
TTAAATCCAGTCTCAGCCAGAAAACCTGAAACCACTCCACCACGGCCACGGAGCCTCCATGAAAAGATA
TTAGAAGAAATTAAGCAGAAAAGAAAGCTGCGGCCTGTATCACCAGAGGAGATTAGAGCTAGCAGATTA
GATGTGACTACCCCTGAATCTACAAAGAATCTTGTGGAGTCATCTATGGTGAATGGAGGTTTGACATCA
CAAACAAAAGAAAACGGGTTAAGTACCTCACAGCAGGTGCCTGCACAGCGGAAGAAGCTCCTCAGAGCC
CCAACCTGGCCGAACCTGGACAGCTCTGAGTCTGAGGAAGAAACGCTGCACAAGTCGACCAGCAGCAGC
AGCGTGTCTCCCTCTTTCCCTGAAGAGCCAGTCTGGAGGCCGTGCCACAAGGAAGAAGCCTCCAAAA
TTCTGCCCATATCATCAACACCCAGCCAGAGAGACGGCAGCCACCCAGAGACGACATTCCATTGAA
AAGGAAAACGCCTACTAACGTGAGGCAGTTCCTGCCGCTTCCAGGCAGAGTTCCCGCTCTCTTGAGGAA
TTCTGCTACCCAGTGGAATGCCTCGCTCTTACTGTGGAAGAAGTGATGCATATTCGCCAGGTCCTGGTG
AAGGCAGAGCTGAAAAATACCAACAGTATAAAGACATCTACACCGCCTTAAAAAAGGAAAGCTCTGC
TTTTGTTGCCGAACAGGAGTTTTCTTCTTCACTTGGTCTTATACCTGTCAGTTCTGTAAGAGGCCG
GTGTGCTCACAGTGTGCAAAAAGATGCGGCTGCCCTCAAACCACTCCACTCTTCTATCTTTTCA
TTGGGACCTTCTGCTCTGCAAAGAGGGGAAAGTAGTATGAGGTCAGAAAAACCTCCACTGCCCATCAT
CGGCCACTTCGGAGCATTGCCAGGTTCTCCTCAAATCTAAGTCTATGGACAAATCAGATGAAGAACTC
CAGTTTCCAAAAGAGTTGATGGAGGACTGGAGCACCATGGAGGTGTGTGGACTGCAAGAAGTTCATT
TCGAAATCATCAGTTCAAGCCGGCGCAGTCTGGTGTGGCCAACAAAAGGGCCCGATTGAAAAGGAAA
ACGCAGTCTTTCTACATGTCCTCGCCAGGCCCTCGGAGTACTGCCCTCAGAGAGGACGATCAGTGAG
ATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001128627

Insert Size: 1869 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001128627.1](#)

RefSeq Size: 5100 bp

RefSeq ORF: 1869 bp

Locus ID: 56907

UniProt ID: [Q08AE8](#)

Cytogenetics: 18p11.21

Protein Pathways: Dorso-ventral axis formation

MW: 71.6 kDa

Gene Summary: Spire proteins, such as SPIRE1, are highly conserved between species. They belong to the family of Wiskott-Aldrich homology region-2 (WH2) proteins, which are involved in actin organization (Kerkhoff et al., 2001 [PubMed 11747823]).[supplied by OMIM, Mar 2008]
Transcript Variant: This variant (3) lacks an in-frame exon and differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (c) is shorter at the N-terminus and lacks an alternate segment compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.