

Product datasheet for **SC323160**

GSPT1 (NM_001130006) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GSPT1 (NM_001130006) Human Untagged Clone
Tag:	Tag Free
Symbol:	GSPT1
Synonyms:	551G9.2; eRF3a; ETF3A; GST1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC323160 representing NM_001130006.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATCCGGGCAGTGGCGGGCGGCGGCGGCGGCGGGGAGCAGCAGCGGCAGCAGCAGC
AGCGACTCGGCGCCTGACTGCTGGGACCAGCGGACATGGAAGCCCCGGGCGGGCCCTTGCAGCGGC
GGCGGCTCCCTGGCGGGCGGCGGCGGAGGCCAGCGGGAGAACCTCAGCGGGCCTTCAGCCGGCAACTC
AACGTC AACGCCAAGCCCTTCGTGCCAACGTCCACGCCGCCGAGTTCGTGCCGTCTTCCTGCGGGGC
CCGGCAGCGCCGACCCCCAGCTGGCGGCGCCCAATAAACACGGAGCCGGCAGCGGCGGGAGGC
CGTGCGGCACCTGTGGAATCCTCTCAAGAGGAACAGTCATTGTGTGAAGTTCAAATTCAGCTGTTAGC
ATGGAATTTCAGAACCTATTGAAAATGGAGAGACAGAAATGTCTCCAGAAGAATCATGGGAGCAGAAA
GAAGAAATAAGTGAAGCAGAGCCAGGGGTGGTTCCTTGGGAGATGGAAGGCCGACAGGAAAGTGCC
CATGAAATGATGGAGGAGGAAGAGGAAATCCAAAACCTAAGTCTGTGGTTGCACCGCCAGGTGCTCCT
AAGAAAGAGCATGTAATGTAGTATTTCATTTGGGCACGTAGATGCTGGCAAGTCAACCATTGGAGGACAA
ATAATGATTTGACTGGAATGGTTGACAAAAGGACGCTTGAAAAGTATGAAAGAGAGCTAAAGAGAAA
AACAGAGAAACTTGGTACTTGTCTTGGGCTTAGACACAAATCAGGAAGAACGAGACAAGGGTAAAAACA
GTAGAAGTGGGTCGTGCCTATTTTAAAACCGAAAAGAAGCATTTCACAATTCAGATGCCCTGGCCAC
AAGAGTTTTGTCCAAATATGATTGGTGGTGCCTCTCAAGCTGATTTGGCTGTGCTGGTAATCTCAGCC
AGGAAAGGAGAGTTTTGAAACTGGATTTGAAAAAGGAGGACAGACAAGAGAACATGCAATGTTGGCAAAG
ACAGCAGGTGTA AACACCTAATTTGTGCTAATTAATAAGATGGATGATCCAACAGTAAATGGAGCAAT
GAGAGATGAAGAATGTAAGGAGAACTAGTGCCATTTTAAAAAAGTTGGCTTCAATCCAAAAAAG
GACATTCATTTATGCCCTGCTCAGGACTTACTGGAGCAAACTCAAAGAGCAGTCGGATTCTGTCTCT
TGGTACATTGGATTACCGTTTATTCCATATCTGGATAATTTGCCGAACCTCAATAGATCAGTTGATGGA
CCAATCAGGCTGCCAATTGTGGATAAGTACAAGGATATGGGCACTGTGGTCTGGAAAAGCTGGAATCA
GGATCTATTTGTAAGGCCAGCAGCTTGTGATGATGCCAAACAAGCACAACTGGAAGTTCTTGGAATA
CTTCCGATGATGTAGAGACTGATACCGTAGCCCCAGGTGAAAACCTCAAATCAGACTGAAAGGAATT
GAAGAAGAGGAGATTCTCCAGGGTTTACTTTGTGATCCTAATAATCTTTGTCATTCTGGACGCACA
TTTGATGCCCAGATAGTGATTATAGAGCAGCAATCCATCATCTGCCCAGGCTATAATGCGGTGCTGCAT
ATTCATACCTGTATTGAGGAGGTGAAAATAACAGCCTTAATCTGCTTGGTAGACAAAAAATCAGGAGAA
AAAAGTAAGACCCGACCCGTTTTGTCAAACAAGATCAAGTATGCATTGCTCGCTTAAGGACAGCAGGA
ACCATCTGCCCTTGAGACCTTAAAGACTCCCTCAGATGGGTGCTTTCACCTTAAGAGATGAGGGTAAG
ACCATTGCAATTGAAAAAGTTCTGAAACTGGTTCCAGAGAAAGACTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001130006

Insert Size: 1911 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130006.1](#)

RefSeq Size: 7166 bp

RefSeq ORF: 1911 bp

Locus ID: 2935

UniProt ID: [P15170](#)

Cytogenetics:	16p13.13
Protein Families:	Transmembrane
MW:	68.6 kDa
Gene Summary:	<p>Involved in translation termination in response to the termination codons UAA, UAG and UGA (By similarity). Stimulates the activity of ETF1 (By similarity). Involved in regulation of mammalian cell growth (PubMed:2511002). Component of the transient SURF complex which recruits UPF1 to stalled ribosomes in the context of nonsense-mediated decay (NMD) of mRNAs containing premature stop codons (PubMed:24486019).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site in an internal exon, compared to variant 1, resulting in an isoform (2) that is 1 aa shorter than isoform 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The extent of this transcript is supported by transcript alignments.</p>