

Product datasheet for **SC323194**

PHTF2 (NM_001127357) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PHTF2 (NM_001127357) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHTF2
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001127357, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGTCCAAAGTCACAGATGCTATAGTCTGGTATCAAAAGAAGATTGGAGCATATGAT
CAACAAATATGGGAAAAATCTGTTGAACAGAGAGAAATCAAGTTTATTAACCTGGGGCTA
AGGAATAAACCAAGAAAAACAGCACATGTGAAACCAGACCTCATAGATGTTGATCTTGTA
AGAGGGTCTGCATTTGCAAAGGCAAAGCCTGAAAGTCTTGACTTCTCTGACCAGAAAG
GGAATTTTCGAGTTGTATTTTTCCCCTTTTTCTCCGGTGGTGGTTACAAGTAACATCA
AAGGTCATCTTTTTCTGGCTTCTTGTCCTTTATCTTCTTCAAGTTGCTGCAATAGTATTA
TTCTGCTCCACTTCTAGCCACACAGCATACCTCTGACAGAGGTGATTGGCCGATATGG
CTGATGCTGCTCCTGGAACTGTGCATTGCCAGATTGTTCCACAAGAACACCCAAACCT
CCTCTAAGTACAGGGGTAAAAGAAGAAGGAAATTAAGAAAAGCAGCCATTTGGAAGTA
CATAGGGAAGGAGATGGTTCTAGTACCACAGATAACACACAAGAGGGAGCAGTTCAGAAC
CACGGTACAAGCACCTCTCACAGCGTTGGCACTGTCTCAGAGATCTCTGGCATGCTGCT
TTCTTTTTATCAGGATCAAAGAAAGCAAAGAATCAATTGATAAATCAACTGAAACTGAC
AATGGCTATGTATCCCTTGATGGGAAGAAGACTGTTAAAAGCGGTGAAGATGGAATACAA
AACCATGAACCTCAGTGTGAAACTATTCGACCAGAAGAGACAGCCTGGAACACAGGAACA
CTGAGGAATGGTCTAGCAAAGATACCCAAAGGACAATAACAAATGTCTCTGATGAAGTC
TCCAGTGAGGAAGGTCTGAAACAGGATACTCATTACGTCGTCATGTGGACAGGACTTCT
GAAGGTGTTCTTCGGAATAGAAAGTCACACCATTATAAGAAACATTACCTAATGAGGAC
GCCCTAAATCGGGTACTAGTTGCAGCTCTCGTGTTCAAGTTCAGACAGGATTCTGAG
AGTGCAAGGCCAGAATCTGAAACAGAAGATGTGTTATGGGAAGACTTGTTACATTGTGCA
GAATGCCATTCATCTTGTACCAGTGAGACAGATGTGAAAATCATCAGATTAATCCATGT
GTGAAAAAAGAATATAGAGATGACCCTTTTCATCAGAGTCATTTGCCCTGGCTCCATGT
TCCCACCCAGGATTAGAAAAAATAAGTGCTATAGTATGGGAAGGTAATGATTGTAAGAAA
GCAGACATGTCTGTACTTGAAATCAGTGAATGATAATGAACAGAGTGAACAGCCATATA
CCAGGAATAGGATACCAGATTTTTGAAAATGCAGTCTCTCATACTGGGTTAACTCCA
TTTGTTTTCCGACTTTCTCAAGCTACAGACTTGGAACTCACAGCACATTCTGCTTCA
GAACCTTATGTGATTGCATTTGGTTCTAATGAAGATGCATAGTTCTTTCTATGGTTATA
ATAAGTTTTGGGTTCCGCTGTCTCTTGTGTGGATTTCTTTTTTTGCTCTGTGTAGCA
GAAAGAAGTTATAAACAGCGATTACTTTTTGCAAACCTTTGGACATTTAACATCTGCA
AGGAGGGCTCGAAAATCTGAGGTTCCCTATTTCCGGTTGAAGAAAGTACAGAATATAAAA
ATGTGGCTATCTCCGTTCCCTATCTTAAGCGTCGAGGTCCTCAGCGATCAGTTGATGTA
ATAGTTTCATCTGCTTTCTTATTGACTATCTCAGTTGATTTATCTGTTGTGCCAGCTA
CTTCATGTACACGAGATCTTCTTGATTGCTACTACAATTGGGAATTGGTAATCTGGTGC
ATCTCGTTAACCTTTTTCTCCTAAGATTTGTTACCCTTGATCAGAAACAAGTAAAAAA
TATAGTAATACCTCAATATTACTTACTGAACAGATAAACCTCTACTTGAAAATGGAGAAA
AAACCTAACAAAAAGGAGGAACTGACACTAGTGAATAATGTTTTAAAACCTGGCTACTAAA
CTGCTAAAGGAGTTGGACAGTCCTTTTAGATTATATGGGCTTACAATGAATCCGCTGCTT
TATAACATCACCCAGGTTGTTATCCTGTCAGCTGTTTCTGGTGTATCAGTGACTTGCTT
GGATTTAATTTAAAGCTATGGAAGATTAAGTCA
    
```

Restriction Sites: Please inquire

ACCN: NM_001127357

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001127357.1</u> , <u>NP_001120829.1</u>
RefSeq Size:	5025 bp
RefSeq ORF:	2256 bp
Locus ID:	57157
UniProt ID:	<u>Q8N3S3</u>
Cytogenetics:	7q11.23-q21.11
Protein Families:	Transcription Factors, Transmembrane
Gene Summary:	May play a role in transcription regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).