

Product datasheet for **SC323197**

SPIRE1 (NM_001128626) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPIRE1 (NM_001128626) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPIRE1
Synonyms:	Spir-1
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001128626, the custom clone sequence may differ by one or more nucleotides

```

ATGGCTCAGGCGGCTGGCCCGGCGGGCGGGGAGCCGCGGACTGAGGCAGTGGGCGGC
GAGGGGCCGCGGGAGCCCGGGCAGCCGGCGCGGGCCGGGGCTCCCGGACGCGCTG
AGCCTGGAGGAGATCCTGCGGCTGTACAACCAGCCATCAACGAGGAGCAGGCGTGGGCC
GTGTGCTACCAGTGTCTGCGTTCCCTGCGCGCCGCCCGCCCGCCAGCCCCGCCAC
CGTGTGCGCTCGGCCGCGCAGATCCGCGTCTGGAGGGACGGCGCCGTACCCTGGCGCC
GCGGCCGACGACGCGGGAGAGCCGCCAGTTGCGGGTAAATTGGGATATTCACAATGT
ATGAAAACAGAGGTCATTGAATCTTTGGGAATTATTATTATAAAGCACTGGACTATGGT
TTGAAGGAGAATGAAGAAAGGAATTAAGCCCTCCCTAGAGCAGCTTATCGATCACATG
GCCAACACGGTGAAGCTGACGGTAGCAATGATGAGGGCTATGAGGCTGCAGAAGAAGGC
CTGGGAGATGAAGATGAAAAGAGAAAAATCTCAGCTATTCGGTCAATAGAGATGTCATG
AAGTTGTGTGCTGCTCATCTCCCTACTGAATCAGATGCACCAAATCATTATCAGGCAGTA
TGTCGTGCACTGTTTGCAGAAAATGGAGCTCCATACATTTCTGACAAAATTAAGAGT
GCGAAAGAGAATCTTAAGAAGATTCAAGAAATGAAAAGAGCGATGAATCTAGCACAGAC
TTGGAAGAGCTGAAAACGCTGACTGGGACGATTCTGGGTACAGGTGATGAGGGATTTG
AGGAATGGGGTAAAACCTAAGAAGGTCCAAGAGCGGCAGTACAACCCCTTGCCCATTTGAA
TATCAGCTCACCCCTTATGAGATGTTAATGGATGACATTCGCTGCAAAAGATACACCTTG
CGAAAAGTGATGGTGAATGGTGAATTTCCCTCGGTTAAAAAGAGTGTCTATGAAATC
ATCCTCGACTTCATCAGATCCAGACCTCCTTAAATCCAGTCTCAGCCAGAAAATGAAA
CCAACCTCACACGCGCCACGGAGCCTCCATGAAAGAATATTAGAAGAAATTAAGCAGAA
AGAAAGCTGCGGCTGTATCACCAGAGGAGATTAGACGTAGCAGATTAGCAATGCGGCCA
CTTAGCATGTCTTACAGTTTTGACTTGTGACATGACTACCCCTGAATCTACAAAAGAAAT
CTTGTGGAGTCATCTATGGTGAATGGAGTTTGACATCACAAAACAAAAGAAAACGGGTTA
AGTACCTCACAGCAGGTGCTGCACAGCGGAAGAAGCTCCTCAGAGCCCCAACTCTGGCC
GAACTGGACAGCTCTGAGTCTGAGGAAGAAACGCTGCACAAGTCGACCAGCAGCAGCAGC
GTGTCTCCCTCTTCCCTGAAGAGCCAGTCTGGAGGCGGTGCCACAAGGAAGAAGCCT
CCAAAATCCTGCCCATATCATCAACACCCAGCCAGAGAGACGGCAGCCACCCAGAGA
CGACATTCATTGAAAAGGAAACGCCTACTAACGTGAGGCAGTTCCTGCCGCTTCCAGG
CAGAGTCCCGCTCTCTTGAAGAAATCTGCTACCCAGTGAATGCCTCGCTCTTACTGTG
GAAGAAGTGATGCATATTCGCCAGTCTGGTGAAGGCAGAGCTGAAAAATACCAACAG
TATAAAGACATCTACACCGCCTTGAAAAAGGAAAGCTCTGCTTTTGTGCCGAACCAGG
AGGTTTTCTTCTTCACTTGGTCTTATACCTGTCAAGTCTGTAAGAGGCCGGTGTGCTCA
CAGTGTGCAAAAAGATGCGGCTGCCCTCAAACCATACTCCACTCTTCCTATCTTTTCA
TTGGGACCTTCTGCTCTGCAAAGAGGGGAAAGTAGTATGAGGTGCAAAAACCCCTCCACT
GCCCATCATCGGCCACTTCGGAGCATTGCCAGGTTCTCTCAAATCTAAGTCTATGGAC
AAATCAGATGAAGAACTCCAGTTTCCCAAAGAGTTGATGGAGGACTGGAGCACCATGGAG
GTGTGTGGACTGCAAGAAGTTCAATTCGAAATCATCAGTTCAAGCCGGCGCAGTCTG
GTGTTGGCCAACAAAAGGCCCCGATTGAAAAGGAAAACGCAGTCTTCTACATGTCCTCG
CCAGGCCCTCGGAGTACTGCCCTCAGAGAGGACGATCAGTGAGATC

```

Restriction Sites: Please inquire

ACCN: NM_001128626

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001128626.1 , NP_001122098.1
RefSeq Size:	5444 bp
RefSeq ORF:	2271 bp
Locus ID:	56907
UniProt ID:	Q08AE8
Cytogenetics:	18p11.21
Protein Pathways:	Dorso-ventral axis formation
Gene Summary:	<p>Spire proteins, such as SPIRE1, are highly conserved between species. They belong to the family of Wiskott-Aldrich homology region-2 (WH2) proteins, which are involved in actin organization (Kerkhoff et al., 2001 [PubMed 11747823]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.</p>