

Product datasheet for **SC323202**

CD18 (ITGB2) (NM_001127491) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CD18 (ITGB2) (NM_001127491) Human Untagged Clone
Tag:	Tag Free
Symbol:	ITGB2
Synonyms:	CD18; LAD; LCAMB; LFA-1; MAC-1; MF17; MFI7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC323202 representing NM_001127491.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTGGGCTGCGCCCCACTGCTCGCCCTGGTGGGGTCTCTCCCTCGGGTGCCTCTCTCAG
GAGTGCACGAAGTTCAAGGTCAGCAGCTGCCGGGAATGCATCGAGTCGGGGCCCCGGCTGCACCTGGTGC
CAGAAGCTGAACTTACAGGGCCGGGGATCCTGACTCCATTGCTGCGACACCCGGCCACAGCTGCTC
ATGAGGGGCTGTGCGGCTGACGACATCATGGACCCACAAGCCTCGCTGAAACCCAGGAAGACCACAAT
GGGGGCCAGAAGCAGCTGTCCCAAAAAAGTGACGCTTTACCTGCGACCAGGCCAGGCAGCAGCGTTC
AACGTGACCTTCCGGCGGGCCAAGGGCTACCCCATCGACCTGTACTATCTGATGGACCTCTCCTACTCC
ATGCTTGATGACCTCAGGAATGTCAAGAAGCTAGGTGGCGACCTGCTCCGGGCCCTCAACGAGATCACC
GAGTCCGGCCGATTGGCTTCGGGTCTTCGTGGACAAGACCGTCTGCCGTTGTAACACGCACCCT
GATAAGCTGCGAAACCCATGCCCAACAAGGAGAAAGAGTGCCAGCCCCGTTTGCCTTCAGGCACGTG
CTGAAGCTGACCAACAACCTCAACCAAGTTTCAGACCGAGGTTCGGGAAGCAGCTGATTTCCGAAACCTG
GATGCACCCGAGGGTGGGCTGGACGCCATGATGCAGGTGCCCGCTGCCGGAGGAAATCCGCTGGCGC
AACGTCACGCGGCTGCTGGTGTTCGCACTGATGACGGCTTCCATTTCCGGGGCAGCGGAAGCTGGGC
GCCATCTGACCCCAACGACGGCCGCTGTACCTGGAGGACAATTGTACAAGAGGAGCAACGAATTC
GACTACCCATCGGTGGCCAGCTGGCGCACAAAGCTGGCTGAAAACAACATCCAGCCCATCTTCGCGGTG
ACCAAGTAGGATGGTGAAGACCTACGAGAAACTCACCGAGATCATCCCAAGTCAGCCGTGGGGGAGCTG
TCTGAGGACTCCAGCAATGTGGTCCATCTCATTAAAGATGCTTACAATAAACTCTCTCCAGGGTCTTC
CTGGATCACAACGCCCTCCCGACACCTGAAAGTCACCTACGACTCCTTCTGCAGCAATCCAGGTGACG
CACAGGAACCGCCAGAGGTGACTGTGATGGCTGCAGATCAATGTCGGATCACCTTCAGGTGAAAG
GTCACGGCCACAGAGTGATCCAGGAGCAGTCGTTTGTATCCGGGCGCTGGGCTTCACGGACATAGTG
ACCGTGCAGGTTCTCCCAAGTGTGAGTGCCGGTCCCGGACCAGAGCAGAGACCGCAGCCTCTGCCAT
GGCAAGGGCTTCTGGAGTGCAGCATCTGCAGGTGTGACACTGGCTACATTGGGAAAACTGTGAGTGC
CAGACACAGGGCCGAGCAGCCAGGAGCTGGAAGGAAGCTGCCGGAAGGACAACAACCTCCATCATCTGC
TCAGGGCTGGGGACTGTGTCTGCGGGCAGTGCCTGTGCCACACCAGCGACGTCCCGGCAAGCTGATA
TACGGGCAGTACTGCGAGTGTGACACCATCAACTGTGAGCGCTACAACGGCCAGGTCTGCGGCGGCCG
GGGAGGGGCTCTGCTTCTGCGGAAAGTCCGCTGCCACCCGGGCTTTGAGGGCTCAGCGTGCCAGTGC
GAGAGGACCACTGAGGGTGCCTGAACCCGCGGCTGTTGAGTGTAGTGGTCTGCGCGGTGCCGCTGC
AACGTATGCGAGTGCCATTCAGGCTACCAGTGCCTCTGTGCCAGGAGTGCCCGGCTGCCCTCACCC
TGTGGCAAGTACATCTCCTGCGCCGAGTGCCTGAAGTTCGAAAAGGGCCCTTTGGGAAGAACTGCAGC
GCGGCGTGTCCGGGCTGCAGCTGTGCAACAACCCGTGAAGGGCAGGACCTGCAAGGAGAGGGACTCA
GAGGGCTGCTGGGTGGCTACACGCTGGAGCAGCAGGACGGGATGGACCGCTACCTCATATGTGGAT
GAGAGCCGAGAGTGTGTGGCAGGCCCAACATCGCCGCCATCGTCGGGGCACCCTGGCAGGCATCGTG
CTGATCGGCATTTCTCTGCTGGTCTGGAAGCTCTGATCCACCTGAGCGACCTCCGGGAGTACAGG
CGCTTTGAGAAGGAGAAGCTCAAGTCCCAGTGGAAACAATGATAATCCCTTTTCAAGAGCGCCACCAG
ACGGTCATGAACCCCAAGTTTGTGAGAGTTAG
ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001127491

Insert Size: 2310 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001127491.2](#)

RefSeq Size: 2963 bp

RefSeq ORF: 2310 bp

Locus ID: 3689

UniProt ID: [P05107](#)

Cytogenetics: 21q22.3

Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Natural killer cell mediated cytotoxicity, Regulation of actin cytoskeleton, Viral myocarditis
MW:	84.8 kDa
Gene Summary:	<p>This gene encodes an integrin beta chain, which combines with multiple different alpha chains to form different integrin heterodimers. Integrins are integral cell-surface proteins that participate in cell adhesion as well as cell-surface mediated signalling. The encoded protein plays an important role in immune response and defects in this gene cause leukocyte adhesion deficiency. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2014]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode the same isoform (1).</p>