

Product datasheet for **SC323220**

FSTL5 (NM_001128427) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FSTL5 (NM_001128427) Human Untagged Clone
Tag:	Tag Free
Symbol:	FSTL5
Mammalian Cell Selection:	Neomycin
Vector:	<u>PCMV6-Neo</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001128427, the custom clone sequence may differ by one or more nucleotides

```

ATGTTTAAAGTGCCTGGTCAGTTGTCTTGGTTCTCGGATTCATTTTTCTGGAGTCGGAAGGA
AGGCCAACCAAGAAGGAGGATATGGCCTTAAATCCTATCAGCCTCTAATGAGATTGCGA
CATAAGGAAAAAATCAAGAAAGTTCAAGAGTCAAAGGATTTATGATTCAGGATGGCCCT
TTTGGATCTTGTGAAAAAAGTACTGTGGTTTGGGAAGACACTGTGTTACCAGCAGAGAG
ACAGGGCAAGCAGAATGTGCCTGTATGGACCTTTGCAAACGTCCTACAAACCTGTGTGT
GGATCTGACGGAGAATTCTATGAAAACCACTGTGAAGTGACAGAGCTGCTTGCCTGAAA
AAACAAAAGATTACCATTGTTCACAATGAAGACTGCTTCTTTAAAGGAGATAAGTGAAG
ACTACTGAATACAGCAAGATGAAAAATATGCTATTAGATTTACAAAATCAAAAATATATT
ATGCAAGAAAAATGAAAATCCTAATGGCGACGACATATCTCGGAAGAAGCTATTGGTGGAT
CAAATGTTTAAATATTTTGTGACAGACAGTAATGGACTGTAGATATTAATGAACTAACT
CAGGTGATAAAACAGGAAGAAGTGGCAAGGATCTCTTTGATTGTACTTTGTATGTTCTA
TTGAAATATGATGATTTTAAATGCTGACAAGCACCTGGCTCTTGAAGAATTTATAGAGCA
TTCCAAGTGATCCAGTTGAGTCTGCCAGAAGATCAGAACTAAGCATCACTGCAGCAACT
GTGGGACAAAAGTCTGTTCTGAGCTGTGCCATTCAAGGAACCCTGAGACCTCCCATTATC
TGGAAAAGGAACAATATTATTCTAAATAATTTAGATTTGGAAGACATCAATGACTTTGGA
GATGATGGGTCTTGTATATTACTAAGGTTACCACAACCTCACGTTGGCAATTACACCTGC
TATGCAGATGGCTATGAACAAGTCTATCAGACTCACATCTTCCAAGTGAATGTTCTCCTCA
GTCATCCGGGTGTATCCAGAGAGTCAGGCTAGAGAGCCTGGGGTAACTGCCAGTCTTAGG
TGCCATGCAGAGGGCATAACCAAGCCTCAGCTTGGCTGGTTGAAGAATGGAATTGATATT
ACACCAAAGCTTTCCAAACAACCTCACGCTTCAAGCAAATGGCAGTGAGGTTACATAAGC
AATGTGCGCTATGAAGATACTGGAGCATACACTGTATCGCAAAGAATGAAGCAGGAGTG
GATGAAGACATCTTCTCTTTTTTGTGGAAGACTCTGCTAGAAAAGACCCTAGCTAACATA
TTATGGAGAGAAGAAGGCTGGGAATTGGGAACATGTTCTATGTTTTTTATGAAGATGGA
ATCAAAGTGATACAACCCATAGAATGTGAATTTAGAGGCACATTAAGCCTAGTGAAGAG
CTCCTTGGATTTAGGATGAAGTCTGTCCCAAAGCTGAGGGAGATGAAGTTCAGAGGTGT
GTGTGGGCATCAGCTGTTAATGTCAAAGACAAGTTCATTTATGTTGCACAGCCAACCTTG
GACAGAGTCCCTTATTGTTGATGTGCAGTCCCAAAAAGTGTTCAGGCAGTGAGCACAGAC
CCTGTCCCAGTTAAATTACACTATGACAAATCACATGATCAGGTCTGGGTGCTAAGCTGG
GGTACCTTGGAGAAGACATCACCAACACTACAGGTAATTACCCTGGCCAGTGGGAATGTG
CCTCACACACGATCCACACCAACCAAGTGGGAAAGCAATTTGACAGAGTGGATGATTTT
TTCATTCCCACCACAACACTCATTATCACCCATATGAGGTTTGGATTTATTCTTCATAAA
GATGAAGCTGCACTACAAAAAATGATCTTGAAACCATGTCATACATCAAGACAATTAAC
TTGAAGGACTATAAGTGCCTTCCCTCAGTCAATGGCATATACACTTGGGAGGCTACTAC
TTCATTGGCTGCAAACTGACAGCACCGGAGCAGTTTCCCAACAGGTCATGGTGGACGGT
GTAAGTACTCAGTCAATGGGTTCAATAGTGATGTGACGGGCACTCCATATGTCTCTCCA
GATGGCCACTACCTTGTGACGATTAATGATGTGAAAGGCTTGTAAAGGTTTCAAGTACATT
ACCATCAGAGGAGAAAATACAGGAGGCTTTTGTATTTACACAAATCTGCACATATCTGAT
CTGGCATTTCACCATCCTTTACTGAAGCCCACCAATATAACATCTACGGTAGTTCAAGC
ACACAACTGATGTGCTCTTTGTGGAGCTCTTCTGGGAAGGTCAAGATGATAAAGAGT
CTCAAGGAACCACTCAAGGCAGAAGAATGGCCTTGGAAACCGGAAAAACAGGCAAATCCAG
GACAGTGGCTTGTGGTCAATACCTGTGACACCTTCCAAGGACTCTCTCTTATCCTA
GATGGACGACTCAATAAATTAACCTGTGAGATCACTGAAGTTGAAAAAGGAAATACAGTC
ATTTGGGTTGGAGATGCC
    
```

Restriction Sites: Please inquire
ACCN: NM_001128427
Insert Size: 4832 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001128427.1](#), [NP_001121899.1](#)

RefSeq Size: 4832 bp

RefSeq ORF: 2541 bp

Locus ID: 56884

UniProt ID: [Q8N475](#)

Cytogenetics: 4q32.2

Protein Families: Secreted Protein