

## Product datasheet for **SC323228**

### **PMS1 (NM\_001128143) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PMS1 (NM_001128143) Human Untagged Clone
Tag:	Tag Free
Symbol:	PMS1
Synonyms:	HNPCC3; hPMS1; MLH2; PMSL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC323228 representing NM\_001128143.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAAACAATTGCCTGCGGCAACAGTTCGACTCCTTTCAAGTTCTCAGATCATCACTTCGGTGGTCAGT
GTTGTAAAAGAGCTTATTGAAAACCTCTTGGATGCTGGTCCACAAAGCGTAGATGTTAAACTGGAGAAC
TATGGATTTGATAAAAATTGAGGTGCGAGATAACGGGGAGGGTATCAAGGCTGTTGATGCACCTGTAAAG
GCAATGAAGTACTACACCTCAAAAATAAATAGTCATGAAGATCTTGAAAAATTTGACAACCTACGGTTTT
CGTGGAGAAGCCTTGGGGTCAATTTGTTGTATAGCTGAGGTTTTAATTACAACAAGAACGGCTGCTGAT
AATTTTAGCACCCAGTATGTTTTAGATGGCAGTGGCCACATACTTTCTCAGAAACCTTCACATCTTGGT
CAAGGTACAACCTGTAACCTGCTTTAAGATTATTTAAGAATCTACCTGTAAGAAAGCAGTTTTACTCAACT
GCAAAAAAATGTAAGATGAAATAAAAAAGATCCAAGATCTCCTCATGAGCTTTGGTATCCTTAAACCT
GACTTAAGGATTGCTTTGTACATAACAAGATTTATCTCAGTGGATTTCTTCAAAGTGTGATGCAGAC
CACTCTTTCAGTCTTTCAACACCAGAAAGAAGTTTCATCTTCATAAACAGTCGACCAGTACATCAA
AAAGATATCTTAAAGTTAATCCGACATCATTACAATCTGAAATGCCTAAAGGAATCTACTCGTTTGTAT
CCTGTTTTCTTCTGAAAATCGATGTTCTACAGCTGATGTTGATGTAATTTAACACCAGATAAAAAGC
CAAGTATTATTACAAAATAAGGAATCTGTTTTAATTGCTCTTGAAAATCTGATGACGACTTGTATGGA
CCATTACCTAGTACAAATCTTATGAAAATAATAAAACAGATGTTCCCGCAGCTGACATCGTTCTTAGT
AAAACAGCAGAAACAGATGTGCTTTTTAATAAAGTGAATCATCTGAAAAGAATTATTCAAATGTTGAT
ACTTCAGTCATTCCATTCAAAAATGATATGCATAATGATGAATCTGAAAAAACACTGATGATTGTTA
AATCACCAGATAAGTATTGGTGACTTTGGTTATGGTCATTGTAGTAGTGAATTTCTAACATTGATAAAA
AACACTAAGAATGCATTTCAAGACATTTCAATGAGTAATGTATCATGGGAGAACTCTCAGACGGAATAT
AGTAAAACTTGTTTTATAAGTTCCGTTAAGCACACCCAGTCAGAAAAATGGCAATAAAGACCATATAGAT
GAGAGTGGGAAAATGAGGAAGAAGCAGGTCTTGAAAACCTCTCGAAAATTTCTGCAGATGAGTGGAGC
AGGGGAAATATACTTAAAAATTCAGTGGGAGAGAATATTGAACCTGTGAAAATTTTAGTGCCTGAAAAA
AGTTTACCATGTAAGTAAGTAATAATAATTATCCAATCCCTGAACAAATGAATCTTAATGAAGATTCA
TGTAACAAAAAATCAAATGTAATAGATAATAAATCTGAAAAGTTACAGCTTATGATTTACTTAGCAAT
CGAGTAATCAAGAAACCCATGTCAGCAAGTCTTTTTGTTCAAGATCATCGTCTCAGTTTCTCATA
GAAAATCCTAAGACTAGTTTAGAGGATGCAACACTACAAATTGAAGAAGTGTGGAAGACATTGAGTGAA
GAGGAAAAACTGAAATATGAAGAGAAGGCTACTAAAGACTTGGAACGATACAATAGTCAAATGAAGAGA
GCCATTGAACAGGAGTCACAAATGTCACTAAAAGATGGCAGAAAAAAGATAAAACCCACCAGCGCATGG
AATTTGGCCCAGAACACAAGTTAAAAACCTCATTATCTAATCAACCAAAACTTGATGAACCTCCTCAG
TCCCAAATGAAAAAAGAAGGAGTCAAAATATTAATAATGGTACAGATCCCCTTTTCTATGAAAACTTA
AAAATAAATTTTAAAGAAACAAAACAAAGTTGACTTAGAAGAGAAGGATGAACCTTGCTTGATCCACAAT
CTCAGGTTTCTGATGCATGGCTAATGACATCCAAAACAGAGGTAATGTTATTAATCCATATAGAGTA
GAAGAAGCCCTGCTATTTAAAAGACTTCTTGAGAAATCAATAAATCTCTGCAGAGCCACTGGAAAAGCCA
ATTATGTTAACAGAGAGTCTTTTTAATGGATCTCATTATTTAGACGTTTTATATAAATGACAGCAGAT
GACCAAAGATACAGTGGATCAACTTACCTGTCTGATCCTCGTCTTACAGCGAATGGTTTTCAAGATAAAA
TTGATACCAGGAGTTTCAATTAAGTAACTGAAAATTAAGTGAAGATAAGGAATGGCTAATTGTCTCCCATTC
TATGGAGTAGCAGATTTAAAAGAAATCTTAATGCTATATTAACAGAAATGCAAAGGAAGTTTATGAA
TGTAGACCTCGCAAAGTGATAAGTTATTTAGAGGGAGAAGCAGTGCCTATCCAGACAATTACCCATG
TACTTATCAAAGAGGACATCCAAGACATTATCTACAGAATGAAGCACCAGTTTGGAAATGAAATTTAAA
GAGTGTGTTTATGGTGGCCATTTTTTCATCATTTAACCTATCTTCCAGAAACTACATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001128143

**Insert Size:** 2682 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001128143.1](#)

**RefSeq Size:** 3421 bp

**RefSeq ORF:** 2682 bp

**Locus ID:** 5378

**UniProt ID:** [P54277](#)

**Cytogenetics:** 2q32.2

**Protein Families:** Druggable Genome, Transcription Factors

**MW:** 101.4 kDa

**Gene Summary:** This gene encodes a protein belonging to the DNA mismatch repair mutL/hexB family. This protein is thought to be involved in the repair of DNA mismatches, and it can form heterodimers with MLH1, a known DNA mismatch repair protein. Mutations in this gene cause hereditary nonpolyposis colorectal cancer type 3 (HNPCC3) either alone or in combination with mutations in other genes involved in the HNPCC phenotype, which is also known as Lynch syndrome. [provided by RefSeq, Jul 2008]  
Transcript Variant: This variant (2) is also known as PMS1 nirs variant 7.