

Product datasheet for **SC323234**

SLFN14 (NM_001129820) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLFN14 (NM_001129820) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLFN14
Synonyms:	BDPLT20
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC323234 representing NM_001129820.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGAGTCTCAAGACTGATACTGAAATGCCGTATCCTGAGGTAATAGTAGATGTGGGCAGAGTGATT
TTTGGAGAAGAAAACAGGAAGAAGATGACCAACAGCTGTTTGAAGAGATCTGAGAATTCTAGAATTATC
CGGGCTATATGTGCACTGTTAAATTCTGGAGGTGGTGTGATCAAAGCAGAGATTGATGATAAAACCTAT
AGTTACCAATGCCATGGGCTGGGACAGGATTTGAAAACCTCTTTTCAAAGCTCCTTCCTTCAGGTTCA
CAGAAATACCTTGACTACATGCAGCAGGGGCACAATCTCCTGATTTTTGTGAAGTCATGGAGCCAGAT
GTTTTAGCCTTCCACTAAGGATTTGCAGCTTGCCTCAATTTGTATCGGAGAGATGTGACTTCTGCT
ATCAACTTGAGTGTAGCAGTGCCTGGAGCTTCTCAGAGAGAAGGGGTTTAGAGCCCAAAGAGGAAGA
CCAAGGGTGAAGAAGTTGCATCCTCAGCAGGTTCTCAATAGATGCATTGAGGAGAGGAAGATATGAGG
ATATTGGCCTCAGAATTTTTTAAAAAGGACAAACTCATGTATAAGGAGAACTCAACTTTACTGAGTCA
ACACATGTTGAATTTAAAGGTTACCACCAAAAAGTCATACCTCGGATTAAGGAAATGCTGCCTCAT
TATGTTTCTGCATTTGCCAACACTCAAGGGGGATATGTCCTATTGGGGTGGATGATAAGAGCAAAGAA
GTGGTTGGATGTAAGTGGGAAAAAGTGAATCCTGACTTACTAAAAAAGAAATCGAAAATGCATAGAA
AAATTGCCTACATTCCACTTCTGCTGTGAGAAGCCAAAGGTAATTTCACTACAAAAATCCTGAATGTG
TACCAAAAAGATGTCTGGATGGTTATGTCTGTGTGATTCAAGTGGAGCCCTTCTGTTGCGTGGTGTTC
GCAGAGGCCCCAGATTCTGGATCATGAAAACAATCTGTACACCGGCTGACAGCTGAGCAGTGGGTG
GTCATGATGCTGGATACTCAGTCAGCTCCTCCAGTTTGGTCCACAGACTACAACCTTGCCTGATTTCA
TCAGCTTCATCTGCACGAAAAAGTCCCGGATATCCATAAAAAGTCCACAAATTTAAGGAGCTGTCGA
CGACATTTGTTCCAGTGACACAGGAAGGTACAATTTAAACCAGAATCCCTCTGTAAGAAGCTGTTTC
TCAGATCATAAAGAAGTGGAGGGTTAATGAAGACCTGATACATCCTTGTCTCAGGGGATTGTGATA
TTTTCTAGAAGCTGGGCTGGTGTGTTGGCTTCAGGAAAGAACAAGTGAATGCTGTGTGATGCTCCTG
ATAGCAGTTAACAGCCCGTGGTACTCTATACAATCTTAATAGACCCCAATTGGCCTGGAGGACTTGAA
TATGCCGAAACACAGCTCATCAGTTAAAGCAGAACTGCAAACTGTTGGTGGTTACACAGGGAAAGTG
TGCATCATTCCAAGGCTGATACACCTGAGCAGCACACAGAGTAGACCTGGTGGAGATCCCCCTGCGTTAC
CCCAGGCTCACAGGCTGTGTGATGAGGAGGAAATGGAAGACTTGTGCAGGCTCTGGTGGTGGTCTCC
CTGTCTCCAGATCTCTCTGAGTGACCAGATGGGCTGTGAATTTTTCAACTTGCTCATAATGGAGCAG
AGCCAGTTGCTTCTGAGAGTCTTCAGAAGACAGTGAATTTTCACTACTGCTTCCAGGAGTCAGG
AAGACAGCCCTAGCCATAAAGATCATGGAGAAAATTAAGGACTTGTTCCTGCAAAACAAAAGAGATC
CTCTATGTTTGTGAAAGCGACTCCCTAAAGGATTTTGTGACCCAACAACCACTGCCAAGCTGTGACC
AGGAAAATTTTCAAGGGGAGTTTCTAAAGATTAACACATAGTGATGGATGAGACTGAGAATTTTC
TGCAGCAAATATGGCAATTGGTACATGAAGGCTAAGAACATCACCCATCCAAAGGCGAAGGGGACTGGA
AGTGAAAACCTTACCATGGGATTTCTGGCTTTTCTTGACCCTTTTCAAATCCATCACGCAGATGTC
AATGGCCTTCCCCTCCATCTGCTCAGTTTCTCGAAAAACAATCACAGTGGATCCACTGTGCTCTG
GAAATAGCGAAGGTTATGAAAGAAGAAATGAAGAGGATCAAAGAAAATCCTCCCTCAACATGTCTCCA
GACACATTAGCATTGTTGAGCAGACTGCCTATGAGGAAGCAACGTGTGCCAGGCTCTGCTGGGTTG
TGTGAGACAAAGACTAATCTGACAAACAGAAACAATAGCTAACTATGTGGCAAGAAAATGTACAGCCTG
TTCCAGTGTGGCTATCTGCCAAAAGATATAGCAATTTCTGTGAGGAGAGGGGAGGACAGAGGACGCTAT
AGGCTTGCACTACTCAAAGCAATGGAATTAATTGAGACCCACAGACCATCAGAGGTTGTGTTAGCCCG
GCCACCGGTGTTGGGGAAGTCACATTGTTTTAGACAGTATTCAGCAATTTTCAGGCCTGGAGAGGACT
GTCGTGTTGGGCTTAGTCCAGAATGTGACCAGTCAGAGGAAATTCATAAGCTCTGCTTTGCTTCAAGA
GCCATTAACACCTCTACCTGCTTATGAAAAGAGGGCAGCCTACTGA
ACGGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001129820

Insert Size: 2739 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_001129820.1</u>
RefSeq Size:	2882 bp
RefSeq ORF:	2739 bp
Locus ID:	342618
UniProt ID:	<u>P0C7P3</u>
Cytogenetics:	17q12
MW:	103.9 kDa
Gene Summary:	The protein encoded by this gene plays an important role in platelet formation and function. Defects in this gene are a cause of thrombocytopenia with excessive bleeding. [provided by RefSeq, Jul 2016]