

Product datasheet for **SC323235**

ZNF717 (NM_001128223) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF717 (NM_001128223) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF717
Synonyms:	OB1; X17; ZNF838
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001128223, the custom clone sequence may differ by one or more nucleotides

```

ATGTTTCCAGTGTCTCTGGCTGTTTCCAAGAGCTACAAGAAAAGAATAAATCTCTGGAG
TTGGTGTCTTTGAGGAGGTAGCTGTGCACTTACCTGGGAGGAGTGGCAGGACCTGGAT
GATGCTCAGAGGACCTGTACAGGGACGTGATGCTGGAGACCTACAGCAGCCTGGTATCA
TTGGGGCATTACATTACCAAACTGAGATGATCTTCAAGCTAGAGCAAGGAGCAGAGCCA
TGGATAGTAGAAGAAACCCAAACCTGAGACTTTTCACTGTCCAGATCATTGATGACCTT
ATTGAAAGGAGCCATGAAAGTCATGATAGATTTTTCTGGCAAATTGTAATCACCACAGC
AACACATCAACTCAGGAGAGAGTTGAATTAGGAAAAACATTTAATTTGAACTCAAACCAT
GTTTTAAATCTGATTATAAATAATGAAACAGTTCAGGAATGAAGCCTGGGCAGTTAAT
GATTGCCAGAACATGCTTTTCCCTATTAAGCCTGGGGAGACACAGTCTGGAGAGAAACCT
CATGTCTGTGATATAACCAGGAGATCCACAGACATCATGAACATCTTACTCAGCATCAC
AAGATTCAAACCTGTCTGCAGACTTTTCAATGTAATGAACAAGGGAAAACCTTCAACACG
GAGGCAATGTTCTTTATACATAAGAGGGTTCATATAGTACAGACCTTTGGTAAATATAAT
GAATATGAGAAAGCCTGTAATAACTCAGCTGTTATTGTCCAAGTGATAACTCAGGTAGGA
CAGCCAACTTGCTGTAGAAAGTCTGACTTCACTAAACATCAGCAGACACACAGGAGAG
AAACCCTATGAATGTGTTGAATGTGAGAAACCTCCATTAGCAAAATCAGACCTTATGCTA
CAGTGCAAGATGCCTACTGAGGAAAAACCTTATGCCTGTAAGTGGTGTGAAAAATTGTTT
AGCTATAAGTCCAGCCTCATTATCCATCAGAGAATTACACAGGGGAAAAAGCCCTATGGA
TGCAATGAATGTGGTAAAACCTTTCGCCGTAAGTCATTCTCACTTTACATGAGAGAACT
CACACAGGGGATAAACCTTACAAATGATTGAATGTGGAAAACTTTTCACTGTAAGTCA
CTTCTCACTTTACATCAGAACTCACTCAGGGGAAAAAGCCCTATCAGTGTAGTGAATGT
GGAAAACTTTAGCCAGAAGTCATACCTCACAAATACATCATAGAACTCACACAGGGGAA
AAGCCCTATGCATGTGACCATTGTGAAGAAGCATTTAGCCATAAGTCAAGGCTTACTGTC
CATCAGAGAACACACAGGGGAAAAACCGTATGAATGTAATGAGTGTGAAAAACCTTTT
ATCAATAAGTCAAACCTCAGGTTACATCAGAGAACTCACACAGGGGAAAAACCTATGAA
TGCAATGAATGTGGGAAAAACGTTTACCCTAAGTCATTCTCACTATCCATCAATGGACT
CACACAGGGGAAAAACCTACGAATGCAATGAATGTGGGAAAAACCTTTCGCTGTAAGTCA
TTCCTCACTGTCCATCAGAGAACTCATGCTGGGAAAAACCATACGCATGTAACGAATGT
GGAAAAACATATAGCCACAAGTCATACCTTACAGTACATCAGAGAACTCACACAGGGGAA
AAACCCTATGAATGTAATGAATGTGAAAAATCCTTCACTGTAAGTCATTCTAACTATA
CATCAGAGAACTCATGCTGGCAAAAAACCTATGAATGTAATGAATGTGAGAAAAACCTTT
ATCAATAAGTTAAACCTTGGGATACACAAGAGAACTCACACAGGGGAAAGACCCTATGAA
TGTAATGAATGTGAAAAACCTTTCGTCAGAAGTCAAATCTCAGACCCATCAGGGAACCT
CACACAGGAGAGAAACCTTACGTATGTAATGAATGTGAAAAACCTTTCATCGCAAGTCA
TTCCTCACCATACACCAGAGAACTCACACGGGGAAAAACCGTATGGATGTAATGAATGTG
GAAAAACCTTTTGTGAGAAATCATACCTTATTATACATCAGAGAACTCACACGGGGGAAA
AGCCCTATGAATGTAATGAATGTGAAAAATCCTTTCATCAGAAGGCAATCTTCAGAAGC
ATCAAGGTATTCACACGGGGGAGAAACCTATGAATGTAGCAAATGTGAAAAACCTTGT
CAGAAAGTCAGTCCTCACTGTACATCATAGAACCATAACCGAGAAAAACCTTATGAATGT
AATGAGTGTGGGAAAAACCTTTTGTGACAAGTCAAACCTCAGTACGCATCAGGGGACTCAC
TCAGGAGAGAAACCTATGAATGTGATGAATGTAGGAAAACTTTTACGATAAGACAGTT
CTCACCATACATCAGAGAACTCACACAGGTGAGAAGCCATTTGAATGTAAGAATGTAGG
AAAACCTTCTCCAGAAAGTCAAACTCTTTGTACATCAGAGAACTCACACAGGGGAAAAA
CCCTTTAGATGTAATGAATGTAGGAAAACTTTCTCCAGAAAGTCAGGCCTCAGTATACAT
CAGAGAACACACAGGAGAAAAACCTTATGAATGTAAGGAATGTGGGAAAAACCTTTTGC
CAGAAAGTCACACCTCAGCAGGCATCAACAAACCCATATAGGAGAGAAATCTGACGTAGCT
GAGGCAGGCTATGTGTTCCCTCAAAATCACTCTTTTTCCCT

```

Restriction Sites: Please inquire
ACCN: NM_001128223

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001128223.1</u> , <u>NP_001121695.1</u>
RefSeq Size:	3068 bp
RefSeq ORF:	2745 bp
Locus ID:	100131827
Cytogenetics:	3p12.3
Gene Summary:	This gene encodes a Kruppel-associated box (KRAB) zinc-finger protein, which belongs to a large group of transcriptional regulators in mammals. These proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation and apoptosis, and in regulating viral replication and transcription. A pseudogene of this gene was identified on chromosome 1. [provided by RefSeq, May 2016] Transcript Variant: This variant (1) encodes the longest isoform (a). Variants 1 and 2 encode the same isoform. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on paralogous alignments.