

## **Product datasheet for SC323403**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Casein Kinase 1 delta (CSNK1D) (NM\_139062) Human Untagged Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: Casein Kinase 1 delta (CSNK1D) (NM\_139062) Human Untagged Clone

Tag: Tag Free

Symbol: Casein Kinase 1 delta

Synonyms: ASPS; CKI-delta; CKId; CKIdelta; FASPS2; HCKID

Mammalian Cell None

Selection:

Vector: pCMV6-XL5

E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene ORF within SC323403 sequence for NM\_139062 edited (data generated by NextGen

Sequencing)

ATGGAGCTGAGAGTCGGGAACAGGTACCGGCTGGGCCGGAAGATCGGCAGCGGCTCCTTC GGAGACATCTATCTCGGTACGGACATTGCTGCAGGAGAAGAGGTTGCCATCATGCTTGAA TGTGTCAAAACCAAACACCCTCAGCTCCACATTGAGAGCAAAATCTACAAGATGATGCAG GGAGGAGTGGGCATCCCCACCATCAGATGGTGCGGGGCAGAGGGGGACTACAACGTCATG GTGATGGAGCTGCTGGGGCCAAGCCTGGAGGACCTCTTCAACTTCTGCTCCAGGAAATTC AGCCTCAAAACCGTCCTGCTGCTTGCTGACCAAATGATCAGTCGCATCGAATACATTCAT TCAAAGAACTTCATCCACCGGGATGTGAAGCCAGACAACTTCCTCATGGGCCTGGGGAAG AAGGGCAACCTGGTGTACATCATCGACTTCGGGCTGGCCAAGAAGTACCGGGATGCACGC ACCCACCAGCACATCCCCTATCGTGAGAACAAGAACCTCACGGGGACGGCGCGGTACGCC TCCATCAACACGCACCTTGGAATTGAACAATCCCGAAGAGATGACTTGGAGTCTCTGGGC TACGTGCTAATGTACTTCAACCTGGGCTCTCTCCCCTGGCAGGGGCTGAAGGCTGCCACC AAGAGACAGAAATACGAAAGGATTAGCGAGAAGAAAATGTCCACCCCATCGAAGTGTTG TGTAAAGGCTACCCTTCCGAATTTGCCACATACCTGAATTTCTGCCGTTCCTTGCGTTTT GACGACAAGCCTGACTACTCGTACCTGCGGCAGCTTTTCCGGAATCTGTTCCATCGCCAG GGCTTCTCCTATGACTACGTGTTCGACTGGAACATGCTCAAATTTGGTGCCAGCCGGGCC GCCGATGACGCCGAGCGGGAGCGCAGGGACCGAGAGGAGCGGCTGAGACACTCGCGGAAC CCGGCTACCCGCGGCCTCCCTTCCACAGCCTCCGGCCGCCTGCGGGGGACGCAGGAAGTG GCTCCCCCACACCCCTCACCCCTACCTCACACACGGCTAACACCTCCCCCGGCCCGTC TCCGGCATGGAGAGAGAGGGGAAAGTGAGTATGCGGCTGCACCGCGGGGCCCCCGTCAAC ATCTCCTCGTCCGACCTCACAGGCCGACAAGATACCTCTCGCATGTCCACCTCACAGAAT AGCATTCCTTTCGAACACCACGGCAAGTAG

Clone variation with respect to NM\_139062.2 113 a=>t





5' Read Nucleotide Sequence: Kinase Domain Sequence:

**OTI Disclaimer:** 

>SC323403 kinase domain raw sequence. By performing <u>BLASTX</u> analysis with this sequence against NCBI refernce protein database, you can confirm the presence of the kinase-

deficient mutation

CSAATGMGCAATGGGCGGTAGGCGTGTACGGTGGAGGTCTATATAAGCAGAGCTCGTTTAGTGAACCGT CAGAATTTTGTAATACGACTCACTATAGGGCGGCCGCGAATTCGGCACGAGGCGGCTCCTTTAGGCAGCT GAAAGGGGATTTAGGCCCGGAAGATCCGAGTCCATCCGCGGCGGGGAGAGGGCAAGCGGGACCGGTAGGG

GCCGGAGCAGCGGCGCGCGCTCGGACTGTCCCATCCGCCCCGT

Restriction Sites: Please inquire
ACCN: NM\_139062
Insert Size: 2500 bp

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:customer.com">customer.com</a> are team at <a href="mailto:customer.com">customer.com</a> or by

calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

This kinase-deficient mutant clone was generated by created by site-directed mutagenesis from the corresponding wild-type clone. See details in "Application of active and kinase-deficient kinome collection for identification of kinases regulating hedgehog signaling." <u>Cell.</u>

2008 May p536-548.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Components

**OTI Annotation:** 



**Reconstitution Method:** 

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM 139062.1, NP 620693.1

 RefSeq Size:
 2106 bp

 RefSeq ORF:
 1230 bp

 Locus ID:
 1453

 UniProt ID:
 P48730

 Cytogenetics:
 17q25.3

**Domains:** pkinase, TyrKc, S\_TKc

**Protein Families:** Druggable Genome, Protein Kinase

Protein Pathways: Circadian rhythm - mammal, Gap junction, Hedgehog signaling pathway

**Gene Summary:** This gene is a member of the casein kinase I (CKI) gene family whose members have been

implicated in the control of cytoplasmic and nuclear processes, including DNA replication and repair. The encoded protein may also be involved in the regulation of apoptosis, circadian rhythm, microtubule dynamics, chromosome segregation, and p53-mediated effects on growth. The encoded protein is highly similar to the mouse and rat CK1 delta homologs. Three transcript variants encoding different isoforms have been found for this gene.

[provided by RefSeg, Feb 2014]

Transcript Variant: This variant (2) has an additional exon at the 3' end compared to transcript

variant 1. This results in a shorter isoform (2) with a different C-terminus compared to

isoform 1.