

Product datasheet for **SC323655**

MSK1 (RPS6KA5) (NM_004755) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MSK1 (RPS6KA5) (NM_004755) Human Untagged Clone
Tag:	Tag Free
Symbol:	MSK1
Synonyms:	MSK1; MSPK1; RLPK
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_004755 edited
GAGCGGAAAGAGACTCGTCTTTGCGTCCGAGTTCTGGAGCCCGCCGACCCCGACTCCTG
GGGCCGCGGCAGCGGCTGCGAGGGGACGGGCGTCCGCTGTCTCCTGGGTTCCCTCGTAG
CGACCCGCGGGATCGGAAAAAAGGAGAAGATGGAGGAGGAGGTGGCAGCAGCGCGGC
GCCGCGGGGACCAGCGCGGACGGCGGACGGAGGAGAGCAGCTCCTCACTGTCAAGCAC
GAGCTGCGGACTGCTAATTTGACAGGACATGCTGAGAAGGTGGGAATAGAAAATTTGAG
CTCCTGAAGGTCCTAGGAACTGGAGCTTATGAAAAAGTATTTCTAGTTTCGTAATAAGT
GGCCATGATACTGGAAAGCTGTATGCCATGATGGTTTTGAAAAAGGCAACAATCGTTCAA
AAGGCCAAAACCACAGAGCATACAAGGACAGAACGACAAGTCTGGAACACATTAGGCAG
TCGCCATTTTGGTAACATTACATTATGCTTTCCAGACAGAAACCAAACCTCATCTCATT
TTAGATTATATAAATGGTGGTGAACTTTTACTCATCTTTCTCAAAGAGAGCGTTTCACA
GAGCATGAGGTGCAGATTTATGTTGGAGAGATTGTGCTTGCCCTTGAACATCTCCACAAG
TTGGGGATTATATATCGTGATATTAAGCTTGAGAATATTCTACTTGATTCTAATGGCCAT
GTGGTGTGACAGATTTTGGTCTGAGTAAGGAGTTTGTGGCTGATGAACTGAAAGAGCA
TATTCCTTTTGTGGAACATTGAATACATGGCACCAGATATTGTCAGAGGGGGAGATTCA
GGACATGACAAGGCAGTTGACTGGTGGAGTTGGGTGTTCTAATGTATGAATTACTAACT
GGAGCATCTCCTTTCAGTGTGATGGAGAAAAAATCCCAAGCTGAGATATCTAGGAGA
ATATTA AAAAGTGAGCCTCCATATCCCCAAGAAATGAGTGCTTTAGCGAAAGACCTAATT
CAGCGCTTTTGGATGAAAGATCCCAAGAAGAGATTGGGATGTGGTCCACGTGATGCAGAT
GAAATCAAAGAACATCTCTTTCCAGAAAAATAAATTGGGATGATTTAGCCGCAAAAAA
GTGCCGTGACCATTTAAGCCAGTCATTCGAGATGAATTAGATGTGAGTAACCTTTCAGAA
GAGTTCACAGAAATGGATCCCACTTATTCTCCCGAGCCCTGCCAGAGTTCTGAGAAG
CTGTTTCAGGGCTATTCTTTGTTGCTCCTCCATCCTATTCAAGCGTAATGCAGCTGTC
ATAGACCCTCTTCAGTTTCACATGGGAGTTGAACGTCCTGGAGTGACAAATGTTGCCAGG
AGTGCAATGATGAAGGACTCTCCATTCTATCAACACTATGACCTAGATTTGAAGGACAAA
CCCCTGGGAGAAGGTAGTTTTTCAATTTGTCGAAAGTGTGTGCATAAAAAAGTAACCAA
GCTTTTGCAGTCATGATAATCAGCAAAAAGGATGGAAGCCAATACTCAAAGGAAATAACA



[View online »](#)

```
GCTCTGAAACTCTGTGAAGGACACCCCAATATTGTGAAGTTGCATGAAGTTTTTCATGAT
CAGCTTCACACGTTTCTAGTGATGGAACCTCTGAATGGAGGAGAACTGTTTGAGCGCATT
AAGAAAAAGAAGCACTTCAGTGAGACGGAAGCCAGCTACATCATGAGGAAGCTTGTTTCA
GCTGTAAGCCACATGCATGATGTTGGAGTGGTGCACAGGGATCTGAAACCTGAGAATTTA
TTGTTCCACCGATGAAAATGACAAATTTGGAAATTTAAAATAATTGATTTTGGATTTGCACGG
CTAAAGCCACCGGATAATCAGCCCCTGAAGACTCCATGCTTACCCTTCATTATGCCGCC
CCAGAGCTCTTGAATCAGAACGGCTACGATGAGTCTGTGACCTGTGGAGCTTGGCGCTC
ATTTTGTACACAATGTTGTCAGGACAGGTTCCCTTCCAATCTCATGACCGAAGTTTGACG
TGTACCAGCGCGGTGAAATCATGAAGAAAATTTAAAAGGGAGATTTCTCCTTTGAAGGA
GAAGCCTGGAAGAATGTATCCCAAGAGGCTAAAGATTTGATCCAAGGACTTCTCACAGTA
GATCCAAACAAAAGGCTTAAAATGTCTGGCTTGAGGTACAATGAATGGCTACAAGATGGA
AGTCAGCTGTCTCCAATCCTCTGATGACTCCGGATATTCTAGGATCTTCCGGAGCTGCC
GTGCATACCTGTGTGAAAGCAACCTTCCACGCCTTTAACAATAACAAGAGAGAGGGGTTT
TGCCCTCAGAATGTTGATAAGGCCCTTTGGCTAAGAGAAGAAAATGAAAAGACTAGC
ACCAAGTACCGAGACGCGCAGCAGTTCCAGTGAGAGTTCCATTCTTCTCTCTCATTCT
CACGGTAAAACACCCACCAAGACACTGCAGCCAGCAATCCTGCCGACAGCAATAAC
CCGGAGACCCTCTTCCAGTTCTCGGACTCAGTAGCTTAGGCATGGTAGGAGTGTATCAGT
GATCCATTGCACCTTTATCCCTCAGCATATGCCTGAGGCGATCTTTATGCTTTTAAAA
ATGTTTCCCGTTGGTCTCATTGGAATCTGCCTCCTAATGATTTTTTTTCAGGAAAACCTGT
TTGGTTATCCTCATTCAAAGCACTGGACAGAGAATGTTACTGTGAATAGAGCACATATT
ACTCTTTTAGCAACCTAGCATGATGCCAACAAGACTATTCTTGAAGAGCAAAGTTCC
TGTAATTTAATTAGGGCTAGATTTGAGCTGCTGTAAGTACAGGTTTTCCAGATGTCT
GCCAACAAGAAATGACTCATACTGTGATGATACCTTTTGCCTTGTGGACAATGTG
GGTTTTTGAAATTTGCACCTTCAAACAATGATTTATCAGAGAAAAGGGGTCTGTTTTCAA
AAAAGATTCTGTAATGAATTTTATGTGTGCATATACTTATTTCTTGAGAGAAGATTTTA
ACTTATTGTTTTTATTTTATGGTTACATAT
```

- Restriction Sites:** Please inquire
- ACCN:** NM_004755
- Insert Size:** 3100 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004755.2](#), [NP_004746.2](#)
RefSeq Size: 3883 bp
RefSeq ORF: 2409 bp
Locus ID: 9252
UniProt ID: [O75582](#)
Cytogenetics: 14q32.11
Domains: pkinase, S_TK_X, TyrKc, S_TKc
Protein Families: Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways: Bladder cancer, MAPK signaling pathway, Neurotrophin signaling pathway

Gene Summary:

Serine/threonine-protein kinase that is required for the mitogen or stress-induced phosphorylation of the transcription factors CREB1 and ATF1 and for the regulation of the transcription factors RELA, STAT3 and ETV1/ER81, and that contributes to gene activation by histone phosphorylation and functions in the regulation of inflammatory genes (PubMed:11909979, PubMed:12569367, PubMed:12763138, PubMed:9687510, PubMed:18511904, PubMed:9873047). Phosphorylates CREB1 and ATF1 in response to mitogenic or stress stimuli such as UV-C irradiation, epidermal growth factor (EGF) and anisomycin (PubMed:11909979, PubMed:9873047). Plays an essential role in the control of RELA transcriptional activity in response to TNF and upon glucocorticoid, associates in the cytoplasm with the glucocorticoid receptor NR3C1 and contributes to RELA inhibition and repression of inflammatory gene expression (PubMed:12628924, PubMed:18511904). In skeletal myoblasts is required for phosphorylation of RELA at 'Ser-276' during oxidative stress (PubMed:12628924). In erythropoietin-stimulated cells, is necessary for the 'Ser-727' phosphorylation of STAT3 and regulation of its transcriptional potential (PubMed:12763138). Phosphorylates ETV1/ER81 at 'Ser-191' and 'Ser-216', and thereby regulates its ability to stimulate transcription, which may be important during development and breast tumor formation (PubMed:12569367). Directly represses transcription via phosphorylation of 'Ser-1' of histone H2A (PubMed:15010469). Phosphorylates 'Ser-10' of histone H3 in response to mitogenics, stress stimuli and EGF, which results in the transcriptional activation of several immediate early genes, including proto-oncogenes c-fos/FOS and c-jun/JUN (PubMed:12773393). May also phosphorylate 'Ser-28' of histone H3 (PubMed:12773393). Mediates the mitogen- and stress-induced phosphorylation of high mobility group protein 1 (HMG1/HMG14) (PubMed:12773393). In lipopolysaccharide-stimulated primary macrophages, acts downstream of the Toll-like receptor TLR4 to limit the production of pro-inflammatory cytokines (By similarity). Functions probably by inducing transcription of the MAP kinase phosphatase DUSP1 and the anti-inflammatory cytokine interleukin 10 (IL10), via CREB1 and ATF1 transcription factors (By similarity). Plays a role in neuronal cell death by mediating the downstream effects of excitotoxic injury (By similarity). Phosphorylates TRIM7 at 'Ser-107' in response to growth factor signaling via the MEK/ERK pathway, thereby stimulating its ubiquitin ligase activity (PubMed:25851810).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (a). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.