

Product datasheet for **SC323721**

PEX26 (NM_017929) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: PEX26 (NM_017929) Human Untagged Clone
Tag: Tag Free
Symbol: PEX26
Synonyms: PBD7A; PBD7B; PEX26M1T; Pex26pM1T
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC (PS100020)
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_017929.3
GGCTGGGCGAGCGCAAAGATGTCGCGCCCGCTGCCGGGAGGCGAGGCCAACTCGGGATA
TCCCGGAGCCTCTGGGGAGGCGGTCACTCCGACGCTGAGGACCTGGGCCTTGGACCCGG
ACTCGTTATGAAGAGCGATTCTTCGACCTCTGCAGCCCCCTCAGGGGGCTCGGGGGACC
CCTGCGCAGCAGCGAGCCGGTGC GCGCGGTCCCGGCCGGGCGCCGGCCTGGACCTTCT
GGAGGAGGCGGCCGACCTCCTGGTGGTGCACCTGGACTCCGGGCGGCGTGGAGACCTG
CGAGCGGGCCTGGCAGAGTCTGGCCAACCACGCCGTGGCAGAGGAACCCGCGGGCACCTC
ATTGGAGGTGAAGTGTCCCTGTGTGTTGTGGGATCCAGGCCCTGGCAGAAATGGATCG
GTGGCAAGAAAGTCTCCTGGGTCCTCAGTATTACCAGGTCCCTGAAAAGCTACCCCC
CAAAGTCTGGAGCTGTGCATTCTTTATACAGCAAAATGCAAGAGCCTGGAGCTGTGCT
GGATGTGGTGGGTGCCTGGCTCCAAGACCCAGCCAATCAAACCTTCCAGAATATGGAGC
CTTGGCAGAAATTCACGTGCAGCGGGTGTGCTGCCTCTGGGCTGCTTATCGGAGGCTGA
GGAGCTAGTGGTGGGCTCTGCAGCCTTTGGTGAGGAGCGCGACTGGATGTACTTCAGGC
CATTACACACAGCGAGGCAGCAGCAGAAACAGGAACACTCAGGCTCTGAGGAGGCCAGAA
GCCAAACCTGGAAGGCTCTGTCTCCACAAGTTCCTGTCACTACCGATGTTGGTTCGCCA
GCTTTGGGACTCTGCGGTGAGCCACTTCTTTCTGCCCCTCAAAGAGTCTCCTGGC
TGCCTTGATCCTCTGTCTCCTGGTGGTGAATTTGATCCAGCTTCCCCTCCTCCCTGCA
CTTCTCTACAAGCTGGCCAGCTCTCCGCTGGATCCGGAAGGCTGCATTTTCTCGCT
CTACCAGTCCGCATCCATGACTGAGGGTCCCTGCGCACACAGCCTCTCTGCTCCTCAC
GTCCGTGGCCACAGAAGCAGAGCGACAGAGCGACACATCCACAGGCGCCCTGGGGAAAT
GGGACCAGCCTAATCTCGCGAGTGCACCTGTGTCTTGCTGCCTGGGTGCCCTCTCCTTTG
CACCCTACTTCGGCTGGTGC GGTAGATGATGTGGAAACAAAGCAGGACCAGCAGGCACA
GCACCTCCAGAACAGTGC CCGGATGACCAGAGGCCCTTGAAAAGGAGGGTTTGGGGA
CAGGGACTGTGTCCATGAAACATTCCATCTTCTGGTGAAGGCAAGGGTTGGTTCTTCA
GGTCAGGATGTTAATGGAGCTGGAAGTTCAGAAAAAGCCTGGTGAAGTGACAAAAAAA
AAAAA



[View online »](#)

Restriction Sites:	ECoRI-NOT
ACCN:	NM_017929
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017929.3 , NP_060399.1
RefSeq Size:	2474 bp
RefSeq ORF:	918 bp
Locus ID:	55670
UniProt ID:	Q7Z412
Cytogenetics:	22q11.21

Gene Summary:

This gene belongs to the peroxin-26 gene family. It is probably required for protein import into peroxisomes. It anchors PEX1 and PEX6 to peroxisome membranes, possibly to form heteromeric AAA ATPase complexes required for the import of proteins into peroxisomes. Defects in this gene are the cause of peroxisome biogenesis disorder complementation group 8 (PBD-CG8). PBD refers to a group of peroxisomal disorders arising from a failure of protein import into the peroxisomal membrane or matrix. The PBD group is comprised of four disorders: Zellweger syndrome (ZWS), neonatal adrenoleukodystrophy (NALD), infantile Refsum disease (IRD), and classical rhizomelic chondrodysplasia punctata (RCDP). Alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Dec 2010]

Transcript Variant: This variant represents transcript variant 1. Both variants 1 and 2 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.